


MOVES[®] SLC[™]

PORTABLE LIFE SUPPORT SYSTEM



Operator's Manual

Contact Information

MANUFACTURER:

✉ **Address:** Thornhill Research Inc.
210 Dundas Street West, Suite 200
Toronto, Ontario
Canada M5G 2E8

☎ **Phone:** +1.416.597.1325 or +1.888.597.1325

☎ **Fax:** +1.416.597.1330

✉ **Email:** info@thornhillresearch.com

🌐 **Website:** <http://www.thornhillresearch.com>

NOTIFIED BODY:

British Standards Institute
Products Services
Maylands Avenue, Hertfordshire
HP2 4SQ United Kingdom
Phone: +44 0.1442.230442

AUTHORIZED REPRESENTATIVES:

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Phone: +31.70.345.8570
Fax: +31.70.346.7299

RQSOLUTIONS
1F/1 Francis Road, Artarmon
NSW 2064 Australia
Phone: + 61.2.9413.2850
Fax: +61.2.9413.2852

THIS PAGE DELIBERATELY LEFT BLANK.

Tale of Contents

1.0	Notices	1
1.1	Patent Notice	1
1.2	Copyright Notice	1
1.3	Trademark Notices	1
1.4	EMC (Electromagnetic Compatibility) Notice	1
1.5	MRI (Magnetic Resonance Imaging) Notice	1
1.6	Regulatory Notice	1
1.7	Classification	2
1.8	Declaration of Conformity Notice	2
1.8.1	Trade Name	2
2.0	MOVES® SLC™ Quick Start Ventilate Guide	3
2.1	Further Information on Procedures	10
3.0	Introduction	13
3.1	General Information	13
3.1.1	Manufacturer's Information	13
3.1.2	Glossary of Terms and Abbreviations	13
3.2	MOVES® SLC™ Intended Use	15
3.2.1	Operating Environment	15
3.2.2	Intended Operator	15
3.2.3	Intended Population	15
4.0	Regulatory Compliance	17
4.1	Regulatory Symbols	17
4.2	Regulatory Standards Compliance	18
5.0	MOVES® SLC™ System Overview	19
5.1	General Overview	19
5.2	System Orientation and External Components	19
5.3	Theory of Operation	20
5.3.1	Breathing Circuit and Oxygen Supplement	20
5.3.2	Oxygen Concentrator	20
5.3.3	Ventilator	20
5.3.4	Respiratory Gas Monitoring	21
5.3.5	Suction	22
5.3.6	Patient Monitoring	23
5.3.7	Power System	24
6.0	The Masimo Rainbow SET® Pulse CO-Oximeter	25
6.1	Overview	25
6.2	Key Features	25
6.3	Indications for Use	25
6.4	Pulse Oximeter Technology Overview	25
6.4.1	Signal Extraction Technology (SET®)	25
6.4.2	General Description for Oxygen Saturation (SpO ₂)	26
6.4.3	Successful Monitoring for SpO ₂ , PR, and PI	26
6.4.4	Functional Oxygen Saturation	26

6.4.5	General Description for Pulse Rate (PR)	26
6.4.6	General Description for Perfusion Index (PI)	26
6.4.7	Rainbow Pulse CO-Oximetry Technology	26
6.4.8	Pulse CO-Oximetry vs. Drawn Whole Blood Measurements	28
6.4.9	General Description for Total Hemoglobin (SpHb)	28
6.4.10	Successful Monitoring for SpHb	28
6.4.11	General Description for SpOC	29
6.4.12	General Description for Carboxyhemoglobin (SpCO)	29
6.4.13	Successful Monitoring for SpCO	29
6.4.14	General Description for Methemoglobin (SpMet)	29
6.4.15	Successful Monitoring for SpMet	29
6.4.16	SpCO, SpMet, and SpHb Measurements During Patient Motion	29
7.0	Battery Storage	31
8.0	Safety Information	33
8.1	Manual Symbols and Messages	33
8.1.1	Notes, Cautions and Warnings	33
8.1.2	Operational Symbols	34
8.1.3	Label Warning Symbols	34
8.1.4	Product Label Symbols	34
8.2	General Warnings	37
8.3	Electrical Warnings	39
8.4	Patient-Specific Warnings	39
8.5	Masimo Rainbow SET [®] Pulse CO-Oximeter Warnings	41
8.6	General Cautions	45
8.7	Electrical Cautions	48
8.8	Patient-Specific Cautions	48
8.9	General Safety	49
8.10	Electrical Safety	50
8.11	Preparing for Emergency Operation	50
8.12	Radio Interference	51
8.13	Battery Handling	51
8.14	Battery Disposal	51
9.0	Getting Started	53
9.1	MOVES [®] SLC [™] System Contents	53
9.1.1	MOVES [®] SLC [™] Accessories Warnings and Cautions	53
9.1.2	System Contents	54
9.2	MOVES [®] SLC [™] Consumables	65
9.3	Attaching the Shoulder Strap to MOVES [®] SLC [™]	68
9.4	Attaching Clamps to the MOVES [®] SLC [™]	70
9.5	Patient Monitoring Accessories	73
9.6	Installing the Ventilator Cartridge and Breathing Circuit	75
9.6.1	About the Ventilator Cartridge	75
9.6.2	Ventilator Cartridge Installation Instructions	77
9.6.3	Installing the Hydrocarbon Filter	79
9.6.4	Installing the Ventilator Breathing Circuit	80
9.7	Delivering Supplementary Oxygen (O ₂)	83

9.8	Installing Suction Accessories	85
9.9	Preparing MOVES [®] SLC [™] for Activation	90
9.9.1	Checking Battery Charge	90
9.9.2	Inspecting the Batteries	91
9.9.3	Installing the Batteries	92
9.9.4	Preparing the Power Supply / Battery Charger	95
9.9.5	Connecting MOVES [®] SLC [™] to the Power Supply / Battery Charger	96
9.9.6	Connecting AC Power	97
10.0	Startup	99
10.1	Adjusting the Screen Display Orientation	99
10.2	User Interface (UI) Controls and Functions	102
10.3	System Visual Indicators (Alarms / System Status)	104
10.4	Changing Settings and Data Views	105
10.4.1	Display of Settings and Views	105
10.4.2	Modifying a Setting	106
10.5	Confirming Display Validity	106
10.6	Startup Sequence	106
11.0	MOVES [®] SLC [™] User Screens	111
11.1	Status Bar	111
11.1.1	Battery Status Icon	113
11.2	System Test Screen	115
11.3	Setup Screen	116
11.3.1	Overview	116
11.3.2	Changing Settings	116
11.3.3	Setup Screen – Ventilate Mode	117
11.3.4	Setup Screen – O2 Supplement	118
11.3.5	Setup Screen – Monitor Only	118
11.3.6	Setup Screen Options	119
11.4	Advanced Screen	122
11.4.1	Accessing the Advanced Screen	122
11.4.2	The Advanced Screen Itself	123
11.4.3	SpO2 Avg (Average) Time	124
11.4.4	SpO2 Sensitivity Mode	125
11.4.5	SpHb Average Window	125
11.4.6	ECG Electromyography (EMG) Filter	125
11.4.7	Additional Information Located on Advanced Screen	126
11.5	Main Screen	126
11.5.1	Selection Order for Options	128
11.5.2	Main Screen Items	130
11.5.3	PC, PEEP and PIP	135
11.5.4	Inverted Display of Patient Monitoring Values	135
11.5.5	Invasive Pressure (IP) Source Buttons	136
11.5.6	Zeroing a Channel	136
11.5.7	Graphs	136
11.6	ECG Screen	136
11.6.1	Overview	136

11.7	Alarm Limits Screen	138
11.7.1	Overview	138
11.7.2	Limit List Active	139
11.8	Alarm On / Off Screen	140
11.8.1	Overview	140
	Alarm On/Off List Active	141
12.0	Using the Remote Screen	143
12.1	Installing the Remote Screen Batteries	143
12.2	Battery Indicator	147
12.3	Connecting the Remote Screen to Wall Power	149
12.4	Connecting the Remote Screen to MOVES® SLC™	151
12.5	The Remote Screen User Interface (UI)	154
12.5.1	When First Connecting	154
12.5.2	No System Test Screen	155
12.5.3	Remote Screen Panel Buttons	157
12.5.4	Graphs Independently Configurable	158
12.5.5	Alarm Indicators	158
12.5.6	Remote Screen Software Version	160
13.0	Connecting the Patient	161
13.1	Connection Overview	161
13.2	Connecting an Intubated Patient	161
13.3	Connecting a Spontaneously Breathing Patient	161
13.4	Attaching the Ear or Finger Clip Pulse CO-Oximeter Sensor	161
13.5	Placement of the Heart Rate Electrodes	163
13.5.1	ECG Cable Color Coding and Naming Conventions	164
13.5.2	Correct Electrode Placement	164
13.5.3	Reducing Artifacts	166
13.5.4	Expiry Date of ECG Electrodes	166
13.6	Zeroing the Pressure in the IP Transducer	167
13.6.1	Zeroing the Transducer Channel	167
13.6.2	Transducer Warnings	167
13.7	Using the Suction Feature	168
13.8	Shutdown Procedures	168
14.0	Using System Graphs / Trends	171
14.1	Overview	171
14.2	Available System Graphs / Trends	172
14.2.1	System Graphs	172
14.2.2	System Trends	173
15.0	Alarms	175
15.1	About Status LEDs	175
15.2	Alarm Priorities and Characteristics	175
15.2.1	Standard Alarms	175
15.2.2	High Priority Communication Failure Alarm	176
15.3	Alarm Queue	176
15.4	Locked Alarms and Messages	176
15.5	Inhibitable and Latching Alarms	177

15.5.1	Inhibitable Alarms	177
15.5.2	Latching Alarms.....	178
15.5.3	Alarms That Have Been Turned Off.....	178
15.6	Alarms Conditions and Causes.....	179
15.7	Safe Gas Mode	205
16.0	Appendix	207
16.1	System Default Settings.....	207
16.2	System Cleaning	208
16.3	System Maintenance.....	208
16.3.1	Replacing Filters.....	208
16.3.2	General Maintenance	210
16.4	Accessories Maintenance	211
16.5	Checking the Accuracy of the Temperature Probe	211
16.6	MOVES [®] SLC [™] Specifications.....	211
16.6.1	Model Number.....	211
16.6.2	Physical Properties.....	211
16.6.3	Oxygen Concentrator Specifications	212
16.6.4	Ventilator Specifications	212
16.6.5	Suction Specifications	213
16.6.6	Electrical Characteristics	213
16.6.7	Environmental Specifications	214
16.7	Patient Monitoring Specifications	214
16.7.1	Heart Rate Monitoring Specifications.....	214
16.7.2	Temperature Monitoring Specifications.....	215
16.7.3	Airflow Monitoring Specifications	215
16.7.4	CO ₂ Monitoring Specifications.....	215
16.7.5	O ₂ Monitoring Specifications	216
16.7.6	ECG Specifications	216
16.7.7	NIBP Specifications.....	217
16.7.8	Invasive Pressure Specifications	217
16.7.9	Pulse Oximetry Specifications.....	218
16.7.10	Drift in Sensing Accuracy	218
16.7.11	Specifications of Masimo Rainbow SET [®] Pulse CO-Oximeter	218
16.8	General Accessories Specifications.....	220
16.8.1	Approved Masimo Pulse Oximeter Accessories	221
16.8.2	Masimo Pulse Oximeter Accessories Specifications.....	222
16.9	Electromagnetic Conformity Information.....	223
16.9.1	IEC 60601-1-2:2007 (Ed 3.0) Table 1 Requirements	223
16.9.2	IEC 60601-1-2:2007 (Ed 3.0) Table 202 Requirements	223
16.9.3	IEC 60601-1-2:200 (Ed 3.0) Table 203 Requirements	225
16.9.4	IEC 60601-1-2:2007 (Ed 3.0) Table 5 Requirements	226

List of Tables

Table 1: Quick Reference to Information on Procedures	10
Table 2: Glossary of Terms and Abbreviations	13
Table 3: Regulatory Symbols Used and Description	17
Table 4: Regulatory Standards Compliance	18
Table 5: Symbols and Messages Used in Manual	33
Table 6: Operational Symbols and Descriptions	34
Table 7: Label Warning Symbols and Descriptions	34
Table 8: Product Label Symbols and Descriptions	34
Table 9: General Warnings	37
Table 10: Electrical Warnings	39
Table 11: Patient-Specific Warnings	39
Table 12: Masimo Rainbow SET [®] Pulse CO-Oximeter Warnings	41
Table 13: General Cautions	45
Table 14: Electrical Cautions	48
Table 15: Patient-Specific Cautions	48
Table 16: MOVES [®] SLC [™] Accessories Warnings and Cautions	53
Table 17: MOVES [®] SLC [™] System Contents	54
Table 18: MOVES [®] SLC [™] Consumables	65
Table 19: Patient Connector Labels and Accessories	74
Table 20: User Interface Controls and Functions	102
Table 21: System Visual Indicator States and Explanations	105
Table 22: Screen Buttons and Descriptions	105
Table 23: MOVES [®] SLC [™] Status Bar Items and Descriptions	111
Table 24: Battery Status Icon Table	114
Table 25: Setup Screen Options and Descriptions	119
Table 26: Advanced Screen Options	124
Table 27: Advanced Screen Additional Info Items	126
Table 28: Main Screen Items and Descriptions	130
Table 29: Alarm Limits and Defaults	139
Table 30: Remote Screen Battery Status Indicators	148
Table 31: Alarm Side Bar States and Explanations	158
Table 32: CO-Oximeter Sensors	161
Table 33: ECG Cable Color Coding and Naming	164
Table 34: Proper Placement of Precordial Leads	166
Table 35: System Graphs and Parameters	172
Table 36: System Trends and Parameters	174
Table 37: Alarm Types and Descriptions	175
Table 38: Alarm Conditions and Causes	179
Table 39: System Default Settings	207
Table 40: Physical Properties of MOVES [®] SLC [™]	211
Table 41: Oxygen Concentrator Specifications of MOVES [®] SLC [™]	212
Table 42: Ventilator Specifications of MOVES [®] SLC [™]	212
Table 43: Suction Specifications of MOVES [®] SLC [™]	213
Table 44: Electrical Specifications of MOVES [®] SLC [™]	213

Table 45: Environmental Specifications of MOVES [®] SLC [™]	214
Table 46: Heart Rate Monitoring Specifications of MOVES [®] SLC [™]	214
Table 47: Temperature Monitoring Specifications of MOVES [®] SLC [™]	215
Table 48: Airflow Monitoring Specifications of MOVES [®] SLC [™]	215
Table 49: CO ₂ Monitoring Specifications of MOVES [®] SLC [™]	215
Table 50: O ₂ Monitoring Specifications of MOVES [®] SLC [™]	216
Table 51: ECG Specifications of MOVES [®] SLC [™]	216
Table 52: NIBP Specifications of MOVES [®] SLC [™]	217
Table 53: Invasive Pressure Specifications of MOVES [®] SLC [™]	217
Table 54: Pulse Oximetry Specifications of MOVES [®] SLC [™]	218
Table 55: Pulse Oximetry Equipment Response Time.....	218
Table 56: General Accessories Specifications of MOVES [®] SLC [™]	220
Table 57: Approved Masimo Pulse Oximeter Accessories	221
Table 58: Masimo Pulse Oximeter Accessories Specifications	222
Table 59: 15.9.1 IEC 6060-1-2:2001 (Ed 2) Table 201 Requirements	223
Table 60: 5.2.2.1f IEC 60601-1-2:2007 (Ed 3.0) Table 2 Requirements.....	223
Table 61: 5.2.2.2 IEC 60601-1-2:2007 (Ed 3.0) Table 3 Requirements	225
Table 62: 5.2.2.2 IEC 60601-1-2:2007 (Ed 3.0) Table 5 Requirements	226

Table of Figures

Figure 5-1: MOVES [®] SLC [™] Upright System Orientation and Components	20
Figure 6-1: Absorption Spectra of Blood Components	27
Figure 6-2: LEDs and Detector	27
Figure 7-1: Battery Compartment Latch Open.....	31
Figure 7-2: Battery Compartment Latch Fully Closed.....	31
Figure 9-1 : Patient Connection Panel	74
Figure 9-2: Ventilator Cartridge – Patient Connection Side (Front)	75
Figure 9-3: Ventilator Cartridge – MOVES [®] SLC [™] Connection Side (Back)	76
Figure 9-4: Ventilator Breathing Circuit.....	80
Figure 9-5: Battery with Condition Indicator Shown.....	90
Figure 9-6: Condition Indicator Showing High Level Charge.....	90
Figure 10-1: MOVES [®] SLC [™] Panel Buttons.....	101
Figure 10-2: MOVES [®] SLC [™] Main Screen	102
Figure 10-3: Power Control Button	106
Figure 10-4: New Patient Screen	107
Figure 10-5: System Test Screen	107
Figure 10-6: Setup Screen	109
Figure 11-1: Status Bar	111
Figure 11-2: System Test Screen	115
Figure 11-3: Configure for a New Patient Screen	116
Figure 11-4: Setup Screen – Ventilate Mode.....	117
Figure 11-5: Setup Screen – O2 Supplement.....	118
Figure 11-6: Setup Screen – Monitor Only	118
Figure 11-7: Accessing the Advanced Screen.....	122
Figure 11-8: The Advanced Screen	123
Figure 11-9: Advanced Screen – SpO2 Avg Time Selected.....	124
Figure 11-10: Screen Items Drawn in Reverse.....	127
Figure 11-11: Screen Items Displaying Dashes.....	127
Figure 11-12: Screen Items Displaying the Fault Icon	127
Figure 11-13: Main Screen with Dashes & Alarms	128
Figure 11-14: IP & Trends Graph Options & Sub-Options.....	129
Figure 11-15: PC & PEEP on Setup Screen.....	135
Figure 11-16: PIP & PEEP on Main Screen.....	135
Figure 11-17: Relations of PEEP to PIP to PC	135
Figure 11-18: Regular Temperature Display.....	136
Figure 11-19: Inverted Display of Temperature (Alarm Condition)	136
Figure 11-20: ECG Screen – Pause Mode	137
Figure 11-21: ECG Screen – Resume Mode	137
Figure 11-22: Alarm Limits Screen – All at Default	138
Figure 11-23: Alarm Limits Screen Active.....	139
Figure 11-24: Alarm On / Off Screen	140
Figure 11-25: Alarm On / Off Screen Active	141
Figure 12-1: MOVES [®] SLC [™] Remote Screen	143
Figure 12-2: Battery Bay 1 on Left Side	143

Figure 12-3: Battery Bay 2 on Right Side	143
Figure 12-4: Remote Screen Battery Indicator	147
Figure 12-5: Initial Screen – No Resources Loaded	154
Figure 12-6: Initial Screen – Resources Already Loaded	155
Figure 12-7: Setup Screen – Initializing.....	156
Figure 12-8: Setup Screen – Monitor Only Mode.....	156
Figure 12-9: Main Screen – On Remote Screen	157
Figure 12-10: Remote Screen – Panel Buttons at Bottom	157
Figure 12-11: Green Bars Indicating No Alarms Active.....	159
Figure 12-12: Red Bars Indicating High-Priority Alarm Active.....	159
Figure 12-13: Yellow Bars Indicating Low or Medium Priority Alarm Active	160
Figure 12-14: Advanced Screen Showing Remote Screen Version	160
Figure 13-1: 1 of 5 Possible Sensor Finger Placements	162
Figure 13-2: Correct Ear Lobe Sensor Placement	162
Figure 13-3: Color-Coded 12-Lead ECG Cable	164
Figure 13-4: Proper Limb Lead Electrode Placement	165
Figure 13-5: Precordial Leads: V1–V6	165
Figure 13-6: Electrode Package	166
Figure 13-7: Electrode Use By Date	166
Figure 13-8: Transducer Stopcock in Running Position (Horizontal)	167
Figure 13-9: Transducer Stopcock in Open-to-Ambient-Air Position (Vertical) Zero Transducer	167
Figure 13-10: Suction Wand Showing Control Vent.....	168
Figure 14-1: Plotting Graphs on Main Screen (detail)	171
Figure 15-1: No Alarms Message in Alarm Queue.....	176
Figure 15-2: Two Locked Alarms.....	177
Figure 15-3: Inhibitible Alarm Able To Be Dismissed.....	177
Figure 15-4: Inhibitible Alarm Cleared.....	177
Figure 15-5: Latched Alarm Able To Be Dismissed	178
Figure 15-6: Latched Alarm Cleared	178
Figure 15-7: Latched Alarm Able To Be Dismissed	178

THIS PAGE DELIBERATELY LEFT BLANK.

1.0 Notices

1.1 Patent Notice

MOVES[®] SLC[™] is covered by one or more US and international patents and patents pending.

1.2 Copyright Notice

Copyright © 2013 Thornhill Research Inc.

No part of this publication may be reproduced, translated into another language, stored in a retrieval system or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise without the prior written consent of Thornhill Research Inc. The information contained herein is subject to change without notice.

1.3 Trademark Notices

MOVES[®] is a registered trademark of Thornhill Research Inc.

SLC[™] is a trademark of Thornhill Research Inc.

All brand and product names mentioned herein are used for identification purposes only and are trademarks or registered trademarks of their respective holders.

1.4 EMC (Electromagnetic Compatibility) Notice

MOVES[®] SLC[™] generates, uses, and can radiate RF (radio frequency) energy. If it is not installed and used in accordance with the instructions in this manual, electromagnetic interference may result. MOVES[®] SLC[™] has been tested and found to comply with the limits set forth in IEC 60601-1-2 (identical to EN 60601-1-2) for medical products. These limits provide reasonable protection against electromagnetic interference in the intended use environments (e.g., hospitals) described in this manual.

1.5 MRI (Magnetic Resonance Imaging) Notice

MOVES[®] SLC[™] contains electromagnetic components whose operation can be affected by intense electromagnetic fields. Do not operate MOVES[®] SLC[™] in an MRI environment or in the vicinity of high-frequency surgical diathermy equipment. Electromagnetic interference could disrupt the operation of the ventilator.

1.6 Regulatory Notice

Federal law restricts the sale of MOVES[®] SLC[™] except by or on the order of a physician.

1.7 Classification

Type of Equipment: Medical Equipment
Electrical Classification: Class I, Type BF (unit)
Registration Classification: Class IIb
Type CF Defibrillation Proof (applied parts)
Adult Lung Ventilator

1.8 Declaration of Conformity Notice

MOVES[®] SLC[™] is declared to conform to the Medical Device Directive of the European community (re: Medical Device Directive 93/42/EEC). This is indicated by the CE Mark shown below.



MOVES[®] SLC[™] also conforms to the following Technical Standards:

IEC 60601-1	IEC 60601-2-49	ASTM E1112-00
IEC 60601-1-2	ISO 80601-2-55	BS EN 794-3 (2009)
ISO 80601-2-13	ISO 80601-2-61	MIL-STD-810G
IEC 60601-2-27	ANSI/AAMI EC13	JECETS
IEC 80601-2-30	ANSI/AAMI SP10	
IEC 60601-2-34		

For more information on these standards, see *Regulatory Standards Compliance on page 18*.

1.8.1 Trade Name

MOVES[®] SLC[™] portable life support system

2.0 MOVES[®] SLC[™] Quick Start Ventilate Guide

Equipment Setup Procedures

1. Open the valve block cover door inside the cartridge cavity and insert the ventilator breathing cartridge.
2. Install the hydrocarbon filter in the REAR panel.
3. Attach the ventilator breathing circuit to the cartridge. Connect the Nafion tubing to the sampling line.
4. Attach the suction canister and tubing to MOVES[®] SLC[™].
5. Connect any required patient monitoring accessories to MOVES[®] SLC[™]. **DO NOT CONNECT TO THE PATIENT AT THIS TIME.**
6. Verify the battery charge levels and insert both batteries.
7. Attach the MOVES[®] SLC[™] AC power unit (includes power supply, battery charger and cables) to the MOVES[®] SLC[™] if required to recharge.

Pre-use Test Instructions

1. Turn the MOVES[®] SLC[™] unit on by pressing the *Power Control Button* (on the REAR panel).



2. The *System Test* screen is displayed.

SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

B A

Vent Check: **PASS**

Open Y-Piece: **PREP REQD**

Closed Y-Piece: **WAITING**

Power Check: **PASS**

Battery Verify: **PASS**

Sensor Check: **PASS**

Sample Line: **WAITING**

Lights+Audio: **INPUT REQD**


Prep for Open Y-Piece test:

Disconnect the Y-piece and ensure that it is unobstructed.

SKIP REMAINING **Prep done** **Skip this test**

Pre-use Test Instructions

3. Some tests are *AUTOMATIC*. Others require preparation of the device into a certain configuration (i.e., *PREP REQD*) or require the user to give feedback (i.e., *INPUT REQD*). The first test selected, the *Open Y-Piece* test, is a *PREP REQD* test.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

B A

Vent Check: **PASS**

Open Y-Piece: **PREP REQD**

Closed Y-Piece: **WAITING**

Power Check: **PASS**


Battery Verify: **PASS**

Sensor Check: **PASS**

Sample Line: **WAITING**

Lights+Audio: **INPUT REQD**

Prep for Open Y-Piece test:




Disconnect the Y-piece and ensure that it is unobstructed.

SKIP REMAINING **Prep done** **Skip this test**

4. Perform the preparation steps indicated if the test is *PREP REQD* (for other tests that are *INPUT REQD*, respond with the *input requested*), then choose *Prep done* and run the test.

NOTE: In order to save the user time, certain tests can be run concurrently. For example, while the *Open Y-Piece* test is still running, the screen will advance to the *Lights + Audio* tests

NOTE: The system should not be used until all tests pass.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

B A

Vent Check: **PASS**

Open Y-Piece: **PREP REQD**

Closed Y-Piece: **WAITING**

Power Check: **PASS**


Battery Verify: **PASS**

Sensor Check: **PASS**

Sample Line: **WAITING**

Lights+Audio: **INPUT REQD**

Prep for Open Y-Piece test:



Disconnect the Y-piece and ensure that it is unobstructed.

SKIP REMAINING **Prep done** **Skip this test**

Pre-use Test Instructions

5. To initiate the *Lights + Audio* tests, user input is required. Start must be confirmed, and then the operator must continue to observe the device and answer the posed questions.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: PASS

Open Y-Piece: 50% DONE

Closed Y-Piece: WAITING

Power Check: PASS


Battery Verify: PASS

Sensor Check: PASS

Sample Line: PREP REQD

Lights+Audio: INPUT REQD

Alarm Lights+Audio Tests:



Press START to initiate tests of the alarm lights and audio.

SKIP REMAINING
Start
Skip this test

6. The first *Lights + Audio* test requires input from the user (*INPUT REQD* test). This test is trying to ensure that the high-priority alarm visual and audio indicators are working properly. The user must respond with Yes or No.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: PASS

Open Y-Piece: PASS

Closed Y-Piece: PREP REQD

Power Check: PASS


Battery Verify: PASS

Sensor Check: PASS

Sample Line: PREP REQD

Lights+Audio: INPUT REQD

Test 1 of 2: RED Lights+Audio



Is the device sounding an alarm melody, and are all four lights blinking red?

SKIP REMAINING
Yes
No

Pre-use Test Instructions

7. The second *Lights & Audio* test also requires input from the user (*INPUT REQD* test). This test is trying to ensure that the medium-priority alarm visual indicator is working properly. The user must respond with Yes or No.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: PASS

Open Y-Piece: PASS

Closed Y-Piece: PREP REQD

Power Check: PASS

Battery Verify: PASS

Sensor Check: PASS

Sample Line: PREP REQD

Lights+Audio: INPUT REQD

Test 2 of 2: AMBER Lights



Are all four lights blinking amber?

SKIP REMAINING

Yes

No

8. The Sample Line test requires user preparation (*PREP REQD* test). Follow the instructions on screen, perform the preparation required, and then choose *Prep done* to start the test.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: PASS

Open Y-Piece: PASS

Closed Y-Piece: PASS

Power Check: PASS


Battery Verify: PASS

Sensor Check: PASS

Sample Line: PREP REQD

Lights+Audio: PASS

Prep for Sample Line test:



Occlude the sample line by using the occluding cap connected to the nafion tube -- ensure the sample line has been unoccluded for 5 seconds first.


SKIP REMAINING

Prep done

Skip this test

Pre-use Test Instructions

9. The *Closed Y-Piece* test requires user preparation (*PREP REQD* test). Follow the instructions on screen then choose *Prep done* to start the test.




SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: **PASS**
 Open Y-Piece: **PASS**
 Closed Y-Piece: **PREP REQD**

Power Check: **PASS**
 Battery Verify: **PASS**
 Sensor Check: **PASS**
 Sample Line: **PREP REQD**
 Lights+Audio: **PASS**


Prep for Closed Y-Piece test:



Occlude the Y-piece by attaching it to the test port in the cartridge. Ensure that the sample port on the Y-piece is blocked.

SKIP REMAINING **Prep done** Skip this test

10. **NOTE:** When testing is in progress, and the MOVES[®] SLC[™] is running, a small white dot can be seen moving up and down at the top right of the screen to indicate that the device is not "frozen". The white "dot" is a way of knowing the device is not "frozen" at any time, not just during startup tests.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen

Vent Check: **PASS**
 Open Y-Piece: **50% DONE**
 Closed Y-Piece: **WAITING**

Power Check: **PASS**
 Battery Verify: **PASS**
 Sensor Check: **PASS**
 Sample Line: **PASS**
 Lights+Audio: **PASS**


Testing in progress...

SKIP REMAINING



Pre-use Test Instructions

11. If all tests are successful, the *System Test* screen displays a *Test Summary* indicating all tests have passed. Choose *Continue*.



SYSTEM TEST

WARNING: **DO NOT CONNECT PATIENT** during this screen


Vent Check: **PASS**
 Open Y-Piece: **PASS**
 Closed Y-Piece: **PASS**
 Power Check: **PASS**
 Battery Verify: **PASS**
 Sensor Check: **PASS**
 Sample Line: **PASS**
 Lights+Audio: **PASS**

TEST SUMMARY: PASSED

The system, as configured at time of testing, passed all tests. Press CONTINUE to acknowledge and proceed.

RETRY ALL TESTS CONTINUE

12. You are taken to the *Setup Screen* and placed in *Monitor Only* mode.



SETUP Monitor Only No Alarms

SYSTEM MODE: Ventilate O2 Supplement Monitor Only

VENTILATOR **CONCENTRATOR**

Vent Mode: IMV A/C SIMV PSV Vent O2: 40%

Breath Sensitivity: Normal Low

Control: Volume Pressure

Tidal Volume: 500 mL

PC: 20 cmH2O

PEEP: 0-3 cmH2O

RR: 10 B/M

I:E Ratio: 1:2

Inspire Time: 1.0 s

Pressure Support: disabled

OTHER

Line Filter: 50Hz 60Hz

ECG Leads: 3 5 12

ECG Rate: slow fast

ECG Range: 3mV 6mV 12mV

NIBP Mode: MANUAL


Suction: 325 mmHg

ADVANCED...


Pre-use Test Instructions

13. Make any changes necessary to the *Ventilate* mode settings and configure any other settings necessary. Also, connect any desired patient monitoring cables to the patient.

NOTE: *Ventilate* and other settings can be changed even though the system is *Monitor Only* mode.



14. Change the system to *Ventilate* mode. Ventilation begins.



15. Attach the patient to the breathing circuit.

Additional Buttons

Image	Name	Use
	Screen Dim Button	Use to adjust screen brightness.
	Alarm Audio Pause Button	Use either to silence ALL alarms permanently (if this feature is permitted) or for a temporary period (120 seconds).

	NIBP Control Button	Use to start or abort a Non-Invasive Blood Pressure (NIBP) reading.
	Suction Control Button	Use to activate MOVES [®] SLC [™] suction functionality. All suction accessories must be connected and ready to use before activating.

2.1 Further Information on Procedures

For further information on procedures, see the sections indicated in the following table.

Table 1: Quick Reference to Information on Procedures

PROCEDURE	Page
1. Attaching the shoulder strap.	68
2. Inserting the ventilator breathing cartridge.	77
3. Installing the hydrocarbon filter.	79
4. Installing the ventilator breathing circuit.	80
5. Inspecting the batteries.	91
6. Installing the batteries.	92
7. Preparing the power supply / battery charger	95
8. Connecting MOVES [®] SLC [™] to the power supply / battery charger.	96
9. Connecting AC power.	97
10. Adjusting the screen display orientation.	99
11. Changing the brightness of the display (see third item in table).	102
12. Performing System Tests.	115
13. Attaching the suction tube and suction canister to the MOVES [®] SLC [™] .	85
14. Making patient monitoring connections between the MOVES [®] SLC [™] unit and the patient.	73
15. Connecting ABP, CVP or ICP equipment to the MOVES [®] SLC [™]	161
16. Selecting the operating mode and associated settings.	116
17. Changing alarm settings.	138

PROCEDURE	Page
18. Connecting to the patient.	161
19. Viewing the Main screen.	126
20. Viewing graphs and trends of patient data.	170
21. Taking NIBP readings manually.	103
22. Using suction.	85, 168
23. Shutting down MOVES [®] SLC [™] .	168

THIS PAGE DELIBERATELY LEFT BLANK.

3.0 Introduction

This Operator's Manual is a reference guide for the MOVES[®] SLC[™] mobile life support system (P/N 122752). This manual includes illustrations, annotated photographs and procedures designed to assist in the operation of various systems, subsystems and components that comprise the unit.

Any operator of the MOVES[®] SLC[™] system must read this manual and accompanying accessory manuals/instructions in their entirety prior to use to safely and effectively operate the system. The operator, or a designated healthcare professional, must disclose the risks and associated mitigation steps to the subject on which the system will be used.

Keep this manual in a dry, convenient location for easy access. All information, illustrations, photographed procedures and specifications in this manual represent the most current product information available at the time of publication.

3.1 General Information

3.1.1 Manufacturer's Information

For warranty, parts, repair or customer service, please contact Thornhill Research Inc. See *Contact Information* on page *iii* for full contact details.

3.1.2 Glossary of Terms and Abbreviations

Table 2: Glossary of Terms and Abbreviations

TERM / ABBREVIATION	DEFINITION
AC	Alternating Current. A type of electrical current in which the direction of the electrical flow switches back and forth regularly.
A/C	Assist / Control Ventilation. In A/C mode, the system delivers a specific tidal volume or PIP at specific intervals based on a patient inspiratory trigger or time trigger.
ABP	Arterial Blood Pressure. ABP is measured in millimeters of mercury (mmHg).
BPM	Beats per Minute
B/M	Breaths per Minute
CFM	Cubic Feet per Minute
CPAP	Continuous Positive Airway Pressure
CVP	Central Venous Pressure. CVP is measured in millimeters of mercury (mmHg).
DC	Direct Current. A type of electrical current in which the electrons always flow in the same direction.
ECG	Electrocardiogram
EtCO ₂	End tidal Carbon Dioxide (mmHg)
FiO ₂	Fraction of Inspired Oxygen by Volume (%)
FRC	Functional Residual Capacity – the volume of air (about 3 liters in an adult) that is present in the lungs at the end of a normal expiration.
Fresh Gas	Gas which has a negligible concentration of CO ₂ .
HC Filter	Hydrocarbon Filter
HR	Heart Rate
I/E ratio	Ratio of Inspiratory Time / Expiratory Time

TERM / ABBREVIATION	DEFINITION
IP	Invasive Pressure. Can refer collectively or individually to ABP, CVP or ICP (which is not a "blood" pressure).
ICP	Intracranial Pressure. ICP is measured in centimeters of water (cmH ₂ O).
IMV	Intermittent Mandatory Ventilation. In IMV mode, the system delivers a specific tidal volume or PIP at specific intervals based on specific time triggers.
LED	Light Emitting Diode
LPM	Liters Per Minute
NATO	North Atlantic Treaty Organization
NIBP	Non-Invasive Blood Pressure. Measured via a blood pressure cuff. NIBP is measured in millimeters of mercury (mmHg).
OEM	Original Equipment Manufacturer. A company that produces hardware to be sold under another company's brand.
Oxygen Supplementation	The provision of gas containing a higher concentration of oxygen than ambient air.
Paw	Airway Pressure
PC	Pressure Control Ventilation
PCB	Printed Circuit Board
pCO ₂	The partial pressure of CO ₂ measured in mmHg.
PEEP	Positive End-Expiratory Pressure in centimeters of water pressure (cmH ₂ O).
PI	Perfusion Index (Perfusion Index, or PI, is a relative assessment of the pulse strength at the monitoring site.)
PIP	Peak Inspiratory Pressure measured in centimeters of water (cmH ₂ O).
PS	Pressure Support. PS assists a patient's inspiratory effort through the application of an additional set level of pressure above PEEP.
PSV	Pressure Support Ventilation. In PSV mode, the system delivers pressure at a specific PS pressure over PEEP during inhalation and a PEEP pressure during exhalation. The triggering and cycling is entirely controlled by the patient.
RR	Respiratory Rate. Breaths per minute (B/M).
SpCO	Saturation percentage of carbon monoxide attached to hemoglobin. CO (carbon monoxide) competes with oxygen for the oxygen-binding sites on hemoglobin. The binding of CO to hemoglobin results in the formation of a compound called Carboxyhemoglobin (COHb). This compound is unable to transport or transfer oxygen.
SpHb	Saturation percentage of hemoglobin. Hemoglobin is the part of a red blood cell that carries oxygen to the body. SpHb measures total hemoglobin and indicates the oxygen carrying capacity of the blood.
SpMet	Saturation percentage of methemoglobin. (Methemoglobin [MetHb] is an oxidized form of hemoglobin that is unable to carry oxygen.)
SpOC	Total oxygen content. (SpHb and SpO ₂ are used together to calculate the actual amount of oxygen in the blood.)
SpO ₂	Oxygen saturation of hemoglobin. Arterial oxygen saturation of hemoglobin as read from a pulse oximeter. It is measured as a percentage (%) of oxy-hemoglobin present in arterial blood in relation to total hemoglobin.

TERM / ABBREVIATION	DEFINITION
SIMV	Synchronized Intermittent Mandatory Ventilation. In SIMV mode, the system delivers breaths synchronized with the patient and ensures that a minimum number of breaths of a specified tidal volume or PIP are delivered. Additionally, any breaths beyond the minimum set number can be supported with a specific level of pressure.
VC	Volume Control Ventilation
Vt	Breath (tidal) volume. Typically measured in milliliters (mL) or liters (L). MOVES [®] SLC [™] measures only in milliliters (mL).

3.2 MOVES[®] SLC[™] Intended Use

MOVES[®] SLC[™] is a portable, computer-controlled, electrically powered emergency-transport ventilator intended to provide continuous or intermittent ventilatory support for the care of adults who require mechanical ventilation. MOVES[®] SLC[™] can also deliver elevated inspired oxygen to ventilated patients.

MOVES[®] SLC[™] is intended to deliver high inspired oxygen concentrations to spontaneously breathing patients who require elevated inspired oxygen.

MOVES[®] SLC[™] is intended to be used in a transport or emergency setting on adult patients who weigh between 40 kg and 120 kg.

MOVES[®] SLC[™] provides the following supplemental functions for the patients it is ventilating or supplying with supplemental oxygen:

a. Suction

The MOVES[®] SLC[™] suction pump is intended for aspiration and removal of fluids, tissue (including bone), gases, bodily fluids or infectious materials from wounds or from a patient's airway or respiratory support system.

b. Patient Monitoring

MOVES[®] SLC[™] is intended to monitor physiological parameters of patients and provide these parameters to a health care provider for interpretation in the form of physiological data and system alarms. Physiological data and system alarms will be available to the care provider from the monitor.



WARNING! MOVES[®] SLC[™] IS INTENDED FOR USE ON ONE PATIENT AT A TIME. IT SHOULD NOT, FOR EXAMPLE, BE USED TO VENTILATE ONE PATIENT WHILE MONITORING ANOTHER.

3.2.1 Operating Environment

MOVES[®] SLC[™] is intended to be operated in a transport or emergency setting. See *Section 16.6.7 Environmental Specifications* on page 214 for a list of environmental operating conditions.

3.2.2 Intended Operator

MOVES[®] SLC[™] is intended to be used by, or under the supervision of, medically qualified and trained personnel.

3.2.3 Intended Population

The intended patient population is adult patients who weigh between 40 kg and 120 kg.


THIS PAGE DELIBERATELY LEFT BLANK

4.0 Regulatory Compliance

4.1 Regulatory Symbols

Regulatory symbols have been added to the labeling on the MOVES[®] SLC[™] unit, power supply and battery charger, and accompanying accessories to indicate regulatory compliance. These symbols, along with a brief description, are shown in the following table.

Table 3: Regulatory Symbols Used and Description

SYMBOL	DESCRIPTION
	CSA (Canadian Standards Association) Mark – This symbol appears on the MOVES [®] SLC [™] Power Supply / Battery Charger. It indicates that approval by the Canadian Standards Association has been granted for use in Canada (C) and the United States (US) under the Certificate Number TBD.
	"Conformité Européen" Mark – The CE mark is a self-declaratory mark that indicates the manufacturer or the importer of record has ensured that all of the applicable European safety and conformity directives and standards have been applied to the product.
The presence of the following symbols indicates compliance with International Electrotechnical Commission (IEC) standard IEC 60601-1 for Medical Electrical Equipment.	
	Consult Accompanying Documents – This symbol appears on the Product Information labels for the MOVES [®] SLC [™] unit and directs the operator to consult the accompanying documents. It also appears frequently in this operator's manual as a general warning and caution symbol.
	Type BF Applied Part – Indicates the device provides an intermediate degree of protection should the patient come in contact with an unintended source of voltage from an external source, but it is not approved for direct cardiac application.
	Type CF Applied Part Defibrillator-Proof - Indicates the applied part provides a high degree of protection should the patient come in contact with an unintended source of voltage from an external source, and that it is approved for direct cardiac application.
	ESD (Electrostatic Discharge) Susceptibility Symbol – The ESD Susceptibility Symbol indicates that handling or use of the item to which the symbol is applied may result in damage from ESD if proper grounding precautions are not taken.

4.2 Regulatory Standards Compliance

Table 4: Regulatory Standards Compliance

STANDARD #	DESCRIPTION
IEC 60601-1	Medical Electrical Equipment (Ed 3.1, 2012), General Requirements for basic safety and essential performance.
IEC 60601-1-2	Medical Electrical Equipment – Part 1-2: General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests (Ed 3.0 2007)
IEC 60601-1-8	Medical Electrical Equipment – Part 1-8: General requirements for safety – Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems (Ed 2.1 2012)
ISO 80601-2-13	Medical Electrical Equipment – Part 2: Particular Requirements for the safety and essential performance of anaesthetic systems (Ed 1.0 2011)
IEC 60601-2-27	Medical Electrical Equipment - Part 2: Particular Requirements for the Safety, including Essential Performance, of Electrocardiographic Monitoring Equipment (Ed. 3.0, 2011)
IEC 80601-2-30	Particular requirements for the safety, including essential performance, of automatic cycling non-invasive blood pressure monitoring equipment (2009)
IEC 60601-2-34	Medical Electrical Equipment - Part 2: Particular Requirements for The Safety, Including Essential Performance, of Invasive Blood Pressure Monitoring Equipment (Ed. 3.0, 2011)
IEC 60601-2-49	Medical Electrical Equipment - Part 2-49: Particular Requirements For The Safety Of Multifunction Patient Monitoring Equipment (Ed. 2.0, 2011)
ISO 80601-2-61	Medical Electrical Equipment - Particular requirements for the basic safety and essential performance of pulse oximeter equipment for medical use (2011)
ISO 80601-2-55	Medical Electrical Equipment - Particular requirements for the basic safety and essential performance of respiratory gas monitors (2011)
ANSI/AAMI EC13	Cardiac monitors, heart rate meters and alarms (2002 / (R)2007)
ANSI/AAMI SP10	Manual, Electronic or Automated Sphygmomanometers (2002 + A1:2003/®2008)
ASTM E1112-00	Standard Specification for Electronic Thermometer for Intermittent Determination of Patient Temperature (2006)
BS EN 794-3	Lung ventilators. Particular requirements for emergency and transport ventilators (1998+A2:2009)
MIL-STD-810G	Department of Defense Test Method Standard for Environmental Engineering Considerations and Laboratory Tests
JECETS	Joint Enroute Care Equipment Test Standard

5.0 MOVES[®] SLC[™] System Overview

5.1 General Overview

The MOVES[®] SLC[™] system is comprised of six main modules:

- Oxygen Concentrator
- Ventilator
- Suction System
- Patient Monitoring System
- Disposable breathing cartridges and breathing circuits for intubated and non-intubated patients
- Power system, comprised of removable hot swappable batteries and an AC power supply and battery charger



NOTE: Color shown in pictures may differ from the actual system or accessories



NOTE: Used breathing cartridges and breathing circuits should be disposed of in accordance with local biohazard regulations.

MOVES[®] SLC[™] is designed to be attached to NATO standard compliant *folding* stretchers such as the Talon II[™] and Raven[™] manufactured by North American Rescue Products Inc.

5.2 System Orientation and External Components

Throughout this manual there are references to aid in orienting the caregiver with the positioning of the unit, especially when the orientation must change in order to conduct a given operating procedure. References to key external components may be given to further assist in the operator's orientation.

The following orientations should be remembered:

- **FRONT:** End of MOVES[®] SLC[™] that contains the patient connection panel
- **REAR:** End of MOVES[®] SLC[™] that contains the battery rack
- **RIGHT and LEFT Sides:** As seen looking from FRONT to REAR.
- **MOVES[®] SLC[™]** is to be operated in the upright position.




Figure 5-1: MOVES[®] SLC[™] Upright System Orientation and Components

5.3 Theory of Operation

5.3.1 Breathing Circuit and Oxygen Supplement

Intubated patients are ventilated with a circle circuit (called a Ventilator Circuit) used with a Ventilator Cartridge. Oxygen and air enter the circuit from two possible sources: the oxygen concentrator or the air pump. The oxygen concentrator can be set to provide an intubated patient with an inspired O₂ concentration (FiO₂) between 30% and 85%. In Ventilator mode, the air pump serves as a back-up source of air flow if the system concentration of oxygen is too low or carbon dioxide is too high or either is unknown.

The ventilator cartridge is designed to remove CO₂ from the circle circuit for an approximately 2-hour period when the system operates at room temperature. The system provides a warning when the CO₂ level in inspired gas exceeds 6 mmHg, which is indicative of cartridge exhaustion, so the cartridge can be changed. If the cartridge is not changed, the system will enter a safe mode and increase oxygen production and air supply.

5.3.2 Oxygen Concentrator

The concentrator provides up to 87% O₂ into the ventilator breathing circuit or directly during O₂ Supplement mode.

5.3.3 Ventilator

The ventilator is comprised of a blower, a sealed ventilator chamber that houses the ventilator bag (air / oxygen reservoir), and a valve block that interfaces with the breathing cartridges and ventilator breathing circuit. The valve block includes the following components:

- Inspiratory and expiratory flow sensors
- One-way valves which are used to direct gas to and from the patient during inhalation and exhalation
- A high pressure relief valve ensures that the pressure in the breathing circuit will never exceed 100 cmH₂O at 60 LPM of flow. (NOTE: PEEP and PIP values are displayed from 0–70 cmH₂O. The pressure relief valve starts to trigger when the pressure exceeds 70 cmH₂O.)
- A port for measuring airway pressure
- A latch that secures the breathing cartridge to the valve block
- A ring mount for the ventilator bag

During patient inhalation, the blower pressurizes the ventilator chamber which displaces the gas from the ventilator bag into the patient's breathing circuit. The blower draws external air through an inlet filter. During exhalation, the patient's expired gas is directed to the ventilator bag. The ventilator bag also receives oxygen from the concentrator or external O₂ source and air from the air pump.

5.3.4 Respiratory Gas Monitoring



NOTE: The system reports O₂/CO₂ gas measurements as ATPD (Ambient temperature [variable] and pressure [variable], dry [no humidity]). O₂/CO₂ readings are corrected based on water vapor pressure measurements.

The Ventilator breathing circuit is equipped with a filter on the Y-piece plus a sampling line and filter which is attached to the GAS SAMPLE port of the patient-connection panel on the MOVES[®] SLC[™]. To prevent equipment failure, only the sampling line and filter supplied with MOVES[®] SLC[™] should be used. The filter should be attached to the GAS SAMPLE port at all times unless the GAS SAMPLE port is capped. In addition, only the filter supplied with the Y-piece should be used with the Y-piece.

The oxygen sensor detects the concentration of oxygen – between 5 and 100% – that is being supplied to and exhaled from the patient and reports it in a plot as a percentage. The CO₂ sensor detects the concentration of CO₂ – between 0 and 10% – that is being supplied to and exhaled from the patient and then, taking into account ambient barometric pressure, reports it in a plot as a partial pressure in mmHg.

Under normal operating conditions (0°C to 40°C), the O₂ and CO₂ sensor should reach specified operating performance in less than two (2) minutes. Until the O₂ and CO₂ sensor reaches specified operating performance, the User Interface Display will display dashes and read as blank (“ – –”) for O₂ and CO₂ values, and a message “CO₂/O₂ sensors warming up” will be displayed in the alarm queue of the UI. When values begin to be displayed, the O₂ and CO₂ sensor has reached specified operating performance.



WARNING! DO NOT REUSE SAMPLING LINES OR FILTERS. THIS COULD PRESENT A DANGER OF INFECTION.



NOTE: Compensation for barometric pressure and temperature is performed by MOVES[®] SLC[™] internal sensing equipment to maintain accuracy of gas calibration over the MOVES[®] SLC[™] environmental operating range. For information on the MOVES[®] SLC[™] environmental operating range see the Appendix.

About O₂ and CO₂ Sensors and Calibration

O₂ and CO₂ calibration are performed at device startup after the O₂ and CO₂ sensors have initialized and warmed up. “CAL” is displayed in each O₂ and CO₂ parameter (i.e., FiO₂, PetCO₂ and PiCO₂) during warm up and calibration. The calibration is based on measuring the air surrounding the MOVES[®] SLC[™]. If the initial attempt at O₂ calibration fails, on each subsequent O₂ calibration check, if the check is in bounds, a calibration is attempted until O₂ calibration succeeds.

A CO₂ calibration is performed every 30 minutes. O₂ calibration is only performed once, and then an O₂ calibration check is performed every 30 minutes. However, if SpO₂ is at or below 85%, the O₂ calibration check is performed every five (5) minutes until the condition clears.



NOTE: The time between O₂ calibration checks is never less than five (5) minutes, even if SpO₂ repeatedly falls below 85%.



NOTE: If the O₂ calibration check fails high, the system automatically enters Safe Gas Mode (SGM). See Safe Gas Mode on page 205 for further information.

The O₂ calibration check fails if the surrounding air's average O₂ reading $\pm$ one (1) standard deviation (data point taken every 90 ms over a five [5] second period) is not contained within the bounds of 19.4 - 22.4%. The O₂ calibration check fails *high* if the air's average O₂ reading plus one (1) standard deviation is above the upper bound; otherwise, the O₂ calibration check fails *low* if the air's average O₂ reading minus one (1) standard deviation is below the lower bound.

If the O₂ calibration check fails *high*, then the MOVES[®] SLC[™] enters Safe Gas Mode (SGM) until an O₂ calibration check does not fail *high* (when the MOVES[®] SLC[™] O₂ calibration is next checked). While the MOVES[®] SLC[™] is in SGM due to an O₂ calibration check failure *high*, the O₂ calibration recheck occurs every five (5) minutes. If the O₂ calibration check fails *low*, subsequent O₂ calibration checks are scheduled at ten (10) minute intervals while each subsequent O₂ calibration check continues to fail *low*.

A high priority alarm stating "O₂ reading may be biased high" will appear in the alarm queue when the O₂ calibration check fails *high*. A low priority alarm stating "O₂ reading may be biased low" alarm will appear in the alarm queue when the O₂ calibration check fails *low*.

Whenever a CO₂ calibration is scheduled, an O₂ calibration check is performed at the same time. This minimizes the time that patient gas monitoring is suspended. Whenever an O₂ calibration check is performed, and the CO₂ calibration is scheduled to occur within the next 11 minutes, then a CO₂ calibration is performed at the same time as the O₂ calibration check.

Asterisk Beside Sensor Values

An O₂ calibration check and/or CO₂ calibration requires 15 seconds to complete. O₂ (FiO₂), CO₂ (PetCO₂ and PiCO₂) and RR values are displayed but are "frozen" while the O₂ calibration check and/or CO₂ calibration is performed. During this time, an asterisk is displayed beside each numeric value. If a parameter has no value (i.e., dashes, "---", are displayed), then no asterisk is shown beside that parameter.

Adverse Effects

There are no known adverse effects associated with MOVES[®] SLC[™] respiratory gas monitoring in itself. However, sampling lines and filters are not reusable and could present a danger of infection were they to be reused. They should be disposed of in accordance with local biohazard regulations.

5.3.5 Suction

The system provides variable suction between -100 and -325 mmHg with flow rates of 20 L/min. A suction kit, consisting of a wand, two hoses, and an 800 mL suction filtration canister with a canister holder, is attached to the suction port of the MOVES[®] SLC[™]. The aspirated air is vented through the exhaust of the MOVES[®] SLC[™] system. Overfill protection for the suction container is provided by a filter in the lid which blocks flow when canister capacity is reached. The MOVES[®] SLC[™] system also contains a suction safety relief valve in the suction path which opens at approximately -415 mmHg.



CAUTION! OPERATING SUCTION WITHOUT THE SUPPLIED FILTER, CANISTER, AND HOSES WILL RESULT IN PERMANENT FAILURE OF THE SUCTION SYSTEM AND THE OXYGEN CONCENTRATOR.

5.3.6 Patient Monitoring

Non-Invasive Blood Pressure (NIBP)

Using the appropriate cuff, NIBP measurements can be completed on either an arm or thigh. NIBP measurements made with the MOVES[®] SLC[™] for adult patient populations are equivalent to those obtained by trained observers using the cuff/stethoscope auscultatory method within the limits prescribed by ANSI/AAMI SP10: 1992 & 2002 (mean error difference of ± 5 mmHg or less, standard deviation of 8 mmHg or less). In addition, blood pressure measurements determined with the MOVES[®] SLC[™] are equivalent to those obtained by a trained observer using the cuff / stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or automated sphygmomanometers.

Pulse Oximetry

The Masimo Rainbow SET[®] Pulse Co-Oximeter measures the functional oxygen saturation of arterial hemoglobin (% SpO₂). Significant levels of dysfunctional hemoglobin, such as carboxyhemoglobin or methemoglobin, may affect the accuracy of the measurement. Factors that may degrade pulse performance or affect the accuracy of the measurement include the following: excessive ambient light, excessive motion, electrosurgical interference, blood flow restrictors (arterial catheters, blood pressure cuffs, infusing lines, etc.), moisture in the sensor, improperly applied sensor, incorrect sensor type, poor pulse quality, venous pulsations, anemia or low hemoglobin concentrations, cardiogreen or other intravascular dyes, carboxyhemoglobin, methemoglobin, dysfunctional hemoglobin, artificial nails or fingernail polish, or a sensor not at heart level. The pulse oximeter is calibrated by the original manufacturer to display functional oxygen concentration.



NOTE: For more information, see Section 6.0 The Masimo Rainbow SET[®] Pulse CO-Oximeter beginning on page 25.

Invasive Pressure (IP)

Data is available from three types of Invasive Pressure (IP) sensors: Arterial, Central Venous, and Intracranial. Arterial Blood Pressure (ABP) is shown with numeric systolic and diastolic readings. The display is capable of displaying a pressure range of -10–300 mmHg for systolic and diastolic pressure readings. Intracranial Pressure (ICP) displays pressures from -14 to 408 cmH₂O. Central Venous Pressure (CVP) display pressures from -10 to 300 mmHg. For all invasive pressures, if the reading is below the sensor range, '<-10' is shown for ABP and CVP, and '<-14' is shown for ICP. If the reading is above the sensor range, '>300' is shown for ABP and CVP, and '>408' is shown for ICP.

Temperature

When the temperature probe is attached to the patient and connected to the MOVES[®] SLC[™], the patient's body temperature is displayed on the Monitor Screen. The temperature can be displayed in degrees Fahrenheit from 82.4° F to 107.6°F or degrees Celsius from 28°C to 42°C. Changing the temperature display from Celsius to Fahrenheit (or vice versa) may only be done by a qualified MOVES[®] SLC[™] service technician.

Electrocardiogram (ECG)

MOVES[®] SLC[™] uses a standard 12-lead ECG system for monitoring the heart and produces 12 ECG channels (I, II, III, aVL, aVR, aVF, V1, V2, V3, V4, V5 and V6). ECG data can be displayed on the Graphs / Trends section of the Monitoring Screen with a vertical scale (for data) and a horizontal scale (for time). The HR (Heart Rate) area of the Monitoring Screen displays heart rate in BPM (Beats per Minute) and data can be drawn from any of the ECG channel waveforms shown on the ECG monitor. If the ECG monitor is being used, and the readings are erratic, the accuracy

of the heart rate cannot be guaranteed. The ECG heart rate meter's response to irregular rhythms has not been assessed.

5.3.7 Power System

The MOVES[®] SLC[™] operates on either rechargeable batteries (DC) or AC power. The power range of the MOVES[®] SLC[™] power supply is as follows:

Input: 100 - 240 VAC, 50-60 Hz, 5.5 A max.

Output: 28 VDC, 14.3 A max.

The MOVES[®] SLC[™] houses up to two lithium polymer batteries. MOVES[®] SLC[™] will operate on a set of 2 batteries for approximately 2.5 hours minimum. Under typical clinical use (ventilator and monitors running, concentrator on for 30 seconds / off for 90 seconds, assuming no leaks), MOVES[®] SLC[™] should operate at least 4 hours on a set of 2 batteries. Battery run time is highly dependent on the use of the oxygen concentrator or suction.

Power is connected to the MOVES[®] SLC[™] system from an AC source via the MOVES[®] SLC[™] power supply and battery charger. The power supply and battery charger system provides 28 VDC at up to 14.3 A. Batteries are charged by connecting AC power to the MOVES[®] SLC[™] with the batteries installed in the system. The system's batteries should take no more than 2.5 hours to fully charge when the system is idle. Batteries can be charged while the system is running although charge times may be longer.



NOTE: *The battery charge level may not appear to increase for approximately 3 hours. This is normal for the battery's initial charge and after extended periods without use.*

6.0 The Masimo Rainbow SET[®] Pulse CO-Oximeter

6.1 Overview

The MOVES[®] SLC[™] uses the Masimo Rainbow SET[®] Pulse CO-Oximeter to provide non-invasive monitoring that measures arterial oxygen saturation (SpO₂), pulse rate (PR), and perfusion index (PI), along with optional measurements of hemoglobin (SpHb), total oxygen content (SpOC), carboxyhemoglobin (SpCO) and methemoglobin (SpMet).

The Masimo Rainbow SET[®] Pulse CO-Oximeter and accessories have been validated and are indicated for use with adult patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, and transport.

6.2 Key Features

- Masimo SET[®] is clinically proven to satisfy all sensitivity and specificity requirements for pulse oximeter technology.
- Rainbow technology uses 7+ wavelengths of light to continuously and noninvasively measure carboxyhemoglobin (SpCO), methemoglobin (SpMet), and total hemoglobin (SpHb), as well as providing a more reliable probe-off detection.
- Total oxygen content (SpOC) provides a calculated measurement of the amount of oxygen in arterial blood, which may provide useful information about oxygen both dissolved in plasma and combined with hemoglobin.
- Perfusion Index (PI) with trending capability indicates arterial pulse signal strength and may be used as a diagnostic tool during low perfusion.

6.3 Indications for Use

Masimo Rainbow SET[®] Pulse CO-Oximeter and accessories are indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂), pulse rate (PR), carboxyhemoglobin saturation (SpCO), methemoglobin saturation (SpMet) and total hemoglobin concentration (SpHb).

Masimo Rainbow SET[®] Pulse CO-Oximeter and accessories have been validated and are indicated for use with adult patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities and transport.

In addition, the Masimo Rainbow SET[®] Pulse CO-Oximeter and accessories are indicated to provide continuous noninvasive monitoring data of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (PR) to multi-parameter devices for the display of those devices.



WARNING! THE MASIMO RAINBOW SET[®] PULSE CO-OXIMETER IS NOT INTENDED FOR USE AS AN APNEA MONITOR.

6.4 Pulse Oximeter Technology Overview

6.4.1 Signal Extraction Technology (SET[®])

Masimo Signal Extraction Technology's signal processing differs from that of conventional pulse oximeters. Conventional pulse oximeters assume that arterial blood is the only blood moving (pulsating) in the measurement site. During patient motion, however, the venous blood also moves, causing conventional pulse oximeters to read low values, because they cannot distinguish between the arterial and venous blood movement (sometimes referred to as noise).

Masimo SET[®] pulse oximetry utilizes parallel engines and adaptive digital filtering. Adaptive filters are powerful because they are able to adapt to the varying physiologic signals and/or noise and separate them by looking at the whole signal and breaking it down to its fundamental components. The Masimo SET signal processing algorithm, Discrete Saturation Transform[®] (DST[®]), in parallel with Fast Saturation Transform (FST[®]), reliably identifies the noise, isolates it and, using adaptive filters, cancels it. It then reports the true arterial oxygen saturation for display on the monitor.

6.4.2 General Description for Oxygen Saturation (SpO₂)

Pulse oximetry is governed by the following principles:

1. Oxyhemoglobin (oxygenated blood) and deoxyhemoglobin (non-oxygenated blood) differ in their absorption of red and infrared light (spectrophotometry).
2. The amount of arterial blood in tissue changes with your pulse (photoplethysmography). Therefore, the amount of light absorbed by the varying quantities of arterial blood changes as well.

6.4.3 Successful Monitoring for SpO₂, PR, and PI

Stability of the SpO₂ readings may be a good indicator of signal validity. Although stability is a relative term, experience will provide a good feeling for changes that are artifactual or physiological and the speed, timing, and behavior of each.

The stability of the readings over time is affected by the averaging mode being used. The longer the averaging time, the more stable the readings tend to become. This is due to a dampened response as the signal is averaged over a longer period of time than during shorter averaging times. However, longer averaging times delay the response of the oximeter and reduce the measured variations of SpO₂ and pulse rate.

6.4.4 Functional Oxygen Saturation

The Masimo Rainbow SET[®] Pulse CO-Oximeter is calibrated to measure and display functional oxygen saturation (SpO₂): the amount of oxyhemoglobin expressed as a percentage of the hemoglobin that is available to transport oxygen.

Note that carboxyhemoglobin is not capable of transporting oxygen, but is recognized as oxygenated hemoglobin by conventional pulse oximetry.

6.4.5 General Description for Pulse Rate (PR)

Pulse rate (PR), measured in beats per minute (BPM), is based on the optical detection of peripheral flow pulse.

6.4.6 General Description for Perfusion Index (PI)

The Perfusion Index (PI) is the ratio of the pulsatile blood flow to the non-pulsatile or static blood in peripheral tissue. PI thus represents a noninvasive measure of peripheral perfusion that can be continuously and noninvasively obtained from a pulse oximeter.

6.4.7 Rainbow Pulse CO-Oximetry Technology

Rainbow Pulse CO-Oximetry technology is governed by the following principles:

1. Oxyhemoglobin (oxygenated blood), deoxyhemoglobin (non-oxygenated blood), carboxyhemoglobin (blood with carbon monoxide content), methemoglobin (blood with oxidized hemoglobin) and blood plasma constituents differ in their absorption of visible and infrared light (using spectrophotometry).

2. The amount of arterial blood in tissue changes with pulse (photoplethysmography). Therefore, the amount of light absorbed by the varying quantities of arterial blood changes as well.

The Masimo Rainbow SET[®] Pulse CO-Oximeter uses a multi-wavelength sensor to distinguish between oxygenated blood, deoxygenated blood, blood with carbon monoxide, oxidized blood and blood plasma.

Absorption Spectra




Figure 6-1: Absorption Spectra of Blood Components



NOTE: The wavelength of blood plasma, which begins near the 1000 nm range and peaks in the 1400 nm range, is omitted from the above graph since most of its wavelength falls outside of the graph parameters.

The Masimo Rainbow SET[®] Pulse CO-Oximeter utilizes a sensor with various light-emitting diodes (LEDs) that pass light through the site to a diode (detector). Signal data is obtained by passing various visible and infrared lights (LEDs, 500 to 1400 nm) through a capillary bed (for example, a fingertip, a hand, a foot) and measuring changes in light absorption during the blood pulsatile cycle. This information may be useful to clinicians. The maximum radiant power of the strongest light is rated at ≤ 25 mW. The detector receives the light, converts it into an electronic signal and sends it to the Masimo Rainbow SET[®] Pulse CO-Oximeter for calculation.




Figure 6-2: LEDs and Detector

Once the Masimo Rainbow SET[®] Pulse CO-Oximeter receives the signal from the sensor, it utilizes proprietary algorithms to calculate the patient's functional oxygen saturation (SpO₂ [%]), blood levels of carboxyhemoglobin (SpCO [%]), methemoglobin (SpMet [%]), total hemoglobin concentration (SpHb [g/dL]) and pulse rate (PR). The SpCO, SpMet and SpHb measurements rely on a multi-wavelength calibration equation to quantify the percentage of carbon monoxide and methemoglobin and the concentration of total hemoglobin in arterial blood. In an ambient temperature of 35° C the maximum skin surface temperature has been measured at less than 106° F (41° C), verified by Masimo sensor skin temperature test procedure.

6.4.8 Pulse CO-Oximetry vs. Drawn Whole Blood Measurements

When SpO₂, SpCO, SpMet, and SpHb measurements obtained from the Masimo Rainbow SET[®] Pulse CO-Oximeter (noninvasive) are compared to drawn whole blood (invasive) measurements by blood gas and/or laboratory CO-Oximetry methods, caution should be taken when evaluating and interpreting the results.

The blood gas and/or laboratory CO-Oximetry measurements may differ from the SpO₂, SpCO, SpMet, SpHb, and SpOC measurements of the Masimo Rainbow SET[®] Pulse CO-Oximeter. Any comparisons should be simultaneous, meaning the measurement on the device should be noted at the exact time that blood is drawn.

In the case of SpO₂, different results are usually obtained from the arterial blood gas sample if the calculated measurement is not appropriately corrected for the effects of variables that shift the relationship between the partial pressure of oxygen (PO₂) and saturation, such as: pH, temperature, the partial pressure of carbon dioxide (PCO₂), 2,3-DPG, and fetal hemoglobin. In the case of SpCO, different results are also expected if concentration of methemoglobin in the blood gas sample is abnormal (greater than 2% for methemoglobin concentration).

High levels of bilirubin may cause erroneous SpO₂, SpMet, SpCO, and SpHb readings. As blood samples are usually taken over a period of 20 seconds (the time it takes to draw the blood) a meaningful comparison can only be achieved if the oxygen saturation, carboxyhemoglobin, and methemoglobin concentration of the patient are stable and not changing over the period of time that the blood gas sample is taken. Subsequently, blood gas and laboratory CO-Oximetry measurements of SpO₂, SpCO, SpMet, SpHb, and SpOC may vary with the rapid administration of fluids and in procedures such as dialysis. Additionally, drawn whole blood testing can be affected by sample handling methods and time elapsed between blood draw and sample testing.

Measurements with Low Signal IQ should not be compared to laboratory measurements.

6.4.9 General Description for Total Hemoglobin (SpHb)

Pulse CO-Oximetry is a continuous and noninvasive method of measuring the levels of total hemoglobin (SpHb) in arterial blood. It relies on the same principles of pulse oximetry to make its SpHb measurement. The measurement is taken by a sensor capable of measuring SpHb, usually on the fingertip for adult patients.

The sensor connects directly to the Pulse CO-Oximeter or with a patient cable. The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data as measurement of total hemoglobin concentration.

6.4.10 Successful Monitoring for SpHb

A stable SpHb reading is associated with correct sensor placement, small physiological changes during the measurement and acceptable levels of arterial perfusion at the measurement site. Physiological changes at the measurement site are mainly caused by fluctuations in the oxygen saturation, blood concentration and perfusion.

6.4.11 General Description for SpOC

The following is the equation for oxygen content via the Pulse CO-Oximeter:

$$SpOC \text{ (ml/dL*)} = 1.31 \text{ (ml O}_2\text{/g Hb)} \times SpHb \text{ (g/dL)} \times SpO_2 + 0.3 \text{ ml/dL}$$

*When ml O₂/g Hb is multiplied by g/dL of SpHb, the gram unit in the denominator of ml/g cancels the gram unit in the numerator of g/dL resulting in ml/dL (ml of oxygen in one dL of blood) as the unit of measure for SpOC.

6.4.12 General Description for Carboxyhemoglobin (SpCO)

Pulse CO-Oximetry is a continuous and noninvasive method of measuring the levels of carboxyhemoglobin concentration (SpCO) in arterial blood. It relies on the same basic principles of pulse oximetry (spectrophotometry) to make its SpCO measurement.

The measurement is obtained by placing a sensor on a patient, usually on the fingertip for adult patients. The sensor connects either directly to the Pulse CO-Oximetry instrument or through an instrument patient cable.

The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data as percentage value for the SpCO, which reflect blood levels of carbon monoxide bound to hemoglobin.

6.4.13 Successful Monitoring for SpCO

A stable SpCO reading is associated with correct sensor placement, small physiological changes during the measurement and acceptable levels of arterial perfusion in the patient's fingertip (measurement site). Physiological changes at the measurement site are mainly caused by fluctuations in the oxygen saturation, blood concentration and perfusion.

6.4.14 General Description for Methemoglobin (SpMet)

Pulse CO-Oximetry is a continuous and noninvasive method of measuring the levels of methemoglobin concentration (SpMet) in arterial blood. It relies on the same basic principles of pulse oximetry (spectrophotometry) to make its SpMet measurement.

The measurement is obtained by placing a sensor on a patient, usually on the fingertip for adults patients. The sensor connects either directly to the Pulse CO-Oximetry instrument or through a patient cable. The sensor collects signal data from the patient and sends it to the instrument. The instrument displays the calculated data as percentage value for the SpMet.

6.4.15 Successful Monitoring for SpMet

A stable SpMet reading is associated with correct sensor placement, small physiological changes during the measurement and acceptable levels of arterial perfusion in the patient's fingertip (measurement site).

Physiological changes at the measurement site are mainly caused by fluctuations in the oxygen saturation, blood concentration and perfusion.

6.4.16 SpCO, SpMet, and SpHb Measurements During Patient Motion

The MOVES[®] SLC[™] displays measurements of SpCO, SpMet, and SpHb during patient motion. However, because of the changes in the physiological parameters such as blood volume, arterial-venous coupling, etc., that occur during patient motion, the accuracy of such measurements may not be reliable during excessive motion. In this case, the measurement value for SpCO, SpMet, or SpHb displays as dashes (---) and a message displays to alert the clinician that the instrument does not have confidence in the value due to poor signal quality caused by excessive motion or other signal interference.

THIS PAGE DELIBERATELY LEFT BLANK.

7.0 Battery Storage



CAUTION! WHEN BATTERIES ARE DISCHARGED AND LEFT IN MOVES[®] SLC[™] FOR A PROLONGED PERIOD, A COMPLETELY DISCHARGED BATTERY CAN RESULT. THE BATTERY CAN STILL BE RECHARGED, BUT IT MAY TAKE MORE THAN THE NORMAL 2.5 HOURS -- AND UP TO 48 HOURS -- TO FULLY CHARGE.

When MOVES[®] SLC[™] is not in operation, its batteries should be disconnected or removed from the unit and stored in a dry area at room temperature. In the picture below, the battery compartment latch is shown as open. When the latch is in this position, the battery is disconnected.



Figure 7-1: Battery Compartment Latch Open

Conversely, make sure the battery compartment latch is fully closed (by pressing it in) to engage the battery. Even when the battery compartment door is closed, the battery will not be connected unless the latch is pressed in fully.



Figure 7-2: Battery Compartment Latch Fully Closed

Also, please note the following important cautions (and final Warning) about storing MOVES[®] SLC[™] batteries:

- Do not store batteries above 60°C (140°F) or below –20°C (–4°F).
- Store batteries in a cool (below 30°C [86°F]), dry area that is subject to little temperature change. Elevated temperatures can result in reduced battery service life.
- Do not place batteries near heating equipment, nor expose to direct sunlight for long periods.



WARNING! BATTERY EXPOSURE TO TEMPERATURES IN EXCESS OF 130°C (266°F) WILL RESULT IN THE BATTERY VENTING FLAMMABLE LIQUID AND GASES.

8.0 Safety Information

8.1 Manual Symbols and Messages

8.1.1 Notes, Cautions and Warnings

This manual contains important messages with symbols labeled **NOTE**, **CAUTION** and **WARNING**. These messages have the following format and meaning:

Table 5: Symbols and Messages Used in Manual

 NOTE:	<i>Supplies additional information that will help complete, offer an alternative to, or explain a portion of a given procedure.</i>
 CAUTION!	PROVIDES INFORMATION TO PREVENT ERRORS OR INDISCRETIONS THAT COULD RESULT IN EQUIPMENT, SYSTEM, OR COMPONENT DAMAGE.
 WARNING!	INDICATES AREAS WHERE INSUFFICIENT KNOWLEDGE OF A GIVEN PROCEDURE, IMPROPER HANDLING, OR LACK OF ATTENTION COULD RESULT IN PERSONAL INJURY OR LOSS OF LIFE!

Read each labeled message carefully, and follow its instructions during operation to reduce the risk of system or component damage and/or personal injury.




WARNING! IMPROPER OPERATION OF THE MOVES[®] SLC[™] COULD ENDANGER A PATIENT!

Since it is virtually impossible to foresee all of the possible consequences resulting from the failure to follow instructions and adhere to safety procedures, the NOTES, CAUTIONS and WARNINGS contained in this manual are not exhaustive. It is the responsibility of the individual operating MOVES[®] SLC[™] to make safety the number one priority during operating procedures.

8.1.2 Operational Symbols

Table 6: Operational Symbols and Descriptions

SYMBOL	DESCRIPTION
	Single Phase Alternating Current
	Direct Current

8.1.3 Label Warning Symbols


Table 7: Label Warning Symbols and Descriptions

SYMBOL	DESCRIPTION
	Fire Hazard: Do not smoke near unit.
	Fire Hazard: Do not operate device near open flame.

8.1.4 Product Label Symbols

Table 8: Product Label Symbols and Descriptions

SYMBOL	DESCRIPTION
	Caution. Read accompanying documentation
	Consult instructions for use
	Fragile. Sensitive to mechanical shock.


SYMBOL	DESCRIPTION
	Do not reuse (single use only)
	No latex used in the manufacture of this product
	Power indicator
	Battery indicator
	Type BF equipment
	Class II equipment
	CE Mark with notified body number
	Authorized Representative in the European Community
	Manufacturer
	Date of manufacture
	Date of expiration (use by)
	Lot number
	Reference or Model Number

SYMBOL	DESCRIPTION
	Serial Number
	Phone
	Fax
	The equipment shall not be disposed of as unsorted municipal waste and shall be collected as electrical and electronic equipment, as applicable, separately as specified by Waste Electrical and Electronic Equipment (WEEE).
	Battery condition indicator
	Non sterile; material can not be guaranteed to be free of contamination.
	Temperature limitation range for usage. Both upper and lower limits are indicated adjacent to horizontal lines.
	Humidity limitation range for usage. Both upper and lower limits are indicated adjacent to horizontal lines.
	Pressure limitation range for usage. Both upper and lower limits are indicated adjacent to horizontal lines.
	Gas sample line port

8.2 General Warnings

Table 9: General Warnings

SYMBOL	GENERAL WARNING
	WARNING! THE POWER SUPPLY / CHARGER CORD IS SUPPLIED WITH A GROUNDING PRONG ON THE MALE CONNECTOR. TO REDUCE THE RISK OF ELECTRICAL SHOCK, THIS PRONG SHOULD NEVER BE REMOVED OR COMPROMISED.
	WARNING! TO AVOID BREATH STACKING WHEN RUNNING VOLUME CONTROLLED IMV, THE RELEASE PRESSURE MUST BE APPROPRIATELY SET.
	WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATE MEANS OF SUPPLYING A HIGH CONCENTRATION OF O ₂ IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE, OR SERIOUS OCCLUSION IN THE CONCENTRATOR SYSTEM.
	WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATIVE MEANS OF VENTILATION CAPABLE OF SUPPLEMENTING A HIGH CONCENTRATION OF O ₂ IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE, OR SERIOUS OCCLUSION IN THE VENTILATOR SYSTEM.
	WARNING! THE OPERATOR SHOULD USE AN ALTERNATIVE MEANS OF VENTILATION UPON EXPERIENCING A PROLONGED APNEA ALARM.
	WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATE MEANS OF SUCTION IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE OR SERIOUS OCCLUSION IN THE SUCTION CIRCUIT.
	WARNING! IT MAY BE NECESSARY TO USE AN ALTERNATIVE MEANS OF OXYGEN SUPPLEMENTATION SHOULD THE HYDROCARBON FILTER REQUIRE REPLACEMENT WHILE TREATING A PATIENT.
	WARNING! THE MOVES [®] SLC [™] SHOULD NOT BE RUN CONTINUOUSLY IN SAFE GAS MODE. SAFE GAS MODE IS INTENDED FOR <u>SHORT TERM USE ONLY</u> TO COMPLETE TRANSPORTS.
	WARNING! IF LEFT IN A HOT ENVIRONMENT OR DIRECT SUNLIGHT FOR A CONSIDERABLE LENGTH OF TIME, THE MOVES [®] SLC [™] ACCESSORIES CASE AND THE ACCESSORIES IN IT CAN BECOME QUITE HOT. MONITORING OF CASE TEMPERATURE IS RECOMMENDED. ALWAYS MAKE SURE THAT ACCESSORIES THAT WILL BE APPLIED DIRECTLY TO THE PATIENT ARE SUITABLE FOR SKIN CONTACT AND WILL NOT CAUSE BURNS.
	WARNING! BATTERY EXPOSURE TO TEMPERATURES IN EXCESS OF 130°C (266°F) WILL RESULT IN THE BATTERY VENTING FLAMMABLE LIQUID AND GASES.

SYMBOL	GENERAL WARNING
	WARNING! THE MOVES[®] SLC[™] OXYGEN CONCENTRATOR DOES NOT FUNCTION WHILE SUCTION IS ON. AN ALTERNATIVE MEANS OF SUPPLYING O₂ WILL BE NECESSARY IF A HIGH PERCENTAGE OF O₂ IS CRITICAL WHILE SUCTIONING.
	WARNING! ONE CHARGED BATTERY MUST BE PRESENT IN THE MOVES[®] SLC[™] UNIT WHENEVER IT IS IN OPERATION, EVEN WHEN IT IS RUNNING ON EXTERNAL POWER. THIS REDUCES THE RISK TO THE PATIENT IN THE EVENT OF A POWER FAILURE.
	WARNING! CHARGING AND DISCHARGING THE BATTERIES WITH OTHER THAN MOVES[®] SLC[™] SYSTEM EQUIPMENT, AND/OR IMPROPER HANDLING, CAN RESULT IN FIRE, EXPLOSION, TOXIC GASES AND SMOKE.
	WARNING! BATTERY TIME SHOWN REMAINING IS APPROXIMATE AND HIGHLY DEPENDENT ON OPERATING CONDITIONS! PUT SAFETY FIRST – ALWAYS CARRY SPARE BATTERIES!
	WARNING! DO NOT OPERATE THE MOVES[®] SLC[™] SYSTEM WITH A DEFECTIVE BATTERY.
	WARNING! DO NOT USE OR CHARGE A DAMAGED BATTERY!
	WARNING! DO NOT OPERATE THE MOVES[®] SLC[™] SYSTEM UNTIL ALL SYSTEM TEST FAILURES HAVE BEEN RESOLVED, AND ALL TESTS HAVE BEEN REPEATED AND PASSED.
	WARNING! ONLY AUTHORIZED SERVICE AND MAINTENANCE PERSONNEL SHOULD REMOVE ANY COVERS FROM MOVES[®] SLC[™]. UNAUTHORIZED REMOVAL OF COVERS FROM MOVES[®] SLC[™] MAY RESULT IN ELECTRIC SHOCK AND POSSIBLY DEATH, AND MAY DAMAGE THE SYSTEM COMPONENTS.
	WARNING! BECAUSE THE MOVES[®] SLC[™] CONTAINS AN OXYGEN CONCENTRATOR, IT SHOULD ONLY BE USED IN A WELL-VENTILATED ENVIRONMENT AWAY FROM POLLUTANTS, FLAMES, SPARKS, OR FUMES.
	WARNING! LEAKS IN THE SAMPLING LINE CAN CAUSE LOW PCO₂ AND/OR LOW O₂ LEVELS.
	WARNING! WHEN MOVES[®] SLC[™] IS NOT IN OPERATION, BATTERIES SHOULD BE UNLATCHED OR REMOVED FROM THE UNIT AND STORED IN A DRY AREA AT ROOM TEMPERATURE. LEAVING BATTERIES INSTALLED IN A NON-OPERATIONAL UNIT MAY CAUSE THEM TO DRAIN TO AN UNRECHARGEABLE LEVEL.
	WARNING! THE CLAMPS HAVE NUMEROUS MOVING PARTS THAT MAY PRESENT A PINCHING OR CRUSHING HAZARD. ALWAYS USE CAUTION WHEN HANDLING BOTH THE FRONT AND BACK CLAMPS.

8.3 Electrical Warnings

Table 10: Electrical Warnings

SYMBOL	ELECTRICAL WARNING
	WARNING! THE MOVES® SLC™ SYSTEM IS NOT ELECTROSURGERY COMPATIBLE.
	WARNING! THE USE OF MULTIPLE (NON-) MEDICAL ELECTRICAL EQUIPMENT CONNECTED TO THE SAME PATIENT MAY POSE A SAFETY HAZARD DUE TO THE SUMMATION OF LEAKAGE CURRENTS FROM EACH INSTRUMENT. ANY COMBINATION OF (NON-) MEDICAL ELECTRICAL EQUIPMENT SHOULD BE EVALUATED BY LOCAL SAFETY PERSONNEL BEFORE BEING PUT INTO SERVICE.
	WARNING! CONDUCTIVE PARTS OF ELECTRODES AND ASSOCIATED CONNECTORS FOR APPLIED PART, INCLUDING THE NEUTRAL ELECTRODE, SHOULD NOT CONTACT OTHER CONDUCTIVE PARTS AND EARTH.

8.4 Patient-Specific Warnings

Table 11: Patient-Specific Warnings

SYMBOL	PATIENT-SPECIFIC WARNINGS
	WARNING! WHEN USING O₂ SUPPLEMENTATION, AN O₂ SAT MONITOR MUST BE USED.
	WARNING! WHEN USING OXYGEN SUPPLEMENTATION MODE, THE MOVES® SLC™ GAS SAMPLING PORT MUST BE CONNECTED TO THE OXYGEN DELIVERY CIRCUIT (E.G., O₂ MASK SAMPLE PORT) AND AN OXYGEN SHUT-OFF DEVICE, SUCH AS BPR'S FIRESAFE™ CANNULA VALVE, SHOULD BE USED IN THE OXYGEN SUPPLY LINE.
	WARNING! DO NOT CONNECT A PATIENT TO MOVES® SLC™ UNTIL THE MOVES® SLC™ SYSTEM IS PROPERLY WARMED UP AND O₂ VALUES ARE DISPLAYED.
	WARNING! SPO₂ SENSOR APPLICATION SITES SHOULD BE INSPECTED AT LEAST EVERY 6 TO 8 HOURS TO ENSURE CORRECT SENSOR ALIGNMENT AND SKIN INTEGRITY. PATIENT SENSITIVITY MAY VARY DUE TO MEDICAL STATUS OR SKIN CONDITION. DISCONTINUE THE USE OF ADHESIVE TAPE STRIPS IF THE PATIENT EXHIBITS AN ALLERGIC REACTION TO THE ADHESIVE MATERIAL.
	WARNING! ACCURACY OF ANY BLOOD PRESSURE MEASUREMENT MAY BE AFFECTED BY THE POSITION OF THE SUBJECT, HIS OR HER PHYSICAL CONDITION AND USE OUTSIDE OF THE OPERATING INSTRUCTIONS DETAILED IN THIS MANUAL. INTERPRETATION OF BLOOD PRESSURE MEASUREMENTS SHOULD BE MADE ONLY BY A PHYSICIAN OR TRAINED MEDICAL STAFF.

SYMBOL	PATIENT-SPECIFIC WARNINGS
	WARNING! DO NOT ATTACH THE BLOOD PRESSURE CUFF TO A LIMB BEING USED FOR IV INFUSIONS AS THE CUFF INFLATION CAN BLOCK THE INFUSION, POTENTIALLY CAUSING HARM TO THE PATIENT.
	WARNING! IF THE BLOOD PRESSURE CUFF IS ON THE SAME LIMB AS A PULSE OXIMETER PROBE, THE OXYGEN SATURATION RESULTS WILL BE ALTERED WHEN THE CUFF OCCLUDES THE BRACHIAL ARTERY.
	WARNING! TO OBTAIN ACCURATE BLOOD PRESSURE READINGS, THE BLOOD PRESSURE CUFF MUST BE THE CORRECT SIZE, AND ALSO BE CORRECTLY FITTED TO THE PATIENT. INCORRECT SIZE OR INCORRECT FITTING MAY RESULT IN INCORRECT READINGS.
	WARNING! WHEN A BLOOD PRESSURE CUFF IS TO BE POSITIONED ON A PATIENT FOR AN EXTENDED LENGTH OF TIME, BE SURE TO OCCASIONALLY CHECK THE LIMB FOR PROPER CIRCULATION.
	WARNING! BEFORE VENTILATING A PATIENT, ENSURE THAT A SPARE VENTILATOR BREATHING CIRCUIT IS READILY AVAILABLE.
	WARNING! CARE SHOULD BE TAKEN WHEN MONITORING PATIENTS WITH PACEMAKERS SINCE HEART RATE METERS MAY FALSELY COUNT PACEMAKER PULSES.
	WARNING! IRREGULAR HEART RHYTHMS SUCH AS PREMATURE ATRIAL OR VENTRICULAR BEATS MAY CAUSE THE HEART RATE TO BE UNDERESTIMATED.
	WARNING! DO NOT REUSE SAMPLING LINES OR FILTERS. THIS COULD PRESENT A DANGER OF INFECTION.
	WARNING! WHEN MONITORING PACEMAKER PATIENTS, HEART RATE METERS MAY CONTINUE TO COUNT THE PACEMAKER RATE DURING OCCURRENCES OF CARDIAC ARREST OR SOME ARRHYTHMIAS. DO NOT RELY ENTIRELY UPON HEART RATE METER ALARMS. KEEP PACEMAKER PATIENTS UNDER CLOSE SURVEILLANCE. SEE <u>16.7.1 HEART RATE MONITORING SPECIFICATIONS</u> ON <u>PAGE 214</u> FOR DISCLOSURE OF THE PACEMAKER PULSE REJECTION CAPABILITY OF MOVES[®] SLC[™].
	WARNING! IMPROPER OPERATION OF THE MOVES[®] SLC[™] SYSTEM COULD ENDANGER A PATIENT!
	WARNING! MOVES[®] SLC[™] IS INTENDED FOR USE ON ONE PATIENT AT A TIME. FOR EXAMPLE, IT SHOULD NOT BE USED TO VENTILATE ONE PATIENT WHILE MONITORING ANOTHER.
	WARNING! DO NOT CONNECT ANY SENSORS, MONITORS, OR THE BREATHING CIRCUIT TO THE PATIENT WHILE PERFORMING SYSTEM TESTS! DOING SO COULD ENDANGER THE PATIENT!
	WARNING! FAILURE TO CHANGE THE VENTILATOR CARTRIDGE WHEN INDICATED MAY LEAD TO THE PATIENT'S SUFFERING FROM AN INCREASE IN INSPIRED CO₂.

SYMBOL	PATIENT-SPECIFIC WARNINGS
	WARNING! ALWAYS CARRY ALTERNATE MEANS OF VENTILATING, SUCTIONING, AND OXYGENATING THE PATIENT.
	WARNING! ALWAYS CARRY BACKUPS OF CONSUMABLES SUCH AS CARTRIDGES AND FILTERS.
	WARNING! NEVER LEAVE A PATIENT UNATTENDED WHEN RUNNING MOVES[®] SLC[™].
	EXTREME WARNING! WHEN USING FLUID FILLED PRESSURE TRANSDUCERS TO MONITOR INTRACRANIAL PRESSURE (ICP), MAKE SURE THAT THE TRANSDUCER AND THE LINE CONNECTING TO THE PATIENT'S DRAIN ARE FREE OF ANY AIR BUBBLES!
	EXTREME WARNING! AFTER COMPLETING FILLING THE ICP TRANSDUCER AND THE LINE, DISCONNECT THE FLUID BAG FROM THE TRANSDUCER, AND CAP THE END WITH THE STERILE CAP PRIOR TO CONNECTING THE TRANSDUCER TO THE PATIENT'S BRAIN!
	EXTREME WARNING! NEVER FLUSH THE ICP TRANSDUCER WHILE CONNECTED TO THE PATIENT!
	EXTREME WARNING! FAILURE TO OBSERVE THE PREVIOUS THREE PRECAUTIONS MAY RESULT IN SERIOUS INJURY OR DEATH!

8.5 Masimo Rainbow SET[®] Pulse CO-Oximeter Warnings

Table 12: Masimo Rainbow SET[®] Pulse CO-Oximeter Warnings

SYMBOL	MASIMO RAINBOW SET [®] PULSE CO-OXIMETER WARNINGS
	WARNING! PULSE RATE MEASUREMENT IS BASED ON THE OPTICAL DETECTION OF A PERIPHERAL FLOW PULSE AND THEREFORE MAY NOT DETECT CERTAIN ARRHYTHMIAS. THE PULSE OXIMETER SHOULD NOT BE USED AS A REPLACEMENT OR SUBSTITUTE FOR ECG-BASED ARRHYTHMIA ANALYSIS.
	WARNING! A PULSE CO-OXIMETER SHOULD BE CONSIDERED AN EARLY WARNING DEVICE. AS A TREND TOWARDS PATIENT HYPOXEMIA IS INDICATED, BLOOD SAMPLES SHOULD BE ANALYZED BY LABORATORY INSTRUMENTS TO COMPLETELY UNDERSTAND THE PATIENT'S CONDITION.
	WARNING! FOR MEASUREMENTS OF HIGH OR LOW SPHB READINGS, BLOOD SAMPLES SHOULD BE ANALYZED BY LABORATORY INSTRUMENTS TO COMPLETELY UNDERSTAND THE PATIENT'S CONDITION.

SYMBOL	MASIMO RAINBOW SET [®] PULSE CO-OXIMETER WARNINGS
	WARNING! SPO₂ IS EMPIRICALLY CALIBRATED TO FUNCTIONAL ARTERIAL OXYGEN SATURATION IN HEALTHY ADULT VOLUNTEERS WITH NORMAL LEVELS OF CARBOXYHEMOGLOBIN (COHB) AND METHEMOGLOBIN (METHB). A PULSE OXIMETER CAN NOT MEASURE ELEVATED LEVELS OF COHB OR METHB. INCREASES IN EITHER COHB OR METHB WILL AFFECT THE ACCURACY OF THE SPO₂ MEASUREMENT. • FOR INCREASED COHB: COHB LEVELS ABOVE NORMAL TEND TO INCREASE THE LEVEL OF SPO₂. THE LEVEL OF INCREASE IS APPROXIMATELY EQUAL TO THE AMOUNT OF COHB THAT IS PRESENT. • NOTE: HIGH LEVELS OF COHB MAY OCCUR WITH A SEEMINGLY NORMAL SPO₂. WHEN ELEVATED LEVELS OF COHB ARE SUSPECTED, LABORATORY ANALYSIS (CO-OXIMETRY) OF A BLOOD SAMPLE SHOULD BE PERFORMED. • FOR INCREASED METHB: THE SPO₂ MAY BE DECREASED BY LEVELS OF METHB OF UP TO APPROXIMATELY 10% TO 15%. AT HIGHER LEVELS OF METHB, THE SPO₂ MAY TEND TO READ IN THE LOW TO MID 80S. WHEN ELEVATED LEVELS OF METHB ARE SUSPECTED, LABORATORY ANALYSIS (CO-OXIMETRY) OF A BLOOD SAMPLE SHOULD BE PERFORMED.
	WARNING! INTERFERING SUBSTANCES: DYES, OR ANY SUBSTANCE CONTAINING DYES, THAT CHANGE USUAL BLOOD PIGMENTATION MAY CAUSE ERRONEOUS READINGS.
	WARNING! HEMOGLOBIN SYNTHESIS DISORDERS MAY CAUSE ERRONEOUS SPHB READINGS.
	WARNING! ELEVATED LEVELS OF TOTAL BILIRUBIN MAY LEAD TO INACCURATE SPO₂, SPMET, SPCO, SPHB, AND SPOC MEASUREMENTS.
	WARNING! MOTION ARTIFACT MAY LEAD TO INACCURATE SPMET, SPCO, SPHB, AND SPOC MEASUREMENTS.
	WARNING! SEVERE ANEMIA MAY CAUSE ERRONEOUS SPO₂ READINGS.
	WARNING! VERY LOW ARTERIAL OXYGEN SATURATION (SPO₂) LEVELS MAY CAUSE INACCURATE SPCO AND SPMET MEASUREMENTS.
	WARNING! WITH VERY LOW PERFUSION AT THE MONITORED SITE, THE READINGS MAY READ LOWER THAN CORE ARTERIAL OXYGEN SATURATION.
	WARNING! DO NOT USE TAPE TO SECURE THE SENSOR TO THE SITE; THIS CAN RESTRICT BLOOD FLOW AND CAUSE INACCURATE READINGS. USE OF ADDITIONAL TAPE CAN CAUSE SKIN DAMAGE OR DAMAGE THE SENSOR.
	WARNING! IF THE SENSOR IS WRAPPED TOO TIGHTLY, OR SUPPLEMENTAL TAPE IS USED, VENOUS CONGESTION / PULSATIONS MAY OCCUR, CAUSING ERRONEOUS READINGS.

SYMBOL	MASIMO RAINBOW SET [®] PULSE CO-OXIMETER WARNINGS
	WARNING! VENOUS CONGESTION MAY CAUSE UNDER READING OF ACTUAL ARTERIAL OXYGEN SATURATION. THEREFORE, ASSURE PROPER VENOUS OUTFLOW FROM MONITORED SITE. SENSOR SHOULD NOT BE BELOW HEART LEVEL (E.G., SENSOR ON HAND OF A PATIENT IN A BED WITH ARM DANGLING TO THE FLOOR).
	WARNING! VENOUS PULSATIONS MAY CAUSE ERRONEOUS LOW READINGS (E.G., TRICUSPID VALVE REGURGITATION).
	WARNING! LOSS OF PULSE SIGNAL CAN OCCUR WHEN: • THE SENSOR IS TOO TIGHT. • THE PATIENT HAS HYPOTENSION, SEVERE VASOCONSTRICTION, SEVERE ANEMIA, OR HYPOTHERMIA. • THERE IS ARTERIAL OCCLUSION PROXIMAL TO THE SENSOR. • THE PATIENT IS IN CARDIAC ARREST OR IS IN SHOCK.
	WARNING! THE PULSATIONS FROM INTRA-AORTIC BALLOON SUPPORT CAN BE ADDITIVE TO THE PULSE RATE ON THE OXIMETER PULSE RATE DISPLAY. VERIFY PATIENT'S PULSE RATE AGAINST THE ECG HEART RATE.
	WARNING! MISAPPLIED SENSORS OR SENSORS THAT BECOME PARTIALLY DISLODGED MAY CAUSE EITHER OVER OR UNDER READING OF ACTUAL ARTERIAL OXYGEN SATURATION.
	WARNING! AVOID PLACING THE SENSOR ON ANY EXTREMITY WITH AN ARTERIAL CATHETER OR BLOOD PRESSURE CUFF.
	WARNING! HIGH INTENSITY EXTREME LIGHTS (INCLUDING PULSATING STROBE LIGHTS) DIRECTED ON THE SENSOR MAY NOT ALLOW THE PULSE CO-OXIMETER TO OBTAIN READINGS.
	WARNING! THE PULSE CO-OXIMETER CAN BE USED DURING DEFIBRILLATION, BUT THE READINGS MAY BE INACCURATE FOR UP TO 20 SECONDS.
	WARNING! BEFORE USE, CAREFULLY READ THE MASIMO SENSOR <i>DIRECTIONS FOR USE</i>.
	WARNING! TISSUE DAMAGE CAN BE CAUSED BY INCORRECT APPLICATION OR USE OF A SENSOR, FOR EXAMPLE BY WRAPPING THE SENSOR TOO TIGHTLY. INSPECT THE SENSOR SITE AS DIRECTED IN THE SENSOR'S <i>DIRECTIONS FOR USE</i> TO ENSURE SKIN INTEGRITY AND CORRECT POSITIONING AND ADHESION OF THE SENSOR.
	WARNING! THE PULSE CO-OXIMETER IS NOT INTENDED FOR USE AS AN APNEA MONITOR.
	WARNING! TO AVOID CROSS CONTAMINATION USE ONLY MASIMO SINGLE USE SENSORS ON THE SAME PATIENT.

SYMBOL	MASIMO RAINBOW SET [®] PULSE CO-OXIMETER WARNINGS
	WARNING! UNLESS OTHERWISE SPECIFIED, DO NOT STERILIZE SENSORS OR PATIENT CABLES BY IRRADIATION, STEAM, AUTOCLAVE OR ETHYLENE OXIDE. SEE THE CLEANING INSTRUCTIONS IN THE <i>DIRECTIONS FOR USE</i> FOR THE MASIMO RE-USEABLE SENSORS.
	WARNING! DO NOT USE THE PULSE CO-OXIMETER OR OXIMETRY SENSORS DURING MAGNETIC RESONANCE IMAGING (MRI) SCANNING. INDUCED CURRENT COULD POTENTIALLY CAUSE BURNS. THE PULSE CO-OXIMETER MAY AFFECT THE MRI IMAGE, AND THE MRI UNIT MAY AFFECT THE ACCURACY OF THE OXIMETRY MEASUREMENTS.
	WARNING! IF USING PULSE CO-OXIMETRY DURING FULL BODY IRRADIATION, KEEP THE SENSOR OUT OF THE IRRADIATION FIELD. IF THE SENSOR IS EXPOSED TO THE IRRADIATION, THE READING MIGHT BE INACCURATE, OR THE UNIT MIGHT READ ZERO FOR THE DURATION OF THE ACTIVE IRRADIATION PERIOD.
	WARNING! EXERCISE CAUTION WHEN APPLYING A SENSOR TO A SITE WITH COMPROMISED SKIN INTEGRITY. APPLYING TAPE OR PRESSURE TO SUCH A SITE MAY REDUCE CIRCULATION AND/OR CAUSE FURTHER SKIN DETERIORATION.
	WARNING! CIRCULATION DISTAL TO THE SENSOR SITE SHOULD BE CHECKED ROUTINELY.
	WARNING! A FUNCTIONAL TESTER CANNOT BE UTILIZED TO ASSESS THE ACCURACY OF THE PULSE CO-OXIMETER OR ANY SENSORS.
	WARNING! DO NOT MODIFY OR ALTER A PULSE CO-OXIMETER SENSOR IN ANY WAY. ALTERATIONS OR MODIFICATION MAY AFFECT PERFORMANCE AND/OR ACCURACY.
	WARNING! DO NOT USE DAMAGED SENSORS OR PATIENT CABLES. DO NOT USE A SENSOR OR PATIENT CABLE WITH EXPOSED OPTICAL OR ELECTRICAL COMPONENTS.
	WARNING! DO NOT ATTEMPT TO REPROCESS, RECONDITION OR RECYCLE ANY MASIMO SENSORS OR PATIENT CABLES AS THESE PROCESSES MAY DAMAGE THE ELECTRICAL COMPONENTS, POTENTIALLY LEADING TO HARM.
	WARNING! EXPLOSION HAZARD. DO NOT USE THE PULSE CO-OXIMETER IN THE PRESENCE OF FLAMMABLE ANESTHETICS OR OTHER FLAMMABLE SUBSTANCE IN COMBINATION WITH AIR, OXYGEN-ENRICHED ENVIRONMENTS, OR NITROUS OXIDE.
	WARNING! AS WITH ALL MEDICAL EQUIPMENT, CAREFULLY ROUTE PATIENT CABLING TO REDUCE THE POSSIBILITY OF PATIENT ENTANGLEMENT OR STRANGULATION.
	WARNING! ALWAYS REMOVE THE SENSOR FROM THE PATIENT AND COMPLETELY DISCONNECT THE PATIENT FROM THE PULSE CO-OXIMETER BEFORE BATHING THE PATIENT.
	WARNING! DO NOT IMMERSE THE SENSOR OR PATIENT CABLE IN WATER, SOLVENTS, OR CLEANING SOLUTIONS (THE SENSORS AND CONNECTORS ARE NOT WATERPROOF).

SYMBOL	MASIMO RAINBOW SET [®] PULSE CO-OXIMETER WARNINGS
	WARNING! INTRAVASCULAR DYES OR EXTERNALLY APPLIED COLORING (SUCH AS NAIL POLISH) MAY LEAD TO INACCURATE SPO₂ MEASUREMENTS.
	WARNING! PATIENT SAFETY – IF A SENSOR IS DAMAGED IN ANY WAY, DISCONTINUE USE IMMEDIATELY.
	WARNING! FAILURE TO APPLY THE SENSOR PROPERLY MAY CAUSE INCORRECT MEASUREMENTS.
	WARNING! USE ONLY MASIMO SENSORS FOR PULSE OXIMETRY OR PULSE CO-OXIMETRY MEASUREMENTS.

8.6 General Cautions

Table 13: General Cautions

SYMBOL	GENERAL CAUTION
	CAUTION! OPERATION OF MOVES[®] SLC[™] OUTSIDE OF SPECIFIED LIMITS MAY CAUSE INACCURATE RESULTS.
	CAUTION! OPERATING SUCTION WITHOUT THE SUPPLIED FILTER, CANISTER, AND HOSES WILL RESULT IN PERMANENT FAILURE OF THE SUCTION SYSTEM AND THE OXYGEN CONCENTRATOR.
	CAUTION! MOVES[®] SLC[™] SHOULD NOT BE USED ADJACENT TO OR STACKED WITH OTHER EQUIPMENT. IF ADJACENT OR STACKED USE IS NECESSARY, MOVES[®] SLC[™] SHOULD BE OBSERVED TO VERIFY NORMAL OPERATION IN THE CONFIGURATION IN WHICH IT WILL BE USED.
	CAUTION! THE USE OF ACCESSORIES AND CABLES OTHER THAN THOSE SPECIFIED, WITH THE EXCEPTION OF ACCESSORIES AND CABLES QUALIFIED AND SOLD BY THORNHILL RESEARCH INC., MAY RESULT IN INCREASED EMISSIONS OR DECREASED IMMUNITY OF MOVES[®] SLC[™] AND MAY CAUSE THE SYSTEM TO BE NON-COMPLIANT WITH THE REQUIREMENTS OF IEC 60601-1-2:2007 (ED 3.0).
	CAUTION! SOME PATIENT MONITORING ACCESSORIES MAY NOT FUNCTION PROPERLY OR MAY RELAY INACCURATE READINGS IF OPERATED OUTSIDE OF THEIR NORMAL OPERATING CONDITIONS.
	CAUTION! USE ONLY THE ACCESSORIES THAT ARE PROVIDED WITH THE MOVES[®] SLC[™] SYSTEM!
	CAUTION! THE MOVES[®] SLC[™] ACCESSORIES PROVIDED WITH THE MOVES[®] SLC[™] SYSTEM ARE FOR USE ONLY WITH THE MOVES[®] SLC[™] SYSTEM.

SYMBOL	GENERAL CAUTION
	CAUTION! BEFORE INSTALLING A HYDROCARBON FILTER, CHECK THE FOUR-DIGIT DATE CODE PRINTED ON THE CARTRIDGE. THE CARTRIDGE LABEL IS STAMPED WITH FOUR CHARACTERS "XXYY", WHERE "XX" IS THE WEEK OF THE YEAR AND "YY" IS THE YEAR. A CARTRIDGE MORE THAN ONE YEAR OLD SHOULD BE DISCARDED SINCE IT MAY DEGRADE THE PERFORMANCE OF OR CAUSE DAMAGE TO THE MOVES [®] SLC [™] OXYGEN CONCENTRATOR.
	CAUTION! OPERATING MOVES [®] SLC [™] WITHOUT A HYDROCARBON FILTER WILL DAMAGE THE UNIT. DO NOT OPERATE WITHOUT A HYDROCARBON FILTER!
	CAUTION! A MOVES [®] SLC [™] UNIT SHOULD NEVER BE PUT INTO TRANSPORT SERVICE WITH LESS THAN A 95% CHARGE IN BOTH BATTERIES.
	CAUTION! NEVER CHARGE BATTERIES IN AMBIENT TEMPERATURES BELOW 32°F (0°C) OR ABOVE 104°F (40°C).
	CAUTION! IF INTENDING TO RUN ON BATTERIES, ENSURE THAT THERE IS SUFFICIENT POWER FOR THE LENGTH OF TIME REQUIRED, OR REPLACE THE BATTERIES.
	CAUTION! ONLY SELECT "SKIP TESTS" IF A SYSTEM TEST HAS PREVIOUSLY BEEN PERFORMED AND PASSED, AND THE SYSTEM HAS JUST BEEN RESTARTED DUE TO POWER FAILURE (E.G., LOSS OF BATTERY OR EXTERNAL POWER).
	CAUTION! THE SUCTION WAND SUPPLIED WITH THE MOVES [®] SLC [™] IS DESIGNED TO MINIMIZE TOTAL OCCLUSION. IT IS STRONGLY ADVISED THAT ONLY THE WAND SUPPLIED BE USED WITH THE MOVES [®] SLC [™] SYSTEM.
	CAUTION! ACCURACY OF ANY BLOOD PRESSURE MEASUREMENT MAY BE AFFECTED IF MOVES [®] SLC [™] IS USED OR STORED OUTSIDE THE RELEVANT TEMPERATURE OR HUMIDITY RANGES DESCRIBED IN THE ENVIRONMENTAL SPECIFICATIONS (SEE SECTION 16.6.7 ON PAGE 214).
	CAUTION! SUBSTITUTION OF COMPONENTS DIFFERENT FROM THOSE SUPPLIED WITH THE MOVES [®] SLC [™] MAY RESULT IN MEASUREMENT ERROR.
	CAUTION! THE SUCTION CANISTER AND SUCTION FILTER ARE INTENDED FOR SINGLE USE ONLY AND SHOULD BE DISPOSED OF IN ACCORDANCE WITH LOCAL BIOHAZARD REGULATIONS.
	CAUTION! BREATHING CIRCUITS, SAMPLE LINES (BUT NOT THE NAFION TUBE) AND FILTERS, BREATHING CARTRIDGES, ECG ADHESIVE SENSOR PADS, ABP/CVP/ICP TRANSDUCER AND SUCTION WAND AND TUBING ARE DISPOSABLE AND SHOULD BE DISPOSED OF IN ACCORDANCE WITH LOCAL BIOHAZARD REGULATIONS.
	CAUTION! ALL MOVES [®] SLC [™] ACCESSORIES' PACKAGING AND DISPOSABLE ACCESSORIES SHOULD BE DISPOSED OF RESPONSIBLY IN ACCORDANCE WITH LOCAL WASTE DISPOSAL STANDARDS AND / OR LOCAL BIOHAZARD REGULATIONS.

SYMBOL	GENERAL CAUTION
	CAUTION! DO NOT SUBMERGE THE MOVES [®] SLC [™] OR POUR CLEANING LIQUIDS OVER OR INTO THE MOVES [®] SLC [™] .
	CAUTION! THE LABEL ON THE PACKAGE OF THE VENTILATOR CARTRIDGE CONTAINS AN EXPIRY DATE. ALWAYS CHECK THE EXPIRY DATE ON THE VENTILATOR CARTRIDGE BEFORE USING IT TO MAKE SURE THAT THE VENTILATOR CARTRIDGE HAS NOT EXPIRED. AS WELL, MONITOR SPARE CARTRIDGES WITH REGARD TO THEIR REMAINING "SHELF LIFE".
	CAUTION! THE SURFACE OF THE MOVES [®] SLC [™] SYSTEM CAN BECOME HOT, ESPECIALLY IF IT IS BEING OPERATED IN DIRECT SUNLIGHT. CARE SHOULD BE TAKEN WHEN TOUCHING OR CONTACTING THE SURFACE OF THE MOVES [®] SLC [™] SYSTEM.
	CAUTION! ECG CABLES SHOULD BE DISCARDED AND REPLACED AFTER TWO (2) YEARS OF CONTINUOUS USE. CHECK CABLE USE BY RECORDING THE DATE THE CABLE WAS FIRST USED.
	CAUTION! CHECK THE EXPIRY DATE ON THE ECG ELECTRODES PACKAGE BEFORE USING. ELECTRODES ARE GOOD FOR 45 DAYS ONCE PACKAGE IS OPENED.
	CAUTION! WHEN YOU DETACH THE MOVES [®] SLC [™] CLAMPS, AND RETURN THEM TO THE ACCESSORIES CASE, BE CAREFUL NOT TO PINCH THE WIRES ATTACHED TO THE PINS IN THE CLAMP APPARATUS AS THIS CAN CAUSE WEAR, ABRADING, AND EVENTUAL BREAKAGE OF THE WIRES.
	CAUTION! WHEN BATTERIES ARE DISCHARGED AND LEFT IN MOVES [®] SLC [™] FOR A PROLONGED PERIOD, A <u>COMPLETELY</u> DISCHARGED BATTERY (NO LED LIGHTS) CAN RESULT. THE BATTERY CAN STILL BE RECHARGED, BUT IT MAY TAKE MORE THAN THE NORMAL 2.5 HOURS. IT HAS BEEN OBSERVED TO TAKE ANYWHERE FROM 6 TO 48 HOURS TO FULLY CHARGE.
	CAUTION! IF MOVES [®] SLC [™] IS EXPOSED TO SIGNIFICANT AMOUNTS OF SAND OR DUST, IT SHOULD BE CLEANED BY AN AUTHORIZED TECHNICIAN IN ACCORDANCE WITH THE MAINTENANCE MANUAL.
	CAUTION! THE VENTILATOR BREATHING CIRCUIT AND THE SAMPLING FILTER CONNECTED TO THE NAFION TUBING SHOULD BE INSPECTED EVERY FOUR (4) HOURS FOR CONDENSATION. THE CIRCUIT AND FILTER SHOULD BE CHANGED AFTER 12 HOURS OF CONTINUOUS USE.
	CAUTION! NO LUBRICANTS OTHER THAN THOSE RECOMMENDED BY THE MANUFACTURER SHALL BE USED ON THE MOVES [®] SLC [™] .

8.7 Electrical Cautions

Table 14: Electrical Cautions

SYMBOL	ELECTRICAL CAUTION
	CAUTION! MOVES [®] SLC [™] IS INTENDED FOR USE BY HEALTHCARE PROFESSIONALS ONLY. MOVES [®] SLC [™] MAY CAUSE RADIO INTERFERENCE OR MAY DISRUPT THE OPERATION OF NEARBY EQUIPMENT. IT MAY BE NECESSARY TO TAKE MITIGATION MEASURES, SUCH AS REORIENTING OR RELOCATING MOVES [®] SLC [™] OR SHIELDING THE LOCATION.
	CAUTION! MEDICAL ELECTRICAL EQUIPMENT NEEDS SPECIAL PRECAUTIONS REGARDING EMC (ELECTROMAGNETIC COMPATIBILITY) AND NEEDS TO BE INSTALLED AND PUT INTO SERVICE ACCORDING TO THE EMC INFORMATION PROVIDED IN THE MOVES [®] SLC [™] OPERATOR'S MANUAL.
	CAUTION! THE POWER SUPPLY/CHARGER CORD IS A SPECIAL MEDICAL-GRADE POWER CORD AND SHOULD NOT BE REPLACED WITH A NON-MOVES [®] SLC [™] SUPPLIED PART.
	CAUTION! PORTABLE AND MOBILE RF (RADIO FREQUENCY) COMMUNICATIONS EQUIPMENT CAN AFFECT MEDICAL ELECTRICAL EQUIPMENT.
	CAUTION! INTERFERENCE MAY OCCUR IN THE VICINITY OF KNOWN RADIO FREQUENCY (RF) TRANSMITTING DEVICES AND EQUIPMENT MARKED WITH THE FOLLOWING SYMBOL: 
	CAUTION! OXIMETER READINGS MAY BE AFFECTED BY THE USE OF AN ELECTROSURGICAL UNIT (ESU).

8.8 Patient-Specific Cautions

Table 15: Patient-Specific Cautions

SYMBOL	PATIENT-SPECIFIC CAUTION
	CAUTION! ALL COMPONENTS IN MOVES [®] SLC [™] THAT COME IN CONTACT WITH RESPIRATORY GASES ARE NON-STERILE.

8.9 General Safety

1. This equipment should be operated only by a trained practitioner in emergency transport medicine.
2. The operator must inspect the MOVES[®] SLC[™] unit and all accessories for visible physical damage (cracks, holes, leaks; missing components, structural hardware and protective covers; signs of tampering, etc.) prior to each use. All system tests must be completed prior to connecting the MOVES[®] SLC[™] unit to a patient.
3. The operator must be fully familiar with the contents of this manual before operating the MOVES[®] SLC[™] system.
4. The MOVES[®] SLC[™] unit must be serviced only by qualified personnel. There are no user serviceable parts inside the MOVES[®] SLC[™] unit.
5. The MOVES[®] SLC[™] system must not be used for any purpose other than as stated in section 3.2 *MOVES[®] SLC[™] Intended Use* on page 15 .
6. The MOVES[®] SLC[™] unit must always be used and stored in accordance with the environmental specifications listed in subsection 16.6.7 *Environmental Specifications* on page 214. Any storage or use outside of these conditions may cause system degradation and harm to patients. After any such storage or use, the MOVES[®] SLC[™] unit must be serviced by qualified personnel. Do not use the MOVES[®] SLC[™] unit in contaminated environments.
7. The operator must ensure that the MOVES[®] SLC[™] Unit and Accessories Cases contain all the necessary components for successful use. See *Section 9.1: MOVES[®] SLC[™] System Contents* on page 53 .
8. Always ensure that system batteries are fully charged, and always maintain alternative methods of ventilation, O₂ supplementation and suction in case of failure of part or all of the MOVES[®] SLC[™] system.
9. If the MOVES[®] SLC[™] unit and accessories are placed in storage, they must be inspected and serviced every 12 months at a minimum.
10. The MOVES[®] SLC[™] system must be operated using only accessories specified in this manual, and supplied by the manufacturer, or their equivalent.
11. The accessories provided with the MOVES[®] SLC[™] system are for use only with the MOVES[®] SLC[™] system. Use of MOVES[®] SLC[™] accessories with another system may result in failure of the accessories and present a risk to and/or harm a patient, operator or bystander.
12. The MOVES[®] SLC[™] system must not be connected to a patient unless the operator has thoroughly read and understood this manual, the operator is qualified or is under qualified supervision, all system tests have been performed, and all of the Safety Procedures contained within this manual have been read and adhered to.
13. The MOVES[®] SLC[™] unit should be placed on *Monitor Only* mode when it is operating in 'Ventilator' or 'O₂ Supplement' mode and filters or cartridges need to be changed.
14. The MOVES[®] SLC[™] unit must not be operated without filters and cartridges present and properly installed.
15. Failure to change the ventilator cartridge when indicated will lead to hypercarbia. Alternative ventilation and/or external O₂ supplementation may be required while changing the cartridge.
16. Cartridges are for single patient use only. The Ventilator Cartridge will require replacement after at least every 2 hours of continuous use under standard temperature and pressure conditions. Different environmental and patient conditions may reduce the life of the Ventilator Cartridge. Always have a spare cartridge and an alternative ventilation method available.
17. The Ventilator Cartridge has an expiry date listed on the package label. Before using the Ventilator Cartridge, always check the expiry date to make sure that the cartridge has not expired.
18. Always have available an alternative means of suction in the event of a power failure, mechanical failure or serious occlusion in the suction circuit.
19. The MOVES[®] SLC[™] system is not recommended for use in treating patients with pacemakers or other implanted medical devices.
20. If the MOVES[®] SLC[™] unit should experience ingress of any particulate matter (dirt, dust, sand, etc.) through suction, fresh gas inlet, ventilator driving gas or valve block ports, the unit should be cleaned and serviced by authorized personnel.

21. Should the MOVES[®] SLC[™] unit experience ingress of any liquids or become contaminated by bodily fluids through suction, fresh gas inlet, air intake, ventilator driving gas or valve block ports, the unit should be transferred to authorized personnel for cleaning, sterilization and/or component replacement.
22. Should the MOVES[®] SLC[™] unit experience a sudden shock (by being dropped or roughly handled), or exposed externally to high voltage, the unit should be inspected and serviced by authorized personnel.
23. The MOVES[®] SLC[™] system outputs concentrated oxygen. No smoking or open flame is permitted near the unit.
24. Do not use oil or grease on or near the MOVES[®] SLC[™] unit or its components.
25. The surface of the MOVES[®] SLC[™] system can become hot, especially if it is being operated in direct sunlight. Care should be taken when touching or contacting the surface of the MOVES[®] SLC[™] system.

8.10 Electrical Safety

1. Connect the MOVES[®] SLC[™] power supply / charger only to an AC source from 100–240V, 50–60 Hz, 5.5 A max. Fluctuations in voltage and current can have adverse effects on the performance of the MOVES[®] SLC[™] system.
2. Do not modify the power supply provided with the MOVES[®] SLC[™] system with additional voltage regulators or similar equipment. The MOVES[®] SLC[™] power supply/charger detects the supply automatically.
3. Do not route power cables immediately adjacent to patient connection cables. Power cables can produce voltage transients (“crosstalk”) that seriously affect data collection cables, especially ECG sensor connections. Voltage transients can resemble ECG readings.
4. Do not operate the MOVES[®] SLC[™] unit if the main AC cord or the power-supply cord to MOVES[®] SLC[™] unit shows any sign of damage such as frayed insulation, or if there are cracked or damaged plugs or receptacles, or a missing grounding plug, or if the ‘PWR’ light on the power supply/charger does not illuminate.
5. Do not operate the MOVES[®] SLC[™] system on battery power if:
 - a. A battery shows any sign of damage such as cracks, holes, or leakage.
 - b. Any battery is known to be defective.
 - c. The battery discharge indicator on either battery fails to illuminate.
 - d. The power source indicator on the startup screen of the MOVES[®] SLC[™] unit fails to detect a battery when a battery is installed, or it indicates a charge that is not in agreement with the indicators on the batteries.
 - e. The power source indicator on the startup screen of the MOVES[®] SLC[™] fails to show a decrease in charge while running.
 - f. The ‘PWR’ light on the power supply/charger does not illuminate when attempting to charge the batteries.
6. The main external AC cord is supplied with a grounding prong on the male connector to reduce the risk of electrical shock. Never remove or compromise this prong.
7. Use only the MOVES[®] SLC[™] power supply P/N 111422, batteries P/N 122985, and supplied external power cables to operate the MOVES[®] SLC[™] system.

8.11 Preparing for Emergency Operation

The operator must always be prepared for emergency situations such as power failure, mechanical failure or serious occlusions in the MOVES[®] SLC[™] or its accessories that may require alternative means of treating a patient.

The following guidelines should always be adhered to:

1. The operator must always have prepared two completely charged battery packs.
2. The operator must always have prepared an alternative means of ventilating a patient.

3. The operator must always have prepared an alternative means of providing a high concentration of oxygen to a patient.
4. The operator must always have prepared an alternative means of monitoring a patient's vital signs (e.g., heart rate, blood pressure, respiratory rate, etc.).
5. The operator must always have prepared an alternative means of applying suction to a patient.
6. When treating a patient and operating the MOVES[®] SLC[™] on external power, the operator must always have at least one charged battery installed in the MOVES[®] SLC[™] unit. In the event of an AC power failure, the MOVES[®] SLC[™] unit will immediately switch to battery power.

8.12 Radio Interference



CAUTION! INTERFERENCE MAY OCCUR IN THE VICINITY OF KNOWN RADIO FREQUENCY (RF) TRANSMITTING DEVICES AND EQUIPMENT MARKED WITH THE FOLLOWING SYMBOL:



8.13 Battery Handling

Please note the following important cautions regarding battery handling:

- The batteries should not be opened, destroyed nor incinerated since they may leak or rupture and release into the environment the ingredients they contain.
- The batteries are designed to be recharged. However, improperly charging a cell or battery may cause the cell or battery to flame.
- Use only the MOVES[®] SLC[™] battery charger (P/N 111422) and MOVES[®] SLC[™] battery charging procedures.
- NEVER disassemble a battery or bypass any safety device.
- Do not crush, pierce, or short battery terminals with conductive (i.e. metal) goods.
- Do not directly heat or solder.
- Do not throw into fire.

8.14 Battery Disposal

MOVES[®] SLC[™] batteries are based on Lithium Polymer (LiPo) chemistry. Always dispose of batteries in accordance with local municipal, state, and federal regulations.



WARNING! CHARGING AND DISCHARGING THE LITHIUM POLYMER BATTERIES WITH OTHER THAN MOVES[®] SLC[™] SYSTEM EQUIPMENT, AND/OR IMPROPER HANDLING, CAN RESULT IN FIRE, EXPLOSION, TOXIC GASES AND SMOKE.



WARNING! DO NOT OPERATE THE MOVES[®] SLC[™] SYSTEM WITH A DEFECTIVE BATTERY.

THIS PAGE DELIBERATELY LEFT BLANK.

9.0 Getting Started

The following section provides information and instructions on installing and connecting various parts and accessories and preparing the MOVES® SLC™ for activation.



CAUTION! THESE ACTIVITIES SHOULD BE CARRIED OUT ONLY BY AUTHORIZED / TRAINED PERSONNEL.

9.1 MOVES® SLC™ System Contents

9.1.1 MOVES® SLC™ Accessories Warnings and Cautions

Table 16: MOVES® SLC™ Accessories Warnings and Cautions

	WARNING! ALWAYS CARRY BACKUPS OF ACCESSORIES SUCH AS CARTRIDGES AND FILTERS.
	WARNING! IF LEFT IN A HOT ENVIRONMENT OR DIRECT SUNLIGHT FOR A CONSIDERABLE LENGTH OF TIME, THE MOVES® SLC™ ACCESSORIES CAN BECOME QUITE HOT. MONITORING OF TEMPERATURE IS RECOMMENDED. ALWAYS MAKE SURE THAT ACCESSORIES THAT WILL BE APPLIED DIRECTLY TO THE PATIENT ARE SUITABLE FOR SKIN CONTACT AND WILL NOT CAUSE BURNS.
	CAUTION! SOME PATIENT MONITORING ACCESSORIES MAY NOT FUNCTION PROPERLY, OR MAY RELAY INACCURATE READINGS, IF OPERATED OUTSIDE OF THEIR NORMAL OPERATING CONDITIONS.
	CAUTION! USE ONLY ACCESSORIES THAT ARE PROVIDED OR AUTHORIZED BY THORNHILL RESEARCH INC!
	CAUTION! THE MOVES® SLC™ ACCESSORIES PROVIDED BY THORNHILL RESEARCH INC. ARE FOR USE ONLY WITH THE MOVES® SLC™ SYSTEM.
	CAUTION! ALL MOVES® SLC™ ACCESSORIES' PACKAGING SHOULD BE DISPOSED OF RESPONSIBLY IN ACCORDANCE WITH LOCAL WASTE DISPOSAL STANDARDS AND / OR LOCAL BIOHAZARD REGULATIONS.




NOTE: Where applicable, refer to accessories documentation for specific instructions for use, warnings, storage and operating guidelines.



9.1.2 System Contents


Table 17: MOVES[®] SLC[™] System Contents

PICTURE	TRI P/N	DESCRIPTION
	122752	MOVES [®] SLC [™] System Manufacturer: Thornhill Research Inc.
	111422	Power Supply / Battery Charger (both together in single unit) and AC cable. Manufacturer: Thornhill Research Inc.



PICTURE	TRI P/N	DESCRIPTION
	122985	Battery (2) Manufacturer: Thornhill Research Inc.
	101238	Suction Canister Holder Manufacturer: Thornhill Research Inc.


PICTURE	TRI P/N	DESCRIPTION
	111462	Shoulder Strap Manufacturer: Thornhill Research Inc.
 Front Clamp – Rear View	124566	Front Clamp Manufacturer: Thornhill Research Inc. NOTE: The front clamp can be distinguished by its <u>flat</u> bottom portion near the locking mechanism. It is labeled FRONT.

PICTURE	TRI P/N	DESCRIPTION
 Front Clamp – Side View		
 Back Clamp – Rear View	124565	Back Clamp Manufacturer: Thornhill Research Inc. NOTE: The back clamp can be distinguished by its <u>curved</u> bottom portion near the locking mechanism. It is labeled BACK.

PICTURE	TRI P/N	DESCRIPTION
 Back Clamp – Side View		
	100832	Sensor, NIBP Cuff Small Adult 17–25 cm Manufacturer: SunTech Medical, PN # 98-0080-04
	100834	Sensor, NIBP Cuff Adult 23–33 cm Manufacturer: SunTech Medical, PN # 98-0080-06
	100836	Sensor, NIBP Cuff Large Adult 31–40 cm Manufacturer: SunTech Medical, PN # 98-0080-08
	100838	Sensor, NIBP Cuff Thigh 38-50 cm Manufacturer: SunTech Medical, PN # 98-0080-10

PICTURE	TRI P/N	DESCRIPTION
	124848	NIBP Extension Cable – 3 m (with Lemo connector) Manufacturer: SunTech Medical, PN # 91-0028-09
	101048	ABP/CVP/ICP Cable Set Manufacturer: Fogg System Company, PN # 7185-0798(A)

	125350	12 Lead ECG Cable Manufacturer: MedLab Medizinische Diagnosegeräte GmbH, PN # M03 008/668/001
	125352	LNCS Patient Cable Red LNC-10, 10 ft. (for SpO₂ Finger Clip, SpO₂ Ear Clip and SpO₂ Adhesive Sensors) Manufacturer: Masimo Corporation, Prod. Ref: Red LNC-10 PN # 2056
	125351	Adult SpO₂ (only) Finger Clip Manufacturer: Masimo Corporation Prod. Ref: LNCS DCI, 3 ft. PN # 1863 NOTE: Red LNC-10 adapter cable required.

	125355	SpO₂ Ear Clip Manufacturer: Masimo Corporation, Prod. Ref: LNC5 TC-I, 3 ft. PN # 1895 NOTE: Red LNC-10 adapter cable required.
	125357	SpO₂ Adhesive Sensor (adult) Manufacturer: Masimo Corporation, Prod. Ref: LNC5 Adeptx-3, 3 ft. PN # 1859 NOTE: Red LNC-10 adapter cable required
	125353	SpO₂ & SpHb Finger Clip Manufacturer: Masimo Corporation, Prod. Ref: Rainbow[®] DC-3 SC 360, 3 ft. PN # 2241 NOTE: Requires activation of optional SpHb parameter support on MOVES[®] SLC[™] System.


	125354 SpO₂, SpCO & SpMet Finger Clip Manufacturer: Masimo Corporation, Prod. Ref: Rainbow[®] DCI-dc3, 3 ft. PN # 2201 NOTE: Requires activation of optional SpCO and SpMet parameter support on MOVES[®] SLC[™] System.
	125732 Rainbow[®] Patient Cable Rainbow[®] RC-4 (for Adhesive Sensors), 4 ft. Manufacturer: Masimo Corporation, Prod. Ref: Rainbow[®] RC-4 PN # 2406
	125739 Adult SpO₂, SpHb & SpMet Adhesive Sensor Manufacturer: Masimo Corporation, Prod. Ref: Rainbow[®] R1 25 PN # 2416 Box of 10 NOTE: Rainbow[®] RC-4 adapter cable required. NOTE: Requires activation of optional SpHb parameter support on MOVES[®] SLC[™] System.


 A white, flexible adhesive sensor with a long, thin cable. The cable has the word "rainbow" printed on it. The sensor head is white and has a small black label with the Masimo logo and "Rainbow R25" text.	125740	Adult SpO₂, SpCO & SpMet Adhesive Sensor Manufacturer: Masimo Corporation, Prod. Ref: Rainbow[®] R25 PN # 2221 Box of 10 NOTE: Rainbow[®] RC-4 adapter cable required. NOTE: Requires activation of optional SpCO and SpMet parameter support on MOVES[®] SLC[™] System.
 A metal component consisting of a circular ring attached to a rectangular base. The base has a threaded rod passing through it, with a spring and a nut visible. The component appears to be made of brass or a similar metal.	101113	Ancra Single Stud Manufacturer: Avio Diepen Inc., PN # 40340-20
 A coiled black cable with a silver-colored connector at one end. The cable is flexible and has a small label near the connector.	124849	Temperature Probe Manufacturer: Exacon Scientific, D-OS4A (adult with connector type A)


PICTURE	TRI P/N	DESCRIPTION
	124241	Nafion Tube Manufacturer: VacuMed, PN # 17049
	124826	MOVES® SLC™ Operator's Manual Manufacturer: Thornhill Research Inc.

9.2 MOVES® SLC™ Consumables

Table 18: MOVES® SLC™ Consumables

PICTURE	TRI P/N	DESCRIPTION
	124831	Ventilator Cartridge (Package of 6) Manufacturer: Thornhill Research Inc.
	101114	Suction Canister Manufacturer: Bemis Health Care, PN # 424410

PICTURE	TRI P/N	DESCRIPTION
	111458	Wand, Suction, Bulb Tip, Yankauer – FDA Manufacturer: Cardinal Healthcare, PN # K82
	125168	Wand, Suction, Bulb Tip, Yankauer – CE Manufacturer: HITECH Medical Co., Ltd. PN # HTD0106
	100915	Hydrocarbon Filter Manufacturer: Airgas Inc. PN # 815182
	101243	Tube, Suction (short, canister to system) – FDA Manufacturer: Cardinal Health Canada, PN # 52A
	125166	Tube, Suction (short, canister to system) – CE Manufacturer: HITECH Medical Co., Ltd. PN # HTD0207 NOTE: Only one of the above will be included depending on customer location.

PICTURE	TRI P/N	DESCRIPTION
	101244	Tube, Suction (long, canister to patient) – FDA Manufacturer: Cardinal Health Canada, PN # 66A
	125167	Tube, Suction (long canister to patient) – CE Manufacturer: HITECH Medical Co., Ltd. PN # HTD0203 NOTE: Only one of the above will be included depending on customer location.
	101210	Ventilator Breathing Circuit Manufacturer: Thornhill Research Inc. NOTE: Nafion tube is shown attached – this tube is NOT included in the breathing circuit packaging but should be attached prior to use.
	111520	Pressure Transducer for ABP/CVP/ICP Set – FDA Manufacturer: Edwards Lifesciences, PN # PX260
	125164	Pressure Transducer for ABP/CVP/ICP Set – CE Manufacturer: Edwards Lifesciences, PN # T100209A NOTE: Only one of the above will be included depending on customer location.


PICTURE	TRI P/N	DESCRIPTION
	124834	ECG Adhesive Pads Manufacturer: Kendall Medi-Trace[™] 530, P/N # 31013926 Pack of 30 pads

9.3 Attaching the Shoulder Strap to MOVES[®] SLC[™]



NOTE: The shoulder strap supplied with the MOVES[®] SLC[™] is attached using two Ancra Single Stud anchors. To attach the shoulder strap, the two stud anchors supplied with the MOVES[®] SLC[™] must be attached first and then the shoulder strap is clipped to them.


1. To begin, depress the pin in one of the Ancra Single Stud anchors.




2. Slide the anchor along one of the short rails at either end of MOVES[®] SLC[™] until it is in the center position.



3. The anchor will lock into place in the space between the two circular openings in the row (as shown in the illustration at the right).



4. Repeat the above actions for a second stud at the other end of MOVES[®] SLC[™].
5. Attach the clips on the two ends of the shoulder strap to the two anchors now fastened to MOVES[®] SLC[™].



6. Strap attached to MOVES[®] SLC[™].



9.4 Attaching Clamps to the MOVES[®] SLC[™]

Clamps can be attached to the MOVES[®] SLC[™] which can then enable it to be clamped to a gurney, bed frame or stretcher to enable patient transport.



WARNING! THE CLAMPS HAVE NUMEROUS MOVING PARTS THAT MAY PRESENT A PINCHING OR CRUSHING HAZARD. ALWAYS USE CAUTION WHEN HANDLING BOTH THE FRONT AND BACK CLAMPS.

1. To begin, orient the MOVES[®] SLC[™] with the FRONT to the operator's left.

NOTE: The FRONT is where the ventilator cartridge is inserted and the patient connections panel is located.




2. Next, lay the MOVES[®] SLC[™] on its side.




3. Locate the FRONT clamp for the MOVES[®] SLC[™].


NOTE: The word *FRONT* is embossed on it.




4. Undo the large wing bolt to its furthest extent (by turning it counterclockwise) to open the jaw of the clamp as wide as possible.




5. Orient the clamp as shown in the photo at right. It should nestle snugly against the side of the MOVES[®] SLC[™].



6. Insert the pin attached to the left of the clamp in the slot shown in the photo at right. Press the blue pin button while inserting the pin. Make sure the pin inserts fully into its intended slot (see photo following).




7. Now, insert the pin attached to the right of the clamp in the upper slot shown in the photo at right. Press the blue pin button while inserting the pin. Make sure the pin inserts fully into its intended slot.



8. The fully mounted FRONT clamp is shown in the photo at right. Repeat the above procedure to mount the BACK clamp at the other end of the MOVES[®] SLC[™].

NOTE: Do not worry if the clamps feel slightly loose. When they are clamped to a gurney, bed rail or stretcher, the lateral object to which they are clamped adds the required reinforcement.



9.5 Patient Monitoring Accessories

This section briefly describes how the patient monitoring accessories are to be connected to the MOVES[®] SLC[™] unit. The operator must thoroughly read all of the procedures contained within this manual and must fully understand the MOVES[®] SLC[™] system operation before connecting the system to a patient. Prior to connecting the system to a patient, the operator must have alternative methods of treating the patient available should a power or mechanical failure occur.



WARNING! DO NOT CONNECT ANY MONITORS TO PATIENT WHILE PERFORMING SYSTEM TESTS! DOING SO COULD ENDANGER THE PATIENT!



CAUTION! USE ONLY ACCESSORIES PROVIDED WITH THE MOVES® SLC™ SYSTEM!

All MOVES® SLC™ patient monitors are connected to the patient connection panel (pictured below) as indicated by the labels.




Figure 9-1 : Patient Connection Panel



NOTE: The unlabeled port between the Gas Sample port and the Temp 2 cover is the Calibration port.

Table 19: Patient Connector Labels and Accessories

CONNECTOR LABEL	ACCESSORY
GAS SAMPLE	Sampling line Luer connection
ECG	Electrocardiogram – 12 Lead – MedLab Medizinische Diagnosegeräte GmbH
NIBP	Non-invasive blood pressure – SunTech
Temp 1 & Temp 2	Temperature – Measurement Specialties/Meas400
IP 1–3	Ports can be interchangeably used for any of the following: ABP — Arterial Blood Pressure transducer (Edwards TruWave PX Series) CVP — Central Venous Pressure transducer (Edwards TruWave PX Series) ICP — Intracranial Pressure transducer (Edwards TruWave PX Series)
SpO ₂	Pulse CO-Oximeter – Masimo Medical
Barometric Sensing Port	No connectors are fitted to this port. Do Not Block.
Calibration Port	No connectors are fitted to this port. Do Not Block.

9.6 Installing the Ventilator Cartridge and Breathing Circuit



WARNING! LEAKS IN THE SAMPLING LINE CAN CAUSE LOW pCO₂ AND/OR LOW O₂ LEVELS.

9.6.1 About the Ventilator Cartridge

The Ventilator Cartridge is intended to be used with the MOVES[®] SLC[™] system to provide positive pressure ventilation for patients who are *intubated*. MOVES[®] SLC[™] recycles exhaled oxygen. This cartridge contains CO₂ absorbent material to remove CO₂ from any re-breathed gas. Because this material has a *shelf life*, there is an expiry date printed on the package label of the Ventilator Cartridge. Always check the expiry date on the Ventilator Cartridge before using it to make sure that the Ventilator Cartridge has not expired. As well, monitor your supply of spare cartridges with regard to their remaining “shelf life”.



NOTE: The Ventilator Cartridge should be kept in its packaging until use.

The Ventilator Cartridge is:

- For use with the ventilator breathing circuit, which includes the ventilator hoses, patient filter, endotracheal tube connector, and sampling line.
- For single patient use ONLY. The cartridge should be discarded and replaced between patients, or when MOVES[®] SLC[™] triggers an audible or visual alarm indicating that the level of CO₂ in the system is above 6 mmHg on inspiration.



WARNING! FAILURE TO CHANGE THE VENTILATOR CARTRIDGE WHEN INDICATED MAY LEAD TO THE PATIENT'S SUFFERING FROM AN INCREASE IN INSPIRED CO₂.





Figure 9-2: Ventilator Cartridge – Patient Connection Side (Front)




Figure 9-3: Ventilator Cartridge – MOVES[®] SLC[™] Connection Side (Back)

9.6.2 Ventilator Cartridge Installation Instructions

1. Turn the cartridge Release / Lock knob to the unlocked position with the arrow facing toward the user.



2. Lift and open the valve block cover door inside the cartridge cavity.



3. Insert the cartridge into the cavity as shown.



4. Turn the cartridge Release / Lock knob to the locked position with the arrow facing away from the user.



9.6.3 Installing the Hydrocarbon Filter



CAUTION! BEFORE INSTALLING A HYDROCARBON FILTER, CHECK THE FOUR-DIGIT DATE CODE PRINTED ON THE CARTRIDGE. THE CARTRIDGE LABEL IS STAMPED WITH FOUR CHARACTERS “XXYY”, WHERE “XX” IS THE WEEK OF THE YEAR AND “YY” IS THE YEAR OF MANUFACTURE. A CARTRIDGE MORE THAN ONE YEAR OLD SHOULD BE DISCARDED SINCE IT MAY DEGRADE THE PERFORMANCE OF OR CAUSE DAMAGE TO THE MOVES[®] SLC[™] OXYGEN CONCENTRATOR.

1. Always install a hydrocarbon filter before powering up MOVES[®] SLC[™]. MOVES[®] SLC[™] will alert you with an alarm if it becomes clogged.
2. Insert the hydrocarbon filter as shown at right into the REAR panel of the MOVES[®] SLC[™] unit.



NOTE: If the hydrocarbon filter needs to be changed when MOVES[®] SLC[™] is in operation, a replacement filter will need to be readily available for quick insertion after the old one is removed.



CAUTION! OPERATING MOVES[®] SLC[™] WITHOUT A FILTER WILL DAMAGE THE UNIT. DO NOT OPERATE MOVES[®] SLC[™] WITHOUT A FILTER!

3. Rotate the hydrocarbon filter clockwise to install, counterclockwise to remove.



WARNING! IT MAY BE NECESSARY TO USE AN ALTERNATIVE MEANS OF OXYGEN SUPPLEMENTATION SHOULD THE HYDROCARBON FILTER REQUIRE REPLACEMENT WHILE A PATIENT IS BEING TREATED.



NOTE: Use only the filter cartridges supplied with MOVES[®] SLC[™]. These NIOSH-certified filters are specifically selected for use with MOVES[®] SLC[™]. These filter cartridges are not for use in atmospheres posing immediate danger to life or health or atmospheres containing less than 19.5% oxygen by volume.

9.6.4 Installing the Ventilator Breathing Circuit




Figure 9-4: Ventilator Breathing Circuit

Refer to the figure above when assembling the Ventilator Breathing circuit. Remove the **single use** circuit from its sealed package.



CAUTION! THE VENTILATOR BREATHING CIRCUIT AND THE SAMPLING FILTER CONNECTED TO THE NAFION TUBING SHOULD BE INSPECTED EVERY FOUR (4) HOURS FOR CONDENSATION. THE CIRCUIT AND FILTER SHOULD BE CHANGED AFTER 12 HOURS OF CONTINUOUS USE.



NOTE: The Nafion tube is **NOT** a single use item and **DOES NOT** come with the circuit. It should be housed in the MOVES[®] SLC[™] accessory case. **DO NOT DISCARD THE NAFION TUBE AFTER USE.** It has a sample filter to prevent contamination.



NOTE: Replace the Nafion tube at six (6) month intervals or as needed. The Nafion tube should be inspected prior to use for signs of physical damage including cracking and kinking. Premature failure of the Nafion tube can present as either an occluded sample line or, more commonly, as a leaky sample line that will produce a dampened PCO₂ trace and lower PetCO₂ readings.


1. Make sure there is a sample filter connected to the “Y” piece. If there is not, connect one.




NOTE: This sample filter, which is used to keep the Nafion tubing clean between uses, should be discarded after each use.




2. Connect one end of the Nafion tubing to the sample filter.




3. Connect the other end of the Nafion tubing to the Luer adapter, and then the Luer adapter to one end of the sampling line.



4. Connect the other end of the sampling line (with another sample filter on its end) to the 'GAS SAMPLE' port on the patient connections panel. The sample filter should be connected directly to the patient connection panel.



5. Attach the patient breathing circuit hoses to the patient hose ports on the breathing cartridge.



9.7 Delivering Supplementary Oxygen (O₂)



WARNING! WHEN USING O₂ SUPPLEMENTATION, AN O₂ SAT MONITOR MUST BE USED.




WARNING! WHEN USING OXYGEN SUPPLEMENTATION MODE, THE MOVES[®] SLC[™] GAS SAMPLING PORT MUST BE CONNECTED TO THE OXYGEN DELIVERY CIRCUIT (E.G., O₂ MASK SAMPLE PORT) AND AN OXYGEN SHUT-OFF DEVICE, SUCH AS BPR'S FIRESAFE[™] CANNULA VALVE, SHOULD BE USED IN THE OXYGEN SUPPLY LINE.


The concentrator provides up to 87% O₂, at a nominal flow rate of 2.75 LPM, directly to the patient during O₂ Supplement mode. Follow these instructions to supply supplementary oxygen to the patient.

1. On the Setup screen, set the system to run in O₂ Supplement mode.

NOTE: In O₂ Supplement mode, the Vent O₂ setting (under CONCENTRATOR) is grayed out, indicating that it (and its value) does not apply in O₂ Supplement mode. As mentioned above, the concentrator provides up to 87% O₂ directly to the patient.



2. Locate the O₂ Outlet port located on the LEFT side of the MOVES[®] SLC[™] unit below the large, black air intake.




3. Open the protective cover of the O₂ Outlet port (by pulling it forward) to access the O₂ Outlet.



4. Attach the free end of the tube attached to a nasal cannula (pair of nasal prongs) to the O₂ Outlet.





5. Attach the nasal cannula at the other end of the tube to the patient by inserting the ears of the prongs into both nostrils and securing the tubing behind the ears and gently securing under the chin.




9.8 Installing Suction Accessories

1. Attach the suction canister holder to the MOVES[®] SLC[™] unit by depressing the clips and sliding them along the TOP rail into the desired position. Release the clips and ensure that they lock into place.



2. The anchor will lock into place in the space **between** any two of the holes in the row (as shown in the illustration at the right).



3. Insert the suction canister into the holder as shown. Make sure that the large red plug indicated at right is firmly in place.



4. Connect the suction hose to the 'SUCTION INLET' port of the suction canister.

NOTE: The short suction hose connects the canister to the MOVES® SLC™ unit.



5. Open the protective cover of the 'SUCTION' port located on the RIGHT side of the MOVES[®] SLC[™] unit. Connect the other end of the suction hose to the 'SUCTION' port.

Pull protective cover forward to open.




6. Connect the other end of the suction hose to the 'SUCTION' port.



7. Connect the patient suction hose (long hose) to 'PATIENT PORT' on the suction canister.



8. Connect the patient suction wand to the other end of the patient suction hose.



CAUTION! THE SUCTION WAND AVAILABLE FROM THORNHILL RESEARCH IS DESIGNED TO MINIMIZE TOTAL OCCLUSION. ONLY THIS WAND OR EQUIVALENT SHOULD BE USED WITH THE MOVES[®] SLC[™] SYSTEM.



WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATE MEANS OF SUCTION IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE OR SERIOUS OCCLUSION IN THE SUCTION SYSTEM.



WARNING! THE MOVES[®] SLC[™] OXYGEN CONCENTRATOR DOES NOT FUNCTION WHILE SUCTION IS ON. AN ALTERNATIVE MEANS OF SUPPLYING O₂ WILL BE NECESSARY IF A HIGH PERCENTAGE OF O₂ IS CRITICAL.

9.9 Preparing MOVES® SLC™ for Activation

9.9.1 Checking Battery Charge

There is a Battery Condition Indicator on the front of each battery.



Figure 9-5: Battery with Condition Indicator Shown

Pressing and holding the button on the Battery Condition Indicator shows the battery condition. As battery power decreases, the illuminated LEDs extinguish from furthest left (green – meaning the battery power is at a high level), through middle indicators (orange – meaning that the battery power is at an intermediate stage), to furthest right (red – meaning that the battery power is very low or exhausted).




Figure 9-6: Condition Indicator Showing High Level Charge

When the power supply / battery charger is connected to a live AC power supply and to MOVES[®] SLC[™], any installed batteries will charge.



NOTE: The battery charge level may not appear to increase for approximately 3 hours. This is normal for the battery's initial charge and after extended periods without use.



NOTE: It is possible for batteries to be exhausted beyond normal recovery. If pushing a battery's status button does not cause any LEDs to be illuminated, place the battery in a MOVES[®] SLC[™] to charge for up to 12 hours. If this does not recover the battery's charge, the battery is non-functional and should be returned to the manufacturer.



WARNING! TO REDUCE RISK TO THE PATIENT IN THE EVENT OF A POWER FAILURE, WHEN RUNNING ON EXTERNAL POWER, ONE CHARGED BATTERY SHOULD BE PRESENT IN THE MOVES[®] SLC[™] UNIT. IN ADDITION, ALWAYS CARRY ALTERNATE MEANS OF VENTILATING, SUCTIONING, AND OXYGENATING THE PATIENT.



CAUTION! A MOVES[®] SLC[™] UNIT SHOULD NEVER BE PUT INTO TRANSPORT SERVICE WITH LESS THAN A 95% CHARGE IN BOTH BATTERIES.



WARNING! BATTERY TIME SHOWN REMAINING IS APPROXIMATE AND HIGHLY DEPENDENT ON OPERATING CONDITIONS! PUT SAFETY FIRST – ALWAYS CARRY SPARE BATTERIES!



CAUTION! NEVER CHARGE BATTERIES IN AMBIENT TEMPERATURES BELOW 32°F (0°C) OR ABOVE 104°F (40°C).

9.9.2 Inspecting the Batteries


Inspect the battery for physical damage such as cracks, holes, and leaks.



WARNING! DO NOT USE OR CHARGE A DAMAGED BATTERY!

9.9.3 Installing the Batteries

1. Lift the battery compartment latch to open it.



2. Open the battery compartment door.



CAUTION! THERE ARE BATTERY COMPARTMENT DOORS ON EITHER SIDE OF THE MOVES® SLC™. THEREFORE, A BATTERY CAN BE INSERTED FROM EITHER SIDE. HOWEVER, THE BATTERY CONNECTIONS MUST BE PROPERLY ORIENTED WITH THE RECEIVING CONNECTIONS INSIDE THE BATTERY COMPARTMENT (ON THE DOORS). CONNECTIONS ARE NOT IDENTICAL ON EITHER END OF THE BATTERY.



3. Push the battery into the battery compartment until it engages. You should not have to force it. If you feel resistance, the battery is probably not properly oriented. Pull the battery out, then check the location of the connections and reorient the battery if necessary.

NOTE: The battery is keyed to be inserted in only one orientation. If you are experiencing trouble inserting the battery, check that "key" on the top of the battery matches the slots at the top of the battery cavity.




4. Lift up and close the battery compartment door until it clicks into place.



5. Push the battery release lever in until it fully locks into place.

NOTE: Also make sure that the battery latch on the opposite side is closed completely.



6. Repeat procedures 1–5 for battery # 2.



NOTE: It is recommended that both batteries be installed during all operational times. If one battery requires a charge, it can be recharged simultaneously during normal operation when the MOVES[®] SLC[™] is connected to the power supply / charger.



NOTE: The MOVES[®] SLC[™] can operate from one battery to allow uninterrupted operation during battery replacement (since the system uses only one battery at a time). When a battery reaches an exhausted state, remove the battery and replace with a charged battery.



WARNING! WHEN MOVES[®] SLC[™] IS NOT IN OPERATION, BATTERIES SHOULD BE UNLATCHED OR REMOVED FROM THE UNIT AND STORED IN A DRY AREA AT ROOM TEMPERATURE. LEAVING BATTERIES INSTALLED IN A NON-OPERATIONAL UNIT MAY CAUSE THEM TO DRAIN TO AN UN-RECHARGEABLE LEVEL.

9.9.4 Preparing the Power Supply / Battery Charger

1. The power supply / battery charger is shown to the right. Note, there are two cords:
 - A light grey cord - this cord delivers power from a wall socket, or line supply (like a generator), to the power supply / battery charger.




CAUTION! THIS CORD COMES WITH A GROUNDING PRONG ON THE MALE CONNECTOR. TO REDUCE THE RISK OF ELECTRICAL SHOCK, THIS PRONG SHOULD NEVER BE REMOVED OR COMPROMISED.



CAUTION! THIS CORD IS A SPECIAL MEDICAL-GRADE POWER CORD AND SHOULD NOT BE REPLACED WITH A NON-MOVES[®] SLC[™] SUPPLIED PART.

- A black cord permanently affixed at one end to the power supply/charger, and with a special nine-pin female connector on the other end. This cord delivers power from the power supply / battery charger to the MOVES[®] SLC[™] unit. The female connector that attaches to the MOVES[®] SLC[™] unit is keyed to ensure proper connection.
2. Insert the light grey power cord into the front of the power supply / battery charger.




9.9.5 Connecting MOVES® SLC™ to the Power Supply / Battery Charger


1. Remove the protective cap from the receptacle of the REAR (battery rack) panel of the MOVES® SLC™ by turning it counter-clockwise .



2. Insert the nine-pin connector into the receptacle. Note the chevrons which indicate the top of the connector. Make sure the chevrons face straight up when inserting the connector, which is keyed to match the connector on MOVES® SLC™.




3. Rotate the locking collar clockwise to secure the connector.



9.9.6 Connecting AC Power

1. Insert the power supply / battery charger AC connector into a wall socket or line supply.
2. Verify the connection by checking for an illuminated PWR LED on the power supply / battery charger.
3. Note the CAUTION notice attached to the cord:
4. GROUNDING RELIABILITY CAN ONLY BE ACHIEVED WHEN THE EQUIPMENT IS CONNECTED TO AN EQUIVALENT RECEPTACLE MARKED HOSPITAL ONLY OR HOSPITAL GRADE.

NOTE: The photo at right shows a North American power connection. However, interchangeable power cords suitable for other regions are also available.



THIS PAGE DELIBERATELY LEFT BLANK

10.0 Startup



CAUTION! THE SURFACE OF THE MOVES[®] SLC[™] SYSTEM CAN BECOME HOT, ESPECIALLY IF IT IS BEING OPERATED IN DIRECT SUNLIGHT. CARE SHOULD BE TAKEN WHEN TOUCHING OR CONTACTING THE SURFACE OF THE MOVES[®] SLC[™] SYSTEM.

10.1 Adjusting the Screen Display Orientation

The display screen can be locked or adjusted to face either the RIGHT or LEFT side of the MOVES[®] SLC[™] to facilitate operation from either side of the unit.

1. In the picture at the right, the screen is locked. Turning the catch indicated 90 degrees to either the left or the right will unlock the screen and allow it to be raised into a working position.




2. In the picture at the right, the screen is unlocked. With the arrow shown on the catch pointing in this direction, the screen will lift and face toward the RIGHT side of the MOVES[®] SLC[™]. This is shown in the following photograph.



3. In the picture at the right, the screen is facing toward the RIGHT side of the MOVES[®] SLC[™].



4. In the picture at the right, the screen is unlocked. With the arrow shown on the catch pointing in this direction, the screen will lift and face toward the LEFT side of the MOVES[®] SLC[™]. This is shown in the following photograph.



5. In the picture at the right, the screen is facing toward the LEFT side of the MOVES[®] SLC[™].



NOTE: The on-screen display can be flipped by pressing any of the panel buttons on the side of the screen you would like the orientation to switch to.



Figure 10-1: MOVES[®] SLC[™] Panel Buttons

10.2 User Interface (UI) Controls and Functions



Figure 10-2: MOVES[®] SLC[™] Main Screen

The following table describes the physical appearance, function and the effect or use of each of the components used to operate the MOVES[®] SLC[™] unit.

Table 20: User Interface Controls and Functions

Button / Icon	Function
See Figure above.	Main Liquid Crystal Display (LCD) Screen This screen is the main component of the user interface display. The operator interacts with the display to operate the MOVES [®] SLC [™] unit, getting information from the screen while selecting modes, options and numerical values for settings using the panel buttons.
	Power Control Button This green and white push button (which is next to the power connection on the end of the MOVES [®] SLC [™] unit) is used to activate or deactivate the unit. When the button is initially pressed, the system activates, the UI screen illuminates and an audible indicator sounds. When pressed and held for five (5) seconds when the unit is running, the button turns the unit OFF. As well, a message appears on the UI screen indicating this functionality.
	Screen Dim Button This small black button at the far end of the buttons panel is used to cycle the screen display and visual alarm indicator lights through six (6) diminishing levels of brightness.

Button / Icon	Function
	Alarm Audio Pause Button The Alarm Audio Pause button can be used either to silence ALL alarms permanently or for a temporary period (120 seconds). Pressing the Alarm Audio Pause button for one (1) second and then releasing it will silence ALL alarms temporarily. Pressing the Alarm Audio Pause button for at least three (3) seconds and then releasing it will silence ALL alarms permanently (if allowed, see note below). Pressing the Alarm Audio Pause button for one (1) second again then releasing it will re-enable ALL alarm audio. When alarms are paused temporarily, a dashed “X” is shown through the bell image (see image at left). When alarms are paused permanently, a solid “X” is shown through the bell image. If an alarm is turned off on the Alarm On/Off screen, its alarm will not sound. NOTE: Silencing all alarms permanently can be disallowed using the Admin screen. For more information, consult the section of this manual that deals with the Admin screen.
	NIBP Control Button The NIBP Control button is used to start or abort a Non-Invasive Blood Pressure (NIBP) reading. Prior to taking an NIBP reading, the NIBP cuff is fitted to the patient with the free end of the sensor cord attached to the patient-connector panel. The NIBP Control button is pressed to begin taking a reading, usually manually. The data taken from the sensor appears in numerical form on the Monitor screen (if currently shown). Pressing the button a second time aborts the reading in progress. The button icon will be shown on the Status Bar if an NIBP reading is in progress (i.e., cuff is pressurized). Otherwise, the icon will not be shown.  NOTE: Using the Setup screen, you can set the NIBP control to manual, or to automatically update at set intervals: every 1–5 minutes, 10 minutes, 15 minutes or STAT. STAT, an abbreviation for the Latin word Statim which means “immediately”, will set MOVES[®] SLC[™] to take readings as often as is considered safe (i.e., a minimum time of 30 seconds is enforced between the end of one NIBP reading and the beginning of the next reading) for a maximum of 15 minutes. After 15 minutes in STAT mode the MOVES[®] SLC[™] will automatically transition the update interval to 5 minutes. Pressing the NIBP Control Button manually to take a reading, or to abort a reading, temporarily overrides any automatic setting; however, the timer is reset to the time of the next automatic reading.

Button / Icon	Function
	Suction Control Button The Suction Control button is used to activate and deactivate the suction pump. All suction accessories must be connected and ready to use before activating. When the button is pressed, the suction pump activates. When the button is pressed again, the suction pump deactivates. While suction is ON, the Suction icon and suction values (in mmHg) are displayed in the Status Bar  NOTE: When in Suction mode, the concentrator/suction pump operates at a default (maximum) intake pressure level of 325 mmHg. This value can be changed (lowered) on the Setup screen.  WARNING! THE O₂ CONCENTRATOR DOES NOT PRODUCE OXYGEN WHEN THE SUCTION FEATURE IS IN USE. IF SUCTION USE IS PROLONGED, THE PATIENT WILL BECOME HYPOXIC UNLESS EXTERNAL O₂ IS SUPPLIED.
	Screen Button Press the Screen button to move to the next screen in sequence. If there are outstanding queries or parameters that need to be satisfied, the functionality of the Screen button will be inhibited until these queries or parameters are satisfied. NOTE: This button can also function as an <i>Accelerator</i>. Pressing and holding the button will take you from any screen to the <i>Monitoring</i> screen. If you are already on the <i>Monitoring</i> screen, pressing and holding the button will take you to the <i>Setup</i> screen.
	Cancel Button Press this button to reject the current value that is being changed and return to the previous value.
	Check Button Press the Check button to start or end editing a currently selected item, or change the currently selected item. The Check button is also used to acknowledge messages.
	Next Button Selects the next item in a group or increases a numeric value.
	Previous Button Selects the previous item in a group or decreases a numeric value.

10.3 System Visual Indicators (Alarms / System Status)

There are four (4) System Visual Indicators, one at each top corner of the MOVES[®] SLC[™]. These indicators display three (3) different colors: green, red and yellow. The severity of the highest active alarm is shown. If an alarm is turned off, the

display defaults to the severity of the next active alarm. These System Visual Indicators provide the user with the ability to see alarms, or system alarm status, from a greater distance and angle than would be possible with only screen display. The brightness of the System Visual Indicators is controlled by the front-panel brightness control.

System Visual Indicator states are explained in the following table.

Table 21: System Visual Indicator States and Explanations

System Visual Indicator State	Indicates
Off	System off or in startup testing
Solid Green	No alarms active
Solid Yellow	Low Priority Alarm active
Flashing Yellow	Medium Priority Alarm active
Flashing Red	High Priority Alarm active



NOTE: Even when the audio of active alarms is temporarily silenced using the Alarm Audio Pause button, the System Visual Indicators continue to display the corresponding visual alarm signal. If all active alarms are turned OFF (using the Alarm ON/OFF screen), or if all alarm conditions become satisfied, the indicator returns to a steady green illumination.

10.4 Changing Settings and Data Views

MOVES[®] SLC[™] has settings that the user can modify, for example, O₂, Maximum Airway Pressure, etc. As well, there are customizable interactions to review data. All settings/data that can be manipulated are identified the same way, selected the same way, and are changed/viewed the same way. All modifiable settings/data views are shown on the screen as either a *screen button* or a *selectable area* with multiple elements displayed. These modifiable settings/data views can be selected and made active using three (3) of the *physical* buttons located below the screen; specifically, the Next and Previous (arrow) buttons in conjunction with the Check button.




NOTE: The last selected screen element is remembered for each screen. When the user returns to a screen, the most recently selected element on that screen is reselected.

10.4.1 Display of Settings and Views

Screen buttons and areas are displayed in one of three states Selectable, Selected, and Active.

Table 22: Screen Buttons and Descriptions

Screen Button / Area State	Description	Example
Selectable	Shown with a double-line white border. Identifies the button or area as something that can be selected for modification.	
Selected	Shown the same as Selectable but with yellow lines forming the double-line border. Identifies the currently selected button or area.	

Screen Button / Area State	Description	Example
Active	Shown with a single, dashed line. Identifies the currently selected button or area as active and that changing of a setting is in progress.	

10.4.2 Modifying a Setting

To modify a setting:

1. Select the desired button or area by pressing the Next or Previous button successively until the desired selection is reached.
2. Press the Check button to initiate modification of the setting.
3. Use the Next or Previous button to select/modify the setting.
4. Press the Check button again once to accept the modified setting.



NOTE: The last selected screen element is remembered for each screen. When the user returns to a screen, the most recently selected element on that screen is reselected.



NOTE: In the case of the Alarm On/Off settings:

- Pressing the Check button toggles between the two states (i.e., on and off).
- To finalize and exit to the top level navigation, press the Cancel (X) button.

10.5 Confirming Display Validity

Note the following in regard to the display of valid information on the MOVES[®] SLC[™] UI Screen:

The display should NOT show any elements that are partially written.

- The display should ALWAYS show information while MOVES[®] SLC[™] is on.



CAUTION! IF ANY OF THE ABOVE CONDITIONS ARE NOT SATISFIED, THE MOVES[®] SLC[™] UNIT SHOULD BE DELIVERED TO AUTHORIZED PERSONNEL FOR SERVICING.

10.6 Startup Sequence



NOTE: There should be **NO external oxygen** connected to the MOVES[®] SLC[™] system during startup tests. Having O₂ connected creates a flow in the inspiratory limb which causes the open-circuit test to fail.

1. Locate the (green) Power Control button on the rear of the MOVES[®] SLC[™] unit (the end of the unit with the hydrocarbon filter and power-cord connection).



Figure 10-3: Power Control Button

2. Press the Power Control button to activate the unit – an audible alert is sounded and the user interface displays the System Test screen.



NOTE: If the system has been run within the last 30 minutes, a prompt screen (New Patient Screen) will appear asking the user if the system should be configured for a new patient (i.e., begin with default values or restore previously used ones.)




Figure 10-4: New Patient Screen




Figure 10-5: System Test Screen



NOTE: All active alarms in the alarm queue should be reviewed for patient monitoring fault alarms immediately after completing the startup test procedure and before connecting to the patient.



WARNING! ONE CHARGED BATTERY MUST BE PRESENT IN THE MOVES[®] SLC[™] UNIT AT ALL TIMES, EVEN WHEN RUNNING ON EXTERNAL POWER. THIS REDUCES THE RISK TO THE PATIENT IN THE EVENT OF A POWER FAILURE.



CAUTION! IF INTENDING TO RUN ON BATTERIES, ENSURE THAT THERE IS SUFFICIENT POWER FOR THE LENGTH OF TIME REQUIRED, OR REPLACE THE BATTERIES.



NOTE: If a battery appears to be installed but is shown as missing on the display, inspect the battery and replace or reinstall. Also ensure that the battery door and latch are correct on both sides of the machine.

3. If not all initial tests are passed, you will be offered three choices:

- Skip remaining
- Prep done
- Skip this test

4. Use the Check button  to confirm a selection. Use the Next button  to move forward to the next item.

Use the Previous button  to move back to an item in the list.



WARNING! DO NOT CONNECT ANY SENSORS, MONITORS, OR THE BREATHING CIRCUIT TO THE PATIENT WHILE PERFORMING SYSTEM TESTS! DOING SO COULD ENDANGER THE PATIENT!



CAUTION! ONLY SKIP A TEST OR TESTS IF A SYSTEM TEST HAS PREVIOUSLY BEEN PERFORMED AND PASSED, AND THE SYSTEM HAS JUST BEEN RESTARTED DUE TO POWER FAILURE (E.G., LOSS OF BATTERY OR EXTERNAL POWER).



WARNING! DO NOT OPERATE THE MOVES[®] SLC[™] SYSTEM UNTIL ALL SYSTEM TEST FAILURES HAVE BEEN RESOLVED, AND ALL TESTS HAVE BEEN REPEATED AND PASSED.

5. Once prerequisites have been satisfied, the user must select Continue to move to the Setup screen.




Figure 10-6: Setup Screen

6. The initial system mode will be Monitor Only. If a different system mode is desired, first configure the appropriate settings and then choose the new system mode (i.e., Ventilate or O2 Supplement).



NOTE: If an ABP/CVP/ICP transducer is to be connected to the system, it will be necessary to reset or 'zero' the transducer. For more information, see the section of this manual that deals with Zeroing the Pressure in the Transducer.



WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATIVE MEANS OF VENTILATION CAPABLE OF SUPPLEMENTING A HIGH CONCENTRATION OF O₂ IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE, OR SERIOUS OCCLUSION IN THE VENTILATOR SYSTEM.



WARNING! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATE MEANS OF SUPPLYING A HIGH CONCENTRATION OF O₂ IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE, OR SERIOUS OCCLUSION IN THE CONCENTRATOR CIRCUIT.



WARNING! THE OPERATOR SHOULD USE AN ALTERNATIVE MEANS OF VENTILATION UPON EXPERIENCING A PROLONGED APNEA ALARM.

THIS PAGE DELIBERATELY LEFT BLANK.

11.0 MOVES[®] SLC[™] User Screens

MOVES[®] SLC[™] has multiple user screens.

11.1 Status Bar

At the top of all screens is the Status Bar. The Status Bar is used to display system status and alarms and the name of the current screen. The Status Bar contains two buttons – the Screen Select button and the Alarm View button. The Status bar shows the system power status (including battery charge and external power status), the number of alarms turned off, the number of alarm limits changed from their default value, alarm audio status, NIBP status, suction status, alarm status, and the system state.






Figure 11-1: Status Bar

The following table explains the Status Bar items.

Table 23: MOVES[®] SLC[™] Status Bar Items and Descriptions

Label/Name	Items on the Status Bar	Description
Plug icon with battery icons 	Battery charge for each battery, and charge status. Also indicates whether external power is connected.	- Graphic showing charge state of each battery and whether external power is present. The level of a battery's charge is shown in red when battery life is low, orange when partially depleted and green when high. The system power status area is highlighted when one or more alarms is associated with the system power status. (NOTE: See the Battery Status Icon Table that follows for more information.) - If the state of a battery cannot be determined, its outline is shown with the fault (X) icon. - If the entire power state cannot be determined, neither battery is shown, and a fault (X) icon is displayed only.
Alarms OFF icon 	Icon shows that some (or all) alarms are off. Number below indicates how many alarms are off.	- The Alarm OFF icon is shown if any alarms have been turned off by the user on the Alarm ON/OFF screen. - The number of alarms 'XX OFF' is shown in text below the icon. - No icon or text is shown if all alarms are ON.
Number of limits changed icon 	Icon shows that some (or all) alarm limits have changed from their default values. Number below indicates how many alarm limits have been changed.	- The Alarm Limits Changed icon is shown if any alarm limits have been changed from their default values by the user on the Alarm Limits screen. - The text below the icon indicates the number of alarm limits changed: '1 limit chgd' shown if only 1 limit is not set to its default value. 'X lims chgd' shown if 2-9 limits are not set to their default values. 'XX lim chgd' shown if 10 or more limits are not set to their default values. - No icon or text is shown if all alarms limits are set to their default values.

Label/Name	Items on the Status Bar	Description
Audio PAUSE icon 	Icon shows that alarm audio is PAUSED. Numbers display how many minutes and/or seconds remain for audio pause. OR	- The Audio PAUSE icon is shown if audio pause is active. This function is controlled by quickly pressing and releasing the Audio PAUSE button on the front panel. PAUSE time is 120 seconds max. - The time remaining for audio pause is shown in text below the icon as 'm:ss', for example, 1:23 (1 minute 23 seconds); or, 0:03 (3 seconds). - The Audio OFF icon is shown if audio OFF is active. This function is controlled by pressing and holding the Audio PAUSE button on the front panel for three (3) seconds. - Audio PAUSE can only be active if Audio OFF is not active. - No icon or text is shown if Audio PAUSE / OFF is not active.
Audio OFF icon 	Icon shows that alarm audio is OFF. Text describes that audio is OFF.	
NIBP icon 	Icon shows that NIBP is active. Text below indicates 'ON' status.	- The NIBP icon is shown if NIBP is obtaining a reading, which is controlled by the NIBP button on the front panel and/or by the auto NIBP setting on the Setup Screen. - Text below the ICON is simply ON. - Note that the cuff press is displayed with the NIBP reading on the main screen. - No icon or text is shown if NIBP is not active. - The NIBP fault icon is shown if the status of NIBP cannot be determined or is known to be in a fault condition.
Suction Icon 	Icon shows that Suction is active. Text below icon shows the suction pressure.	- The numerical value below the Suction icon indicates the current amount of suction present in mmHg. - The Suction fault icon (X) is shown if the status of Suction cannot be determined or is known to be in a fault condition.
Screen Select Button	Button indicates which screen is active. 	When the Screen Select button is selected (double yellow lines), pressing the Check button makes it active (dashed single yellow line). Then pressing the Next and Previous buttons will cycle through the list of screens. Pressing the Check button a second time will confirm the screen choice and cause the chosen screen to be displayed.
System State	Text indicates which mode is active. 	Text describes the system mode: - Monitor Only - O2 Supplement - When the system mode is Ventilate mode, the system state text will describe the specific ventilation configuration (e.g., VC-SIMV+PS, which indicates the system is in Volume Controlled Synchronized Intermittent Mandatory Ventilation with spontaneous breaths receiving Pressure Support).

Label/Name	Items on the Status Bar	Description
Alarm View Button	Alarm Queue	3 line alarm message queue top right of screen. The alarm queue shows the highest priority alarm.  To view additional alarms (if existing), use the Next or Previous buttons to select the Alarm Queue (double yellow lines). Then, press the Check button to make it active (dashed yellow line). Then, pressing the Next and Previous buttons will cycle through the alarms. Pressing the Cancel button will exit the Alarm Queue and return to top level navigation.

11.1.1 Battery Status Icon

The Battery status icon shows the charge of the two system batteries individually; it also indicates if external power is connected. The batteries are drawn beside each other and labeled A and B. Charge is shown by filling in each chamber accordingly, indicated by 5 charge segments (the same number of segments as the physical battery LED status on the battery). When an alarm associated with system power is active, the area behind the battery status icons is highlighted with a red background if highest associated alarm is High priority. Otherwise, a yellow highlight is used for Medium/Low alarm priority. When the system is connected to external power, a plug icon is drawn above the battery icons.

Note the following:

- Each battery is individually hot swappable.
- If a battery has no charge, no segments will be shown.
- If a battery is not present, no segments and no battery outline will be shown.
- If MOVES[®] SLC[™] is evaluating a battery's status, a question mark will be drawn over that battery's status area.
- If a battery's status cannot be determined, a fault icon will be drawn over the battery's status area.
- If a battery is charging, a 'lightning bolt' will be drawn over the battery icon and its segments.
- If the power status is not known, a large fault icon will be displayed.

The color of the segments for a battery depends on the number of segments shown. If one segment is shown, the color is red; if two or three, the color is orange; if more than three segments are shown, the color is green. This matches the color of the top most LED on the battery charge indicator on the battery itself.

Table 24: Battery Status Icon Table

Battery Status Icon Table		
Battery A Full Charge Battery B Half Charge Running on Battery No associated alarm(s)	Battery A Full Charge Battery B Full Charge Running on Battery No associated alarm(s)	Battery A Unknown Battery B Near Empty Running on Battery High priority alarm
Battery A Missing Battery B Half Charge Running on Battery High priority alarm	Battery A Missing Battery B Missing Running on External Power High priority alarm	Battery A Unknown Battery B Half Charge, Charging Running on External Power High priority alarm
Battery A Full Charge Battery B Full Charge Running on External Power No associated alarm(s)	Power Status Unknown Power Status Unknown Alarm status ignored	Battery A Near Empty Battery B Near Full Running on Battery Medium/Low priority alarm

11.2 System Test Screen



NOTE: The startup tests perform functional testing on components of the system that require operator intervention to validate – other parts of the system are continually tested and will show failure status in the alarm queue.



NOTE: The operator should wait until the O₂/CO₂ sensor has warmed up and O₂/CO₂ values are displayed before connecting a patient.



NOTE: All active alarms in the alarm queue should be reviewed for patient monitoring fault alarms immediately after completing the startup test procedure and before connecting to the patient.

Upon starting MOVES[®] SLC[™] the user, in most cases, will be presented with the System Test screen.




Figure 11-2: System Test Screen

The System Test screen will indicate to the user which tests have passed and which need preparation before they can be conducted. A description of the preparation required is given at the right of the screen. Once preparation has been done,

the user presses the Next button  to select *Prep done* and then the Check button  to initiate the test. The user is also presented with the options of skipping a particular test or all remaining tests.



NOTE: In order to save the user time, certain tests can be run concurrently. For example, while the Open Y-Piece test is still running, the screen will advance to the Battery Verify test, which can be run concurrently. After completing that test, the screen will advance to the Sample Line test, which can also be run concurrently.

The System Test screen will not immediately appear upon startup if:

1. The system has been shut down for a period of time less than or equal to 3 minutes. In this case, the system will auto resume to the point it was at prior to shutdown. (**NOTE:** A temporary loss of power is assumed here.)
2. The system has been shut down for a period of time greater than 3 minutes but less than 30 minutes. In this case, the system will query the user as to whether the patient is new or continuing (and the following screen will be displayed).




Figure 11-3: Configure for a New Patient Screen

Selecting YES will return the system settings to default values (see *Section 16.1 System Default Settings* on page 207 for details). Selecting NO will keep the last system settings used. **NOTE:** NO should only be selected if the operator is aware of the last system settings configuration or the operator should review the settings on the Setup, Alarm Limits, Alarm ON/OFF and Advanced screens.

11.3 Setup Screen

11.3.1 Overview

The Setup Screen is used to view and modify the primary operational settings of MOVES[®] SLC[™]. The color of a field is dependent on the mode selected. If a field is white, it is currently in use; otherwise, the field is gray. This lets the user see all the primary settings, including those that are not used by the current system mode. Current settings are underlined (and are indicated even if an operational mode is not in use).

11.3.2 Changing Settings



NOTE: Settings can be changed even in a mode where the setting is not active.

To change system settings:

1. Use the Next and /or Previous buttons to navigate to the System Mode area (it becomes surrounded by a double yellow line).
2. Press the Check button (the System Mode area becomes surrounded by a dashed yellow line).
3. Use the Next and /or Previous buttons to navigate to the setting you want to change (it becomes surrounded by a single yellow line).
4. Press the Check button (the setting becomes surrounded by a dashed yellow line).
5. Use the Next and /or Previous buttons to select a new setting.
6. Press the Check button to confirm the new setting (the setting becomes surrounded by a single yellow line).
7. Press the Cancel button to exit the System Mode area.

11.3.3 Setup Screen – Ventilate Mode



NOTE: Because there are four (4) Vent Modes (IMV, A/C, SIMV, PSV), certain settings will be grayed out under each. However, these settings can still be changed even though the mode they apply to is not currently selected.

SYSTEM MODE: Ventilate O2 Supplement Monitor Only

VENTILATOR

Vent Mode: IMV A/C SIMV PSV

Breath Sensitivity: Normal Low

Control: Volume Pressure

Tidal Volume: 500 mL

PC: 20 cmH2O

PEEP: 0-3 cmH2O

RR: 10 B/M

I:E Ratio: 1:2

Inspire Time: 1.0 s

Pressure Support: disabled

CONCENTRATOR

Vent O2: 40%

OTHER

Line Filter: 50Hz 60Hz

ECG Leads: 3 5 12

ECG Rate: slow fast

ECG Range: 3mV 6mV 12mV

NIBP Mode: MANUAL

Suction: 325 mmHg

ADVANCED...

Figure 11-4: Setup Screen – Ventilate Mode

11.3.4 Setup Screen – O2 Supplement

SETUP **No Alarms**

O2 Supplement

SYSTEM MODE: Ventilate O2 Supplement Monitor Only

VENTILATOR

Vent Mode: IMV A/C SIMV PSV

Breath Sensitivity: Normal Low

Control: Volume Pressure

Tidal Volume: 500 mL

PC: 20 cmH2O

PEEP: 0-3 cmH2O

RR: 10 B/M

I:E Ratio: 1:2

Inspire Time: 1.0 s

Pressure Support: disabled

CONCENTRATOR

Vent O2: 40%

OTHER

Line Filter: 50Hz 60Hz

ECG Leads: 3 5 12

ECG Rate: slow fast

ECG Range: 3mV 6mV 12mV

NIBP Mode: MANUAL

Suction: 325 mmHg

ADVANCED...

Figure 11-5: Setup Screen – O2 Supplement

11.3.5 Setup Screen – Monitor Only

SETUP **No Alarms**

Monitor Only

SYSTEM MODE: Ventilate O2 Supplement Monitor Only

VENTILATOR

Vent Mode: IMV A/C SIMV PSV

Breath Sensitivity: Normal Low

Control: Volume Pressure

Tidal Volume: 500 mL

PC: 20 cmH2O

PEEP: 0-3 cmH2O

RR: 10 B/M

I:E Ratio: 1:2

Inspire Time: 1.0 s

Pressure Support: disabled

CONCENTRATOR

Vent O2: 40%

OTHER

Line Filter: 50Hz 60Hz

ECG Leads: 3 5 12

ECG Rate: slow fast

ECG Range: 3mV 6mV 12mV

NIBP Mode: MANUAL

Suction: 325 mmHg

ADVANCED...

Figure 11-6: Setup Screen – Monitor Only

11.3.6 Setup Screen Options

Table 25: Setup Screen Options and Descriptions

Setup Screen Options		
Fields on the Setup Screen	Used by System Mode	Description
System Mode	All	- Ventilate Mode, O2 Supplement, Monitor Only. - Default: Monitor Only - Change immediately after change confirmed
Vent Mode	Ventilate Mode	- MV, A/C, SIMV, PSV. - Default: IMV - Change at start of next inhalation after change confirmed.
Breath Sensitivity	Ventilate Modes: A/C, SIMV	- Normal, Low. - Defaults: Normal
Control	All Ventilate modes except PSV	- Pressure, Volume. - Default: Pressure - Change at start of next inhalation after change confirmed.
Tidal Volume	Ventilate Modes: IMV, A/C, SIMV when Control is set to Volume	100 to 750 ml in 25 ml intervals - Default: 500 ml - After a change is accepted, the option selected will be active after the patient's current breath.
PC	Ventilate Modes: IMV, A/C, SIMV when Control is set to Pressure	- Pressure Control 10 to 55 cmH₂O in intervals of 1 cmH₂O - Default: 20 cmH₂O - After a change is accepted, the option selected will be active after the patient's current breath. NOTE: PC (Pressure Control) = PIP (Peak Inspiratory Pressure) – PEEP (Positive End Expiratory Pressure)
PEEP	Ventilate Mode	- Positive End-Expiratory Pressure 0-3, 4 to 20 cmH₂O in intervals of 1 cmH₂O - Default: 0-3 cmH₂O - After a change is accepted, the option selected will be active after the patient's current breath. NOTE: When the PEEP value is increased by 10 cmH₂O or more some additional breath triggers may be observed. This includes starting the ventilator at a PEEP setting of 10 cmH₂O or more.

Setup Screen Options

RR	All Ventilate modes except PSV	• Respiratory Rate • 6 to 40 B/M in steps of 1 B/M • Default: 10 B/M • Change at start of next inhalation after change confirmed.
I/E ratio	Ventilate Modes: IMV, A/C	• Inspiratory/Expiratory Ratio • 1:1.0, 1:1.5, 1:2.0, 1:2.5, 1:3.0 • Default: 1:2.0 • Change at start of next inhalation after change confirmed.
Inspire Time	Ventilate Mode: SIMV	• Inspiratory Time • 0.3 to 3.0 seconds in steps of 0.1 seconds • Default: 1.0 second • Change at start of next inhalation after change confirmed.
Pressure Support	Ventilate Modes: SIMV, PSV	• Disabled, 5 to 40 cmH₂O in intervals of 1 cmH₂O • Default: Disabled • Change at start of next inhalation after change confirmed.
Vent O2	Ventilator Mode Only	• Ventilate Mode: Air, 30%, 40%, 50%, 60%, 70% or 85% • Default: 40% • Change immediate after change confirmed.
Line Filter	All	• 50 Hz, 60 Hz • Last setting remembered (i.e., not defaulted with new patient). • Change immediate after change confirmed.
ECG Leads	All	• 3, 5 or 12 • Default: 12 • Change immediate after change confirmed. • Prescribes the minimum number of lead wires required to be connected before ECG readings become available. ○ 3: Left Arm, Right Arm, and Left Leg ○ 5: Left Arm, Right Arm, Left Leg, Right Leg, and Chest • 12: Left Arm, Right Arm, Left Leg, Right Leg, and Chest 1-6

Setup Screen Options

ECG Rate	All	- Speed of ECG. ECG is displayed on graphs on Monitor screen or dedicated ECG screen. - Slow or Fast. - Default: Fast - Slow is 60 Hz data point display (10.8 mm/sec sweep speed) - Fast is 150 Hz (27 mm/sec sweep speed) - Change immediate after change confirmed.
ECG Range	All	- Maximum amplitude range of ECG. ECG is displayed on graphs. 3.0 mV, 6 mV or 12.0 mV range. - Default: 3.0 mV - Change immediately after change confirmed. NOTE: Range is centered about zero (e.g., 3.0 mV is +/- 1.5 mV).
NIBP	All	- NIBP Measurement Period - Manual, Stat, 1, 2, 3, 4, 5, 10, 15 min - Default: Manual - In <i>Manual</i> mode no automatic NIBP measurements are taken. Press the NIBP button to start a measurement. <i>Stat</i> mode acquires NIBP measurements as fast as possible while guaranteeing 30 seconds between the end of a measurement and the start of the next measurement. - Change immediately after change confirmed. NOTE: After 15 minutes of <i>Stat</i> mode, the mode will automatically switch to 5 minute mode. IEC 80601-2-30 201.105.2 states that <i>Stat</i> mode can only last for 15 minutes, and then must go to manual or long term automatic. NOTE: If NIBP measurement period is set to anything other than <i>Manual</i>, and three automatic NIBP measurements fail to obtain a reading, then the system will automatically switch the mode to <i>Manual</i> and notify the operator of this change.
Suction	All	100 to 325 mmHg in steps of 25 mmHg - Default: 325 mmHg
Advanced	All	Conveys user to the Advanced screen where SpO₂ Average Time, SpO₂ Sensitivity Mode and SpHb Averaging Window options can be set. Additionally, SpO₂ feature availability can be verified and system firmware versions reviewed.

11.4 Advanced Screen

11.4.1 Accessing the Advanced Screen

The Advanced Screen is accessed from the Setup Screen. The user navigates to the System Mode area and presses the Check button to make it active. The user then navigates to the Advanced option and presses the Check button again to proceed to the Advanced Screen.




Figure 11-7: Accessing the Advanced Screen

11.4.2 The Advanced Screen Itself

The selection order for the buttons on the Advanced screen is the following: Screen Button, Alarm View Button, Options Area.




Figure 11-8: The Advanced Screen

The Options Area is selected as on other screens, by navigating to it and pressing the Check button to make it active. When this is done, the first option will be automatically selected.




Figure 11-9: Advanced Screen – SpO2 Avg Time Selected

The following table lists the options and values available on the Advanced Screen.

Table 26: Advanced Screen Options

Advanced Screen Options		
Option	Values Available	Default Value
SpO2 Avg Time	2-4s, 4-6s, 8s, 10s, 12s, 14s, 16s	8s
SpO2 Sensitivity Mode	MAX, NORMAL, APOD	NORMAL
SpHb Avg Window	LONG, MEDIUM, SHORT	LONG
ECG EMG Filter	OFF, ON	ON

11.4.3 SpO2 Avg (Average) Time

The user-selectable SpO₂ averaging feature allows the operator to select the desired level of visibility to subtle variations in the measured value. Eight (8) second averaging is generally considered the most common averaging interval, and it recommended for most patients since it is short enough to provide visibility to subtle desaturations while also being long enough to minimize major changes in SpO₂ due to quick, transitory desaturations.

11.4.4 SpO₂ Sensitivity Mode

The three sensitivity settings allow the operator to adapt the SpO₂ measurement sensitivity to the patient's level of SpO₂ signal strength and quality at the measurement site.

MAX

MAX mode is used to obtain and display data even when the signal is very weak due to impaired perfusion (can be used, for example, during treatment or examination, i.e., when someone is with the patient). If the sensor becomes detached from the patient, this mode makes virtually no provision for measurements that are displayed erroneously.

NORMAL

NORMAL is the recommended mode for typical monitoring, such as intensive care units.

APOD (Adaptive Probe Off Detection)

APOD is the recommended mode when it is not possible to detect immediately if a sensor has become detached. This mode offers better protection against erroneous measurements being displayed, even though the sensor has become detached from the measurement point (e.g., as a result of the patient moving).

11.4.5 SpHb Average Window

This feature provides 3 options of Short, Medium, and Long averaging times. The user-selectable averaging feature allows the operator to select the desired level of visibility to subtle variations in the SpHb value. This allows an operator to fine tune SpHb responsiveness to achieve the desired level of visibility to rapid variations in SpHb values.

11.4.6 ECG Electromyography (EMG) Filter

When enabled (i.e., ON) this feature filters the ECG waveforms for EMG interference in the frequency range of 15-30 Hz.

11.4.7 Additional Information Located on Advanced Screen

Table 27: Advanced Screen Additional Info Items

Masimo Info	• Board S/N (serial number) • Board Type • DSP F/W (digital signal processor firmware version) • MCU F/W (microcontroller firmware version) • Product ID (identification)
Masimo Supported Params (parameters): The screen provides a list of all supported parameters for the <i>Pulse CO-Oximeter (Board)</i> and the currently attached sensor (<i>Current Sensor</i>). A parameter must be supported by the Pulse CO-Oximeter and by the sensor in order for the parameter measurement to be displayed. Listed below are all possible parameters: • SpO2 • PR - Pulse Rate • PI - Perfusion Index • SpCO • SpMet • SpHb • SpOC NOTE: See Section 6.4 Pulse Oximeter Technology Overview beginning on page 25 for a description of each of these parameters.	Board: • SpO2, PR and PI are always supported • The extended parameters SpCO, SpMet, SpHb and SpOC are optional features which may or may not be present depending on the requested configuration when MOVES[®] SLC[™] was manufactured. Current Sensor: • All possible parameters the currently attached sensor can measure will be listed, even if the Pulse CO-Oximeter does not support the parameter. • If no sensor is presently connected to the Pulse CO-Oximeter, then << not available >> is shown
System Firmware Info (version date)	• Monitor (Patient Monitor) • Power (Power Manager) • Ventilator • Concentrator • UI (User Interface) • UI Base • Battery A • Battery B



NOTE: If batteries are not inserted the << info not available >> notification is displayed.

11.5 Main Screen

The Main Screen is used to display patient status and monitored values. It is the primary screen of MOVES[®] SLC[™]. The Main Screen is for status only. There are no functional settings to change.



NOTE: The values displayed on the Monitoring Screen are those that are actually measured from the patient and not how the system is configured (e.g., The Monitoring Screen displays the actual volume, pressure and respiratory rate measured from the ventilator)

The available buttons are for the following:

- Choosing which chart to view
- Choosing the source of heart rate to view
- Zeroing the invasive pressures channels
- Labeling the invasive pressure channels

Any value that has an active alarm associated with it is drawn in reverse color (non black background).



Figure 11-10: Screen Items Drawn in Reverse

For all numeric or text display items, if no valid data is available, dashes are displayed.



Figure 11-11: Screen Items Displaying Dashes

For all numeric or text display items, if a fault prevents data display the fault icon (X) is displayed.




Figure 11-12: Screen Items Displaying the Fault Icon

For all display items, a trusted number is always shown or represented; otherwise, for values that are not possible – and therefore not trusted – dashes are displayed and the systems acts as if the data is not available. An alarm is usually generated for this untrusted number.




Figure 11-13: Main Screen with Dashes & Alarms

11.5.1 Selection Order for Options

The selection order for Main Screen buttons is the following:

1. Screen Button
2. Alarm View Button
3. HR Source
4. Chart 1 – Data Type
 - If Trends, followed by Data Source, Data Type & Time
 - If ECG, followed by Data Source
 - If IP, followed by Data Source
 - All others followed by Chart 2
5. Chart 2 – Data Type
 - If Trends, followed by Data Source, Data Type & Time
 - If ECG, followed by Data Source
 - If IP, followed by Data Source
 - All others followed by IP1
6. IP 1
7. IP 2
8. IP 3




Figure 11-14: IP & Trends Graph Options & Sub-Options

11.5.2 Main Screen Items

Table 28: Main Screen Items and Descriptions

Label	Items on the Main Screen	Description
Vte	Expired Tidal volume	- Numeric 10– 2500 ml in steps of 10 ml, or a series of dashes (- -) for unknown. - Only displayed when ventilating; otherwise, a series of dashes (- -) is shown. ">2.5 L" is shown if 2510 ml or higher.
Vti	Inspired Tidal volume	- Vti is displayed in the location of Vte, alternating with Vte display every 2 seconds. This is only shown if: "High expired tidal volume" alarm is active <u>OR</u> "Leak detected" alarm is active. - Numeric 10–2500 ml in steps of 10 ml, or a series of dashes (- -) for unknown. - Only displayed when ventilating; otherwise, series of dashes (- -) is shown. ">2.5 L" is shown if 2510 ml or higher.
RR	Respiratory Rate, calculated from the CO ₂ monitor	- Numeric 0–99 B/M in steps of 1, or a series of dashes (- -) for unknown value. ">99" is shown if 100 B/M or higher.
PIP	Peak Inspiratory (Airway) Pressure	- Numeric 0–100 cmH₂O in steps of 1, or a series of dashes (- -) for unknown / out of range. - Only displayed when ventilating, otherwise, or a series of dashes (- -) is shown.
PEEP	Positive End Expiratory (Airway) Pressure	Same as PIP
FiO₂	Fraction of Inspired Oxygen (O ₂)	- Numeric 0–100 % in steps of 1 or (- -) for unknown value. - If >100 and ≤105, shows 100, if >105 shows (- -). - If unavailable due to calibration / warming-up shows "CAL".
ConcFlow	Concentrator Flow	- Volumetric flow of oxygen from the concentrator into the ventilator or O₂ Supplementation port (depending on system mode) at a nominal oxygen concentration of 90%. - Numeric 0.0 to 4.0 LPM ">4.0 LPM" shown for any value over 4.0 LPM.
Temp1	Patient Temperature 1	- Numeric 82.4 to 107.6° F (or 28 to 42° C) in steps of 0.1 or (- -) for unknown value. - Values outside of human body temperature range are shown as follows : "<28.0" or ">42.0" (Celsius) "<82.4" or ">108.0" (Fahrenheit)

Label	Items on the Main Screen	Description
Temp2	Patient Temperature 2	Same as Temp1.
PetCO₂	Partial Pressure of End-Tidal Carbon Dioxide	- Numeric 0–99 mmHg in steps of 1, or a series of dashes (- -) for unknown value; otherwise, ">99" is shown. - If unavailable due to calibration / warming-up, shows "CAL".
PiCO₂	Partial Pressure of Inspired Carbon Dioxide	Same as PetCO₂
SpO₂	Saturation percentage of oxygen attached to hemoglobin. (The amount of oxygen being carried by the red blood cells in the blood.)	- Numeric 0–100% in steps of 1 or a series of dashes (- -) for unknown value. - If >100, shows ">100". - If unavailable due to initialization or calibration shows "CAL".
NOTE: The following four (4) "Sp" items are only conditionally displayed when the appropriate Masimo sensor is present. To display these four (4) items, the device must first have the corresponding feature unlocked. Also, the proper Masimo sensor which measures the value must be present. In addition, the values can only be found on the left of the pleth real-time graph.		
SpCO	Saturation percentage of carbon monoxide attached to hemoglobin (i.e., carboxyhemoglobin). (The amount of carbon monoxide being carried by the red blood cells in the blood.)	- Shown as a percentage - If unavailable due to initialization or calibration shows "CAL".
SpMet	Saturation percentage of methemoglobin. (Methemoglobin [MetHb] is an oxidized form of hemoglobin that is unable to carry oxygen.)	- Shown as a percentage - If unavailable due to initialization or calibration shows "CAL".
SpHb	Total hemoglobin (Hb) concentration. (Hemoglobin is the part of a red blood cell that carries oxygen to the body. SpHb, which measures total hemoglobin, indicates the oxygen carrying capacity of the blood.)	- Shown as a ratio: g Hb / dl blood (grams hemoglobin / deciliter blood) - If unavailable due to initialization or calibration shows "CAL".
SpOC	Total oxygen content. (SpHb and SpO ₂ are used together to calculate the actual amount of oxygen in the blood.)	- Shown as a ratio: ml O₂ / dl blood (milliliters oxygen / deciliter blood) - If unavailable due to initialization or calibration shows "CAL".
n/a	Perfusion Index (Perfusion Index, or PI, is a relative assessment of the pulse strength at the monitoring site.)	0–20.00% displayed as bar graph beside SpO₂ - When there is no data, the graph is not drawn.

Label	Items on the Main Screen	Description
NIBP	Non-Invasive Blood Pressure The previous NIBP reading (if any), will be displayed while the new reading is being obtained. The current progress of the new reading will be displayed (if in progress) by showing the cuff pressure in a vertical bar graph beside the previous reading. The cuff pressure graph is from 0 to 260 mmHg (if over 260, 260 is displayed). If the user cancels the current NIBP reading, the previous one will be removed from the screen. If the NIBP reading is invalid and/or cannot be obtained, the previous NIBP reading will be erased.	• If either the systolic or diastolic value is over the range of the sensor, ">300" is shown for that value. • If the NIBP is one (1) or more minutes old, a message states "XXm ago" where XX is the number of minutes since the last reading. A NIBP measurement older than 15 minutes is discarded and no longer shown.  • If the NIBP reading was obtained using a cuff pressure that was close to the systolic value, which may make the reading lower than the actual systolic value of the patient, the systolic value is displayed with a series of dashes (- - -) to indicate the systolic pressure could not be obtained. The NIBP will automatically obtain another reading using a higher cuff pressure immediately after the reading. • If neither the systolic or diastolic pressure is available, a single set of dashes (- - -) is shown and no divider line is drawn. • While an NIBP reading is in progress, the current BP reading, if any, remains. As well, a bar graph representing cuff pressure from 0 to 260 mmHg is displayed to the right of the reading. • The current cuff pressure is displayed below the cuff pressure bar.

Label	Items on the Main Screen	Description
ABP, CVP, or ICP	IP Invasive Pressure There are three IP inputs on MOVES[®] SLC[™]. They can be used to measure one or more of Arterial Blood Pressure (ABP), Central Venous Pressure (CVP) or Inter-Cranial Pressure (ICP) ABP is displayed in a similar fashion to NIBP; that is, it is displayed showing systolic over diastolic in mmHg. CVP or ICP is displayed as a single number, the mean pressure with up to 3 digits. CVP is displayed in mmHg, and ICP is displayed in cmH₂O. The software will be able to detect if a port is in use via hardware. There can be none, one, two, or three IPs in use at any time.	- For ABP numeric, if either the systolic or diastolic value is over the range of the sensor, ">300" is shown for that value. If either the systolic or diastolic value is under the range of the sensor, "<-10" is shown for that value. - For CVP numeric, Numeric -10 to 300 mmHg in steps of 1 is shown. If greater than 300 mmHg, shows ">300". If less than -10 mmHg, "<-10" is shown. For an unknown value, a series of dashes (- - -) is shown. - For ICP numeric, Numeric -14 to 408 cmH₂O in steps of 1 is shown. If greater than 408 cmH₂O, shows ">408". If less than -14, "<-14" is shown. For an unknown value, a series of dashes (- - -) is shown. 

Label	Items on the Main Screen	Description
HR	Heart Rate	- The source of the HR is one of ECG, PulseOx, or IP in the following priority (if set to AUTO): - APB1 - APB2 - APB3 - SPO2 - ECG - HR Source: ABP or ECG - Numeric 0–250 BPM in steps of 1 or a series of dashes (- -) for an unknown value. ">250" is shown for values over 250 BPM - HR Source: SPO2 (Pulse Oximeter) - Numeric 0–239 BPM in steps of 1 or a series of dashes (- -) for an unknown value. ">239" is shown for values over 239 and below 260 BPM  - The user can also set the source by pressing the source button below the HR and selecting the source from a list. 
Chart Label 1	Chart area 1	- Button to select which chart to view. - Other buttons as required by the chart.
Chart Label 2	Chart area 2	Same as Chart 1

11.5.3 PC, PEEP and PIP

It may be noticed that although an option for “PC” is present on the Setup Screen, only the items “PIP” and “PEEP” are displayed on the Main Screen.




Figure 11-15: PC & PEEP on Setup Screen




Figure 11-16: PIP & PEEP on Main Screen

On ventilators with an active PEEP (Positive End Expiratory Pressure), such as MOVES[®] SLC[™], the user sets the PC (Pressure Control level), which is the pressure ABOVE PEEP, that you want maintained. If you adjust the PEEP value, the PIP (Peak Inspiratory Pressure) value goes up along with it.




Figure 11-17: Relations of PEEP to PIP to PC

However, monitoring (as it is shown on the Main Screen) always displays the true airway pressures: PIP and PEEP.

11.5.4 Inverted Display of Patient Monitoring Values

Normally, patient monitoring values on the Main Screen are shown in their representative colors (e.g., Temp 1 in tan, Temp 2 in turquoise) against a black background (see Figure 11-18). However, when an active alarm is associated with monitoring value, the values are shown in either white or black against a background of their representative color (see Figure 11-19).



NOTE: The active state of the alarm will be indicated by inverted display regardless of whether the alarm is on or off. For example, if the patient temperature is high, causing an alarm condition, the patient temperature is drawn in reverse color, even if the patient temperature alarm has been turned off by the user. If an alarm is latched, meaning it must be acknowledged by the user even if the alarm condition no longer active, the monitoring display is only inverted while the alarm condition is active.



Figure 11-18: Regular Temperature Display



Figure 11-19: Inverted Display of Temperature (Alarm Condition)

11.5.5 Invasive Pressure (IP) Source Buttons

There are three (3) physical IP connections on the patient sensor panel. Each one maps to its corresponding IP channel box on the monitoring screen (labeled IP1, IP2 and IP3). Each channel is independently configurable, and can be set for ABP, ICP or CVP mode. The current mode is displayed on the config button and may be changed by selecting the button and using navigation keys. An additional option is available, "ZERO", which does not change the mode, but rather forces an immediate re-zeroing of the channel.

11.5.6 Zeroing a Channel

An IP port will auto detect the insertion of an IP probe. An inserted probe requires zeroing before data is available, and will show ZERO REQD until the channel config button -- now titled "ZERO" -- is clicked, zeroing the channel. The channel box will then immediately display data with format/units appropriate for the selected mode. The user can change modes without re-zeroing. The port will also detect a probe disconnect and raise an alarm to notify the operator. This alarm may be dismissed to acknowledge the condition.

11.5.7 Graphs

Changing the mode will dynamically update the configuration (axes, units) of the corresponding real-time graph, if it is currently displayed. Trend graphs are maintained for each channel and mode combination independently.

11.6 ECG Screen

11.6.1 Overview

The ECG Screen will display up to 12 graphs (if all 10 ECG leads are used). Selecting 3, 5 or 12 on the Set Up Screen sets the minimum number of leads that need to be connected before ECG graph information is displayed. If more leads are connected than the number selected, additional graphs will still be shown (i.e., if 3 or 5 is chosen, but all leads are connected, all 12 graphs will be shown).



NOTE: The timeout to the Monitoring Screen is disabled while the ECG Screen is displayed.




Figure 11-20: ECG Screen – Pause Mode

The only user-adjustable value is Pause / Resume. Using the arrow buttons, the user navigates to it and presses the Check button to toggle between Pause / Resume. When the Check button is pressed when the button reads PAUSE then the ECG graphs will be frozen and can be examined more easily. Press the Check button when the button reads RESUME to restart normal graphing of the ECG waveform.




Figure 11-21: ECG Screen – Resume Mode

11.7 Alarm Limits Screen

11.7.1 Overview

The Alarm Limits Screen allows for control of alarm limit values. The selection order for the buttons on the Alarm Limits screen is the following: Screen Button, Alarm View Button, Limit List (limits are selected from top to bottom with no wrap).

A list of limits is displayed and the threshold name and value, along with the units of each limit, is shown. The list is scrollable once the Limit List part of the screen is selected.

If no limits have been changed from their default values, the phrase "All limits at default" is shown at right below the Status Bar. If any limits have been changed from their default values, an icon and a summary of the number of limits not at their default value is displayed (see *Figure 11-23: Alarm Limits Screen Active* on page 139).




Figure 11-22: Alarm Limits Screen – All at Default

11.7.2 Limit List Active

If the Limit List is active, a single limit is selected, as indicated by a yellow box around the limit (see Figure following). To modify the limit, the user presses the Check button to make the limit active and then uses the arrow buttons to select a value. The Check button is pressed once again to confirm the new value. The Cancel button is used to abandon an actively selected limit and its new value and return to the previously selected limit value.

Limit values that have been changed from their default are shown with an asterisk. Also, an icon in the status bar indicates that limit values have been changed and accompanying text lists how many. As well, a message at the right below the status bar show the number of limits not at their default values.



NOTE: For more information on the number of limits changed icon, see [Table 23: MOVES® SLCTM Status Bar Items and Descriptions](#) on [page 111](#).




Figure 11-23: Alarm Limits Screen Active

The following table lists all the limits, their available values, and their defaults.

Table 29: Alarm Limits and Defaults

Limits	Description
1. High airway pressure warning	20 to 58 cmH ₂ O in steps of 2, default to 32
2. High airway pressure release	40 to 58 cmH ₂ O in steps of 2, default to 40
3. Low expired PCO ₂	15 to 35 mmHg in steps of 1, default to 25
4. High expired PCO ₂	50 to 60 mmHg in steps of 1, default to 55
5. Low RR	5 to 15 B/M in steps of 1, default to 5

Limits	Description
6. High RR	20 to 40 B/M in steps of 1, default to 30
7. Low systolic BP	70 to 140 mmHg in steps of 1, default to 90
8. High systolic BP	140 to 200 mmHg in steps of 1, default to 180
9. Low SpO ₂	85 to 95% in steps of 1, default to 90
10. Low HR	30 to 60 BPM in steps of 1, default to 50
11. High HR	100 to 200 BPM in steps of 1, default to 120
12. Low CVP	-5 to 5 mmHg in steps of 1, default to 0
13. High CVP	10 to 30 mmHg in steps of 1, default to 20
14. Low ICP	-7 to 7 cmH ₂ O in steps of 1, default to 0
15. High ICP	13 to 40 mmHg in steps of 1, default to 27

11.8 Alarm On / Off Screen

11.8.1 Overview

The Alarm On / Off screen allows for alarm On / Off control. The selection order for Alarm ON/OFF screen buttons is Screen Button, Alarm View Button, ON/OFF list (ON/OFF list is selected from top to bottom with no wrap).




Figure 11-24: Alarm On / Off Screen

A summary of the alarms turned off is shown at top right as follows: All Alarms ON, 1 Alarm OFF, 2 Alarms OFF, etc.

The screen below presents a list of all alarms that can be turned off. The list is sorted when the user enters the screen according to the following:

1. First according to last active, with alarms that are active now at the top.
2. Alphabetically according to the alarm message.

Alarm level is indicated to the right. A red "H" indicates a HIGH priority alarm; a yellow "M" indicates a MEDIUM priority alarm; a yellow "L" indicates a LOW priority alarm. The letters "msg" denote a "MESSAGE", not an alarm.



NOTE: Alarms and messages shown with a padlock beside them in the alarm list CANNOT be turned off or altered in any way, even by someone with Administrator privileges.

Alarm On/Off List Active

If the list is active, one alarm in the list is selected, indicated by a yellow select box. The name of the alarm is the message associated with that alarm. Pressing the Next and Previous arrow buttons on the control panel selects the next or previous item in the list, scrolling the list as needed. Note that there is no wrap around.

description	on/off	last active	level
Low expired CO2	ON	NOW	H
Sample line disconnect	ON	NOW	M
Apnea	ON	3 min	H
CO2/O2 sensors warming up	ON	3 min	msg
Air pump failure	ON	never	
Alarm indication and keypad failure	ON	never	
Altitude outside of system limits	ON	never	
Backup power source required	ON	never	
Barometer fault	ON	never	
Battery A Unknown power level	ON	never	
Battery A not charging	ON	never	
Battery B Unknown power level	ON	never	

Figure 11-25: Alarm On / Off Screen Active

Pressing check toggles the ON/OFF state of the alarm and updates the alarm summary. Pressing the Cancel button exits the list active state.

THIS PAGE DELIBERATELY LEFT BLANK.

12.0 Using the Remote Screen

MOVES[®] SLC[™] can also be used with the remote screen (P/N 125725) shown below. The remote screen is an optional accessory.



Figure 12-1: MOVES[®] SLC[™] Remote Screen

12.1 Installing the Remote Screen Batteries

The MOVES[®] SLC[™] remote screen can be run on battery or AC power. The unit contains two batteries that are hot swappable, which means one battery can be replaced at a time without the unit turning off (as long as another battery with sufficient charge is present). There are battery bays on the left and right sides of the device.



Figure 12-2: Battery Bay 1 on Left Side




Figure 12-3: Battery Bay 2 on Right Side

1. To open the battery bay door, push down on the bottom lip of the recessed section in the middle.



2. Door open, bay ready to receive battery.



3. Orient the battery with the short yellow tab uppermost as shown in the photo at right.



4. Insert the battery as shown in the photo at right.



5. Fully seat the battery as shown in the photo at right.



6. Close the door, and push up on the lip to lock it.



7. **NOTE:** The USB port has been disabled and is non-functional.



12.2 Battery Indicator

The battery indicator in the following Figure indicates that both batteries are installed and are fully charged.




Figure 12-4: Remote Screen Battery Indicator

The following table explains the possible battery indications and the status they represent.

Table 30: Remote Screen Battery Status Indicators

Battery Indicator	Battery Status
Not lit	The battery pack is not inserted or not being charged.
Orange	Charging is in progress.
Green	The battery is fully-charged.
Blinking Green	When you close the battery cover with the battery pack already inserted, you can check the battery level by blinking times. Blinking times to Battery Level: • 5 times = 95 % to 100 % • 4 times = 50 % to 94% • 3 times = 25 % to 49% • 2 times = 5 % to 24% • 1 time = 0 % to 4 %
Red	The remaining battery is approximately 9% or less.
Blinking Red	When blinking approx. every 1 second: The battery pack or the charging circuit is not operating properly. When blinking approx. every 4 seconds: The battery cover is open. You can remove the battery pack in this case. When blinking approx. every 0.5 second: The battery cover is open. If you remove the battery pack in this case, the power supply is cut off and the computer is shut down. Close the battery cover immediately.
Blinking Orange	The battery cannot be charged temporarily due to the following reasons: • Its internal temperature is out of the acceptable range. • The power supply is not enough because software applications or peripheral devices are consuming a large amount of power.
Blinking Green and Orange Alternately	The temperature is low and the computer is warming up to prevent the hard disk drive or flash memory from malfunctioning. The computer will start automatically after warming-up.



NOTE: Once the battery is fully charged, the computer performs recharging only when the battery level becomes less than approximately 95%, so overcharging is avoided.

12.3 Connecting the Remote Screen to Wall Power

1. The cord shown at right, which includes an AC to DC adaptor, is used to connect wall power to the remote screen.



2. Locate the door labeled "DC IN 16V" on the right side of the remote screen next to battery bay 2. Push down on the bottom lip of the recessed section in the middle to open the door.



3. Insert the cylindrical end of power cord into receptacle shown at right.



4. Power cord connected.

NOTE: Installed batteries are automatically recharged if they are below a full charge when wall power is connected to the remote screen.




12.4 Connecting the Remote Screen to MOVES[®] SLC[™]

The remote screen comes with a cord one end of which is permanently attached to the screen. The other cylindrical end attaches to the auxiliary communications port on MOVES[®] SLC[™].


1. Remote screen cable connector used to connect to the auxiliary communications port on MOVES[®] SLC[™].
2. Auxiliary communications port on MOVES[®] SLC[™]. To access the connection lift the rubber cover.



3. The auxiliary communications port connection.



4. Align the two slot keys on the remote screen cable connector with the two receiving channels in the auxiliary communications port connection and insert.



5. Remote screen cable connector connected to the auxiliary communications port connection on MOVES[®] SLC[™].



12.5 The Remote Screen User Interface (UI)

12.5.1 When First Connecting

When the remote screen is physically connected to the auxiliary communications port on MOVES[®] SLC[™], and the MOVES[®] SLC[™] is running, once the remote screen is turned on it will interface automatically with MOVES[®] SLC[™] and, after a moment, begin sharing information. No software program needs to be launched by the user. Initially, however, the user will see a Thornhill Research Welcome screen followed by one of two Initial screens. Which Initial screen the user sees depends on whether the remote screen is starting from a total disconnect or from the hibernate/sleep state (i.e., whether software resources are loaded already or not).

If the remote screen is starting from a total disconnect, the user will see the following screen.




Figure 12-5: Initial Screen – No Resources Loaded

A blue progress bar is shown as resources are loaded. A green progress bar is shown while connecting to MOVES[®] SLC[™].

If the remote screen is starting from the hibernate/sleep state (i.e., software resources are already loaded), the user will see the following screen.



Figure 12-6: Initial Screen – Resources Already Loaded

When the remote screen starts from the hibernate/sleep state, resources are already loaded. Only the green progress bar is shown while connecting to MOVES® SLC™.



NOTE: Sometimes the connection can be so quick that this second screen is virtually invisible.

12.5.2 No System Test Screen

The remote screen does not display the System Test screen. If tests have not been completed prior to connecting the remote screen, the remote screen will open on the Setup screen and display the status "Initializing". The Suction and NIBP icons will be shown crossed out at the top of the screen, just left of center, since these two functions are not available until the startup tests have been passed. See Figure following. While the system is initializing (i.e., System Tests are not yet complete), the System Mode cannot be modified on the Setup screen; however, the settings on the Setup screen can be preconfigured in preparation for use after the System Tests are complete.




Figure 12-7: Setup Screen – Initializing

Once the startup tests have been passed, the Setup screen will open in Monitor Only mode and the crossed out Suction and NIBP icons will disappear from the top of the screen.




Figure 12-8: Setup Screen – Monitor Only Mode

As with the display screen physically attached to MOVES[®] SLC[™], the Setup screen “times out” or defaults to the Main screen after one minute of inactivity so that patient monitoring can be maintained.




Figure 12-9: Main Screen – On Remote Screen

12.5.3 Remote Screen Panel Buttons

Navigating and Selecting Via Panel Buttons

It will be immediately noticed that the actual physical panel buttons on MOVES[®] SLC[™] unit have been replicated at the bottom of the remote screen.




Figure 12-10: Remote Screen – Panel Buttons at Bottom

The remote screen is a touch-sensitive screen, and the user can navigate through the displayed user interface and select options with the touch-sensitive on-screen panel buttons in the same way as with the physical ones on the actual MOVES[®] SLC[™].

Touch Navigating

Rather than navigating step-by-step through the screens, choices and options as the user must do with the physical buttons, and can do with the on-screen buttons, when using the remote screen the user can also select initial choices immediately with a finger press. However, then the Next and Previous arrow buttons and the Check button must be used to define options and confirm choices.

12.5.4 Graphs Independently Configurable

Although most of the functionality on the remote screen is linked and analogous to that on the display screen physically attached to MOVES[®] SLC[™], the graphs are an exception. The graphs shown on the remote screen are NOT linked. Changing the graphs on the remote screen will not change the graphs on the display screen. This has been done to allow the user the option of displaying twice as much graph monitoring information.

12.5.5 Alarm Indicators

On the MOVES[®] SLC[™] unit, there are four (4) System Visual Indicators, one at each top corner. These indicators display the alarm system status (See *Section 10.3 System Visual Indicators (Alarms / System Status)* on page 104 for more details). On the remote screen, this functionality has been duplicated by Alarm Status Bars on either side of the screen (see Figures following).

The severity of the highest active alarm is shown. If an alarm is turned off, the display defaults to the severity of the next active alarm. The color-coding indicates the same alarm states.

Table 31: Alarm Side Bar States and Explanations

Status	Indicates
No bars	System off or in Initializing mode (i.e., System Tests incomplete)
Green bars	No alarms active
Yellow bars	Low or Medium Priority Alarm(s) active
Red bars	High Priority Alarm(s) active



Figure 12-11: Green Bars Indicating No Alarms Active



Figure 12-12: Red Bars Indicating High-Priority Alarm Active



Figure 12-13: Yellow Bars Indicating Low or Medium Priority Alarm Active

12.5.6 Remote Screen Software Version

The remote screen version can be determined by accessing the Advanced Screen. For information on accessing the Advanced Screen, see *Section 11.4 Advanced Screen* beginning on page 122.




Figure 12-14: Advanced Screen Showing Remote Screen Version

13.0 Connecting the Patient

13.1 Connection Overview

The patient may be intubated, in which case the patient will be connected to a Ventilator Breathing Circuit. The patient may also be breathing spontaneously but be receiving supplemental oxygen via a nasal cannula (nasal prongs) or oxygen mask.



WARNING! DO NOT CONNECT A PATIENT TO MOVES[®] SLC[™] UNTIL THE MOVES[®] SLC[™] SYSTEM IS PROPERLY WARMED UP AND O₂ VALUES ARE DISPLAYED.



WARNING! BEFORE VENTILATING A PATIENT, ENSURE THAT A SPARE VENTILATOR BREATHING CIRCUIT IS READILY AVAILABLE.

13.2 Connecting an Intubated Patient

1. Ensure that the Ventilator Cartridge and breathing circuit have been connected to MOVES[®] SLC[™].
2. Ensure Startup Test procedures have been completed and all tests passed.
3. Configure MOVES[®] SLC[™] to operate in Ventilate Mode.
4. Attach the elbow to the patient's endotracheal tube.

13.3 Connecting a Spontaneously Breathing Patient

1. Ensure Startup Test procedures have been completed and all tests passed.
2. Configure MOVES[®] SLC[™] to operate in O₂ Supplement mode.
3. Attach the oxygen supply tubing (connected to the nasal cannula or mask) to the MOVES[®] SLC[™] O₂ Outlet port.
4. Attach the nasal cannula or mask to the patient.

13.4 Attaching the Ear or Finger Clip Pulse CO-Oximeter Sensor

The MOVES[®] SLC[™] system comes with some or all of the pulse CO-oximeter sensors listed in the following table. Sensor position and whether the sensor is reusable or disposable is also indicated.

Table 32: CO-Oximeter Sensors

Sensor	Position	Reusable / Disposable
1. Adult SpO ₂ (only) Finger Clip (LNCS DCI)	Any of the five (5) fingers on either hand (but not the toes)	Reusable
2. SpO ₂ Ear Clip (LNCS TC-I)	Either ear lobe	Reusable
3. Adult SpO ₂ Adhesive Sensor (M-LNCS Adtx)	Any of the five (5) fingers on either hand (but not the toes)	Disposable (i.e., single use)
4. SpO ₂ & SpHb Finger Clip (Rainbow [®] DC-3 SC 360)	Any of the five (5) fingers on either hand (but not the toes)	Reusable

Sensor	Position	Reusable / Disposable
5. SpO ₂ , SpCO, & SpMet Finger Clip (Rainbow® DCI-dc3)	Any of the five (5) fingers on either hand (but not the toes)	Reusable
6. Adult SpO ₂ , SpHb and SpMet Adhesive Sensor (Rainbow® R1 25)	Any of the five (5) fingers on either hand (but not the toes)	Disposable (i.e., single use)
7. Adult SpO ₂ , SpCO and SpMet Adhesive Sensor (Rainbow® R25)	Any of the five (5) fingers on either hand (but not the toes)	Disposable (i.e., single use)

The following photographs show correct placement of the finger clip and ear clip sensors.




Figure 13-1: 1 of 5 Possible Sensor Finger Placements



Figure 13-2: Correct Ear Lobe Sensor Placement

AVOIDING INACCURATE READINGS

Pulse oximetry sensors work by transmitting red and infrared light through body tissue. Most light is absorbed by the tissue between the transmitting probe on one side and the actual sensor on the other. The small amount of light that is not absorbed is used to calculate oxygen saturation. Therefore, anything that disturbs the light flow can cause inaccurate readings. Note the following possible problems:

- **Motion Artifact** – The most common cause of inaccurate SpO₂ readings is movement. Movement affects the ability of the light to travel from the light-emitting diode (LED) to the photo detector. Rhythmic movement such as tremors and seizure activity, as well as shivering, and vibrations caused by ground or air transport, can cause problems with detecting saturation and may measure false high pulse readings. To overcome these problems, move the sensor to the ear as it is usually least affected by motion.
- **Ambient Light** – Because pulse oximeters measure the amount of light transmitted through arterial blood, bright light that shines directly on the sensor – whether from the sun or an overhead exam light – can skew the readings. To fix this problem, move the sensor, or cover it with something opaque.
- **Light Absorbent Substances** – Anything that absorbs light within the reading area may cause false-low readings (e.g., dried blood, polish, dyes [including intravenous]).

PULSE CO-OXIMETER SENSOR WARNINGS



WARNING! SPO₂ SENSOR APPLICATION SITES SHOULD BE INSPECTED AT LEAST EVERY 6 TO 8 HOURS TO ENSURE CORRECT SENSOR ALIGNMENT AND SKIN INTEGRITY. PATIENT SENSITIVITY MAY VARY DUE TO MEDICAL STATUS OR SKIN CONDITION. DISCONTINUE THE USE OF ADHESIVE TAPE STRIPS IF THE PATIENT EXHIBITS AN ALLERGIC REACTION TO THE ADHESIVE MATERIAL.



WARNING! OXIMETER READINGS MAY BE AFFECTED BY THE USE OF AN ELECTROSURGICAL UNIT (ESU).



WARNING! CIRCULATION DISTAL TO THE SENSOR SITE SHOULD BE CHECKED ROUTINELY.



WARNING! A FUNCTIONAL TESTER CANNOT BE USED TO ASSESS THE ACCURACY OF A PULSE OXIMETER MONITOR OR SENSOR.



WARNING! DO NOT USE DAMAGED SENSORS OR PATIENT CABLES. DO NOT USE A SENSOR OR PATIENT CABLE WITH EXPOSED OPTICAL OR ELECTRICAL COMPONENTS.

13.5 Placement of the Heart Rate Electrodes



WARNING! CONDUCTIVE PARTS OF ELECTRODES AND ASSOCIATED CONNECTORS FOR THE APPLIED PART, INCLUDING THE NEUTRAL ELECTRODE, SHOULD NOT CONTACT OTHER CONDUCTIVE PARTS AND EARTH.



CAUTION! ECG CABLES SHOULD BE DISCARDED AND REPLACED AFTER TWO (2) YEARS OF CONTINUOUS USE. CHECK CABLE USE BY RECORDING THE DATE THE CABLE WAS FIRST USED.

The snap-on ends of the MOVES[®] SLC[™] ECG cable are color coded (see Figure following).




Figure 13-3: Color-Coded 12-Lead ECG Cable

The electrodes usually consist of a conducting gel, embedded in the middle of a self-adhesive pad onto which cables clip. It is important that the ECG electrodes to which the cable clips attach be placed properly to reduce motion artifact and receive the best signal possible.

13.5.1 ECG Cable Color Coding and Naming Conventions

The clamps of the ECG cables that are attached to the electrodes are color coded to ease connection for medical personnel. There are two definitions for these colors: AHA colors, used in the USA and Canada, and IEC colors used in all other countries (usually Europe). The following table explains the two systems.

Table 33: ECG Cable Color Coding and Naming

Electrode	IEC Abbreviation	IEC Color	AHA Abbreviation	AHA Color
Right Arm	R	Red	RA	White
Left Arm	L	Yellow	LA	Black
Right Leg	N	Black	RL	Green
Left Leg	F	Green	LL	Red
Chest 1	C1	White-Red	V1	Brown-Red
Chest 2	C2	White-Yellow	V2	Brown-Yellow
Chest 3	C3	White-Green	V3	Brown-Green
Chest 4	C4	White-Brown	V4	Brown-Blue
Chest 5	C5	White-Black	V5	Brown-Orange
Chest 6	C6	White Violet	V6	Brown-Violet

13.5.2 Correct Electrode Placement

The following illustrations show correct electrode placement. The first illustration shows correct limb lead placement. The second illustration and accompanying table shows and explains the proper placement of the precordial leads V1–V6 (or C1–C6 in the IEC naming convention).




Figure 13-4: Proper Limb Lead Electrode Placement




Figure 13-5: Precordial Leads: V1–V6

Table 34: Proper Placement of Precordial Leads

Electrode Identifier	Electrode Position on Body Surface
V1 / C1	Fourth intercostal space at right border of sternum
V2 / C2	Fourth intercostal space at left border of sternum
V3 / C3	Fifth rib between V2 / C2 and V4 / C4
V4 / C4	Fifth intercostal space on left mid clavicular line
V5 / C5	Left anterior axillary line at the horizontal level of V4 / C4
V6 / C6	Left midaxillary line at the horizontal level of V4 / C4

13.5.3 Reducing Artifacts

A common problem with 12-lead ECG's is that a lot of artifacts are generated if the patient is moving around and not fully relaxed. Here are a few guidelines on how to reduce artifacts and capture a good ECG:

- Place the patient in a supine or semi-Fowler's position. If the patient cannot tolerate being flat, you can do the ECG in a more upright position.
- Instruct the patient to place their arms down by their sides and to relax their shoulders.
- Make sure the patient's legs are uncrossed.
- Remove any electrical devices, such as cell phones, away from the patient as they may interfere with the machine.
- If you're getting artifacts in the limb leads, try having the patient sit on top of their hands.

13.5.4 Expiry Date of ECG Electrodes



CAUTION! CHECK THE EMBOSSED EXPIRY DATE ON THE ECG ELECTRODES PACKAGE BEFORE USING. ELECTRODES ARE GOOD FOR 45 DAYS ONCE PACKAGE IS OPENED.



Figure 13-6: Electrode Package



Figure 13-7: Electrode Use By Date

13.6 Zeroing the Pressure in the IP Transducer

To zero the pressure in the invasive pressure transducer, turn the stopcock from the horizontal *running* position to the vertical *open-to-ambient-air* position.




Figure 13-8: Transducer Stopcock in Running Position (Horizontal)



Figure 13-9: Transducer Stopcock in Open-to-Ambient-Air Position (Vertical)
Zero Transducer

13.6.1 Zeroing the Transducer Channel

An IP port will auto detect the insertion of an IP probe. An inserted probe requires zeroing before data is available, and will show ZERO REQD until the onscreen channel config button -- now titled "ZERO" -- is clicked, zeroing the channel. The channel box will then immediately display data with format/units appropriate for the selected mode. The user can change modes without re-zeroing.

13.6.2 Transducer Warnings



EXTREME WARNING! WHEN USING FLUID FILLED PRESSURE TRANSDUCERS TO MONITOR INTRACRANIAL PRESSURE (ICP), MAKE SURE THAT THE TRANSDUCER AND THE LINE CONNECTING TO THE PATIENT'S DRAIN ARE FREE OF ANY AIR BUBBLES!



EXTREME WARNING! AFTER COMPLETING FILLING THE TRANSDUCER AND THE LINE, DISCONNECT THE FLUID BAG FROM THE TRANSDUCER, AND CAP THE END WITH THE STERILE CAP PRIOR TO CONNECTING THE TRANSDUCER TO THE PATIENT'S BRAIN!



EXTREME WARNING! NEVER FLUSH THE ICP TRANSDUCER WHILE CONNECTED TO THE PATIENT!



EXTREME WARNING! FAILURE TO OBSERVE THESE PRECAUTIONS MAY RESULT IN SERIOUS INJURY OR DEATH!

13.7 Using the Suction Feature

To use the Suction feature:

1. Ensure that the suction accessories are set up as shown in the section *Installing Suction Accessories* (see page 85).
2. Activate the suction by pressing the Suction Control Button .
3. Deactivate the suction at anytime by pressing the Suction Control Button a second time.



NOTE: The suction wand supplied with MOVES[®] SLC[™] has a control vent that must be occluded by the operator's thumb to allow suction to occur (See Figure following).




Figure 13-10: Suction Wand Showing Control Vent



NOTE: You can adjust suction pressure between 100 and 325 mmHg in increments of 25 mmHg. (For details, see the subsection *Changing Settings* in the section *Setup Screen* that begins on page 116.)



CAUTION! THE OPERATOR SHOULD ALWAYS HAVE AVAILABLE AN ALTERNATE MEANS OF SUCTION IN THE EVENT OF POWER FAILURE, MECHANICAL FAILURE, OR A SERIOUS OCCLUSION IN THE SUCTION SYSTEM.



WARNING! THE O₂ CONCENTRATOR DOES NOT PRODUCE O₂ WHEN THE SUCTION FEATURE IS IN USE. AN ALTERNATIVE O₂ SUPPLY WILL BE NECESSARY IF HIGH O₂ CONCENTRATION IS CRITICAL.

13.8 Shutdown Procedures

1. To deactivate the MOVES[®] SLC[™] system, press and **hold** the Power Control button  for five (5) seconds until the display turns OFF.
2. Detach all patient monitoring accessories from the patient connection panel.
3. Dispose of breathing circuit, sample lines (**but NOT THE NAFION TUBE**) and filters, breathing cartridge, ECG adhesive sensor pads, IP transducer, suction wand and tubing.
4. Clean and sterilize all cables, NIBP cuffs and tubing with a disinfecting spray. Ensure that no moisture enters the tubing. Return items to storage in the accessories case.

5. Carefully remove the suction canister and suction filter and dispose of them and waste in a sanitary manner in accordance with local biohazard regulations.
6. Unlatch the MOVES[®] SLC[™] batteries or remove them from the system and store them in a dry area at room temperature.



WARNING! LEAVING BATTERIES INSTALLED IN A NON-OPERATIONAL UNIT MAY CAUSE THEM TO DRAIN TO AN UNRECHARGEABLE LEVEL.



CAUTION! THE SUCTION CANISTER AND SUCTION FILTER ARE INTENDED FOR SINGLE USE ONLY AND SHOULD BE DISPOSED OF IN ACCORDANCE WITH LOCAL BIOHAZARD REGULATIONS.



CAUTION! BREATHING CIRCUIT, SAMPLE LINES (BUT NOT THE NAFION TUBE) AND FILTERS, BREATHING CARTRIDGES, ECG ADHESIVE SENSOR PADS, IP TRANSDUCER AND SUCTION WAND AND TUBING ARE DISPOSABLE AND SHOULD BE DISPOSED OF IN ACCORDANCE WITH LOCAL BIOHAZARD REGULATIONS.

THIS PAGE DELIBERATELY LEFT BLANK.

14.0 Using System Graphs / Trends

14.1 Overview

Plotting appears on two graphs on the Main Screen. Either graph can be set to TRENDS, which plots historical values over a user-designated time frame (30 minutes, 1, 2, 4, 8, and 16 hours).




Figure 14-1: Plotting Graphs on Main Screen (detail)

The graphs have a vertical scale (from a minimum to a maximum value) for data of interest, and a horizontal time scale. Because plot data is stored, changing plots (e.g., from ECG to Airway) allows the part of the plot to be displayed immediately. Normally, patient data is displayed over time from left to right; however, if there is no data, the sweep continues to erase old data and plots new data when it is available.



NOTE: Data points above the top of the plot are plotted at the top of the plot. Data points below the bottom of the plot are plotted at the bottom of the plot. In other words, if the data is out of range, the display will only show either maximum or minimum values until the data comes back into range.

14.2 Available System Graphs / Trends

14.2.1 System Graphs

The following table lists all available System Graphs.

Table 35: System Graphs and Parameters

GRAPH	PARAMETERS
1. ECG – 3 mV	Minimum value: -1.5 mV Maximum value: 1.5 mV Plotting rates: Fast (27 mm/sec), Slow (10.8 mm/sec) Sampling Rates: 150 Hz (fast); 60 Hz (slow) Reference line: 0 mV Time scale: Real time
2. ECG – 6 mV	Minimum value: -3.0 mV Maximum value: 3.0 mV Plotting rates: Fast (27 mm/sec), Slow (10.8 mm/sec) Sampling Rates: 150 Hz (fast); 60 Hz (slow) Reference line: 0 mV Time scale: Real time
3. ECG – 12 mV	Minimum value: -6.0 mV Maximum value: 6.0 mV Plotting rates: Fast (27 mm/sec), Slow (10.8 mm/sec) Sampling Rates: 150 Hz (fast); 60 Hz (slow) Reference line: 0 mV Time scale: Real time
4. IP–ABP	Minimum value: 60 mmHg Maximum value: 160 mmHg Plotting rate: 7.7 mm/sec Sampling Rate: 20.8 Hz Reference line: 100 mmHg
5. IP–CVP	Minimum value: 0 mmHg Maximum value: 20 mmHg Plotting rate: 7.7 mm/sec Sampling Rate: 20.8 Hz Reference line: 10 mmHg
6. IP–ICP	Minimum value: 0 mmHg Maximum value: 27 mmHg Plotting rate: 7.7 mm/sec Sampling Rate: 20.8 Hz Reference line: 14 mmHg
7. Plethysmograph – Pleth	Plotting rate: 7.7 mm/sec Sampling Rate: 20.8 Hz Reference line: None  NOTE: The pleth waveform is scaled to a fixed size for signal strengths above 10% or below 0.5% (i.e., all waveforms with a signal strength over 10% or less than 0.5% will have the same amplitude; however, the user will still see a oscillating waveform). Between the values indicated, the waveform is scaled according to signal strength.

GRAPH	PARAMETERS
8. Airway	Minimum value: 0 cmH ₂ O Maximum value: 40 cmH ₂ O Plotting rate: 2.66 mm/sec Sampling Rate: 7.4 Hz Reference line: 30 cmH ₂ O
9. pCO ₂	Minimum value: 0 mmHg Maximum value: 50 mmHg Plotting rate: 4.00 mm/sec Sampling Rate: 11.1 Hz Reference line: 40 mmHg Time scale: Real time

14.2.2 System Trends

There are nine (9) trend charts:

- FiO₂ / SpO₂
- PIP
- TEMPS (both TEMP1 and TEMP2)
- HR
- ABP (available under IP1, IP2, or IP3)
- CVP (available under IP1, IP2, or IP3)
- ICP (available under IP1, IP2, or IP3)
- NIBP
- PetCO₂

Trend plots can be set to display trend data for:

- 30 minutes
- 1 hour
- 2 hours
- 4 hours
- 8 hours
- 16 hours

The following table lists all System Trends.

Table 36: System Trends and Parameters

TREND	PARAMETERS
1. Trend Plot – ABP (available under IP1, IP2 or IP3)	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 60–160 mmHg
2. Trend Plot – CVP (available under IP1, IP2 or IP3)	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0–20 mmHg
3. Trend Plot – ICP (available under IP1, IP2 or IP3)	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0–27 cmH ₂ O
4. Trend Plot - NIBP	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 60–160 mmHg
5. Trend Plot – PetCO ₂	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0–50 mmHg
6. Trend Plot - FiO ₂ / SpO ₂	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0 to 100%
7. Trend Plot - PIP	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0 to 40 cmH ₂ O
8. Trend Plot – TEMP	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 89°F–108°F (32°C–42°C)
9. Trend Plot – HR	Time: 30 min, 1, 2, 4, 8, 16 hours Range: 0 to 200 bpm

15.0 Alarms

Alarm conditions which require immediate or priority user attention can be distinguished from normal status alarms at a distance via their highlighted screen colors, audible annunciation, and alarm status LEDs. If a user is too far away from MOVES[®] SLC[™] to see the contents of the display, the user will still be informed of the alarm condition by seeing that the status LEDs are NOT green.

15.1 About Status LEDs

There are four (4) status LEDs, each situated on a top corner of the MOVES[®] SLC[™]. They are placed away from the display so that they can be seen when the display cannot. Note the following:

- When the system is NOT in alarm, the status LEDs are solid green.
- When the system is IN alarm, the status LEDs are solid yellow, flashing yellow, or flashing red (from low to medium to high in order of priority).
- If alarms have been silenced, the status LEDs continue to indicate alarm level.

Alarms and messages are displayed in the top right of the screen's status area. When a new alarm condition becomes active, or an existing alarm becomes inactive, the system will update the visual and audible alarm state with the highest priority active alarm. To see the text of an alarm other than that of the highest-priority alarm, the user can navigate to the Alarm Queue on the status bar, press the Check button, and then use the Next and Previous buttons to page through the alarms.

15.2 Alarm Priorities and Characteristics

15.2.1 Standard Alarms

The following table details the MOVES[®] SLC[™] standard alarms, their audio and visual characteristics, and their priorities.

Table 37: Alarm Types and Descriptions

ALARMS TABLE			
Alarm Priority	Description / Response	Visible Alarm	Audible Alarm
High	Immediate operator response required	Red LED 60% duty cycle flashing at 2.8 Hz	Three short notes of ascending pitch followed by two long notes of ascending pitch. The musical sequence is a common triad (C-E-G) followed by a perfect fourth (G-C).
Medium	Prompt operator response required	Yellow LED 60% duty cycle at 0.8 Hz	Three notes of equal, medium length and ascending pitch. The musical sequence is a common triad (C-E-G).
Low	Operator awareness required	Yellow LED Constant ON	None
Information	Message	Green LED (Indicating no Alarm) Constant ON	None



NOTE: Low priority alarms and/or messages will not be audibly indicated.

15.2.2 High Priority Communication Failure Alarm

The High Priority Communication Failure alarm is triggered when there is no longer communication between the MOVES[®] SLC[™] user interface and its internal systems.

Alarm Priority	Description	Visible Alarm	Audible Alarm
High	Alarm system failure. Immediate operator response required	Red LED 60% duty cycle flashing at 2.8 Hz	A high short note followed by two short low notes an octave lower, then a high long note followed by a low long note an octave lower.

15.3 Alarm Queue

The Alarm Queue is displayed in the top right corner of the screen's Status Bar. Initially, only the text of the highest-priority alarm is displayed. To see the text of an alarm other than that of the highest-priority alarm, the user navigates to the Alarm Queue, presses the Check button, and then uses the Next and Previous buttons to page through the alarms. The Alarm Queue shows all active alarms, non-active latched alarms, and non-alarm messages. If there are no alarms, the message NO ALARMS is displayed.




Figure 15-1: No Alarms Message in Alarm Queue

15.4 Locked Alarms and Messages

Some alarms and messages CANNOT be turned off or altered in any way, even by someone with Administrator privileges. These alarms and messages are shown with a padlock beside them in the alarm list on the Alarm ON/OFF screen. Selecting them and pressing the Check button will have no effect on their status (unlike other alarms which can be toggled ON / OFF or Dismissed).




Figure 15-2: Two Locked Alarms

15.5 Inhibitable and Latching Alarms

Inhibitable and Latching alarms can be distinguished from General alarms by the “dismiss” option initially shown in the top right corner.

15.5.1 Inhibitable Alarms

Some, but not all, MOVES[®] SLC[™] alarms are Inhibitable. An Inhibitable alarm can be dismissed and cleared even though the condition or event that triggered the alarm is still active.

An Inhibitable alarm is shown in two states in the images below.




Figure 15-3: Inhibitable Alarm Able To Be Dismissed




Figure 15-4: Inhibitable Alarm Cleared

To clear the Inhibitable alarm, select it in the Alarm Queue and then press the Check button. Once an alarm is cleared, it will be removed from the Alarm Queue.



NOTE: The alarm will, of course, also be cleared if the condition or event that triggered the alarm is addressed or corrected.

15.5.2 Latching Alarms

Some alarms (known as latching alarms) require operator acknowledgment before they are removed from the Alarm Queue (i.e., they will persist in the Alarm Queue even if the condition or event that triggered the alarm is addressed / corrected). This guarantees that the operator is made aware that the alarm condition occurred **even if the alarm condition has since disappeared**. A latching alarm is visually identical to an Inhibitable alarm. A latching alarm is shown in two states in the images below.




Figure 15-5: Latched Alarm Able To Be Dismissed



Figure 15-6: Latched Alarm Cleared

To clear the Latched alarm, select it in the Alarm Queue and then press the Check button. Once an alarm is cleared, it will be removed from the Alarm Queue. A Latched alarm will also not persist after MOVES[®] SLC[™] is shut down.



NOTE: All alarms are logged in the log file.



NOTE: See the extensive Alarms Table which follows for identification of which alarms are Inhibitable, which are Latched, and which are General and neither.

15.5.3 Alarms That Have Been Turned Off

Alarms that are not Locked can be turned OFF using the Alarm ON/OFF screen. A white triangle with an "X" through it presented on a red hexagonal shape indicates that some alarms have been turned off. This icon is shown in the top left corner of the Status Bar. The number of alarms turned off is indicated below the icon.




Figure 15-7: Latched Alarm Able To Be Dismissed

15.6 Alarms Conditions and Causes

Table 38: Alarm Conditions and Causes

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
1. Internal fault: Service soon	Low	Latching		Low priority fault indicating possible calibration issues with patient monitoring and/or system sensors, concentrator pump temperature sensor fault or internal memory malfunction.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
2. Low Temperature: System warming ...	Low			System temperature is less than 0 C (32 F) or concentrator pump temperature is less than 0 C (32 F) while system temperature is in the range of -18 C (-0.4 F) to 70 C (158 F).	1. Wait for system to warm up or proceed to a higher temperature environment. 2. If erroneous, return MOVES[®] SLC[™] for servicing.
3. System Temperature too high	Low			System temperature exceeds 70 C (158 F).	1. Proceed to a lower temperature environment. 2. If erroneous, return MOVES[®] SLC[™] for servicing.
4. System Temperature too low	Low			System temperature is less than -18 C (-0.4 F).	1. Proceed to a higher temperature environment. 2. If erroneous, return MOVES[®] SLC[™] for servicing.
5. TEMP1: High patient temperature	Low			Patient temperature is over 38 C (100.4 F) and no more than 42 C (107.6 F).	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect temperature patient connection. 3. Check for related alarms and follow actions if required. 4. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
6. TEMP1: Low patient temperature	Low			Patient temperature is no less than 28 C (82.4 F) and under 35 C (95.0 F).	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect temperature probe patient connection. 3. Check for related alarms and follow actions if required. 4. Clinical intervention may be required.
7. TEMP1: Temperature probe disconnect	Low		Inhibitible	Temperature probe patient connection is disconnected from either system or patient after previously being connected.	1. Fully disconnect then reconnect temperature probe patient connection. 2. Check integrity of interface with patient. 3. Replace temperature probe patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
8. TEMP2: High patient	Low			Patient temperature is over 38 C (100.4 F) and no more than 42 C (107.6 F).	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect temperature patient connection. 3. Check for related alarms and follow actions if required. 4. Clinical intervention may be required.
9. TEMP2: Low patient	Low			Patient temperature is no less than 28 C (82.4 F) and under 35 C (95.0 F).	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect temperature probe patient connection. 3. Check for related alarms and follow actions if required. 4. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
10. TEMP1: Temperature probe disconnect	Low		Inhibitible	Temperature probe patient connection is disconnected from either system or patient after previously being connected.	1. Fully disconnect then reconnect temperature probe patient connection. 2. Check integrity of interface with patient. 3. Replace temperature probe patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
11. Concentrator temperature sensor fault	Low			Redundant internal temperature sensor communication failure with concentrator PCB .	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
12. Concentrator failure	Medium			Concentrator PCB communication failure with UI or concentrator PCB reports concentrator pump fault.	1. Turn suction on and off for 2 seconds. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.
13. Concentrator in degraded mode	Medium			Concentrator PCB reports partial failure due to unexpected data, including: product tank flow, product tank pressure and suction pressure. Oxygen flow or concentration may be reduced.	1. Monitor inspired O₂ to ensure adequate oxygenation is being maintained in the ventilation circuit. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.
14. Suction not available	Medium			Concentrator PCB communication failure with UI or concentrator PCB reports concentrator pump fault.	1. Turn suction on and off for 2 seconds. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
15. Ventilator failure: Use backup vent if attempting to ventilate	HIGH, if vent problem while attempting to ventilate LOW otherwise	Latching		Communication failure between ventilator module and UI, erroneous feedback from blower to ventilator module, blower air inlet occlusion, ventilator pressure sensor fault, or low supply voltage to ventilator module while in ventilation mode (1) if ventilator not running or (2) if ventilator is running. NOTE: In the event of a communication failure, the ventilator will continue to operate with the last user-set parameters and will operate the air pump at least 9 LPM. NOTE: If a "Leak alarm" and a "Ventilator failure" alarm occur simultaneously, the leak is only substantial if you see alternating Ve/Vi values displayed.	1. Ensure that the blower inlet is not occluded. 2. Power cycle MOVES® SLC™. 3. Return MOVES® SLC™ for servicing.
16. Vent vol-flow fault: Try PC-IMV mode if venting and not already in PSV or PC-IMV	Medium			Ventilator flow sensor communication failure with ventilator PCB detected during system start-up tests.	1. Power cycle MOVES® SLC™. 2. Repeat system start-up tests. 3. Revert to pressure control ventilation modes if necessary. 4. Return MOVES® SLC™ for servicing.
17. Air pump: failure	Medium			Erroneous feedback from air pump is detected by concentrator PCB.	1. Power cycle MOVES® SLC™. 2. Return MOVES® SLC™ for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
18. Leak detected	Medium			Exhaled tidal volume is more than 100 mL less than the delivered inspired tidal volume, or the internal ventilator bag has collapsed. NOTE: A high expiratory flow rate (over 80 LPM) can be misinterpreted as being lower than actual, resulting in an incorrect leak report. NOTE: If a "Leak alarm" and a "Ventilator failure" alarm occur simultaneously, the leak is only substantial if you see alternating Ve/Vi values displayed.	1. Ensure all breathing circuit connections (including patient tubing, Y-piece and sample line) are tight and no leaks are detected. 2. Ensure ventilator cartridge is properly connected to MOVES[®] SLC[™]. 3. Replace breathing circuit including ventilator cartridge. 4. Leak may be patient-related. Clinical intervention may be required. 5. Return MOVES[®] SLC[™] for servicing.
19. Internal power monitoring fault	Medium			Internal supply voltage mismatches between any or all of the ventilator PCB, UI PCB, concentrator PCB or power manager PCB.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
20. Power manager fault	Medium			Communication failure between power manager PCB and UI.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
21. Battery A: Unknown power level	Low			Power manager PCB detects presence of battery but fails to communicate or power manager fault alarm is active.	1. Replace affected battery or batteries. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.
22. Battery B: Unknown power level	Low			Power manager PCB detects presence of battery but fails to communicate or power manager fault alarm is active.	1. Replace affected battery or batteries. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
23. Battery A: Not charging	Low			Battery charge level has not increased by more than 1% after one hour of charging and is less than 96%.	1. Remove and replace problematic battery. 2. Replace MOVES® SLC™ power supply/charger. 3. Power cycle MOVES® SLC™. 4. Return battery and/or power supply/charger and/or MOVES® SLC™ for servicing.
24. Battery B: Not charging	Low			Battery charge level has not increased by more than 1% after one hour of charging and is less than 96%.	1. Remove and replace problematic battery. 2. Replace MOVES® SLC™ power supply/charger. 3. Power cycle MOVES® SLC™. 4. Return battery and/or power supply/charger and/or MOVES® SLC™ for servicing.
25. Backup power source required	High			Only one valid power source (battery A, battery B or power supply/charger) is detected. Battery charge level must be greater than or equal to 10% to be considered valid.	1. Connect a second valid power source to MOVES® SLC™. 2. If second valid power source is connected and alarm persists, ensure that battery doors are fully closed and/or power supply/charger connection is secure.
26. Low reserve power level	Medium			Battery charge level is less than 20% when two valid power sources (battery A, battery B or power supply/charger) are detected.	• Connect power supply/charger to MOVES® SLC™ to allow battery or batteries to charge.
27. Critically low reserve power	High			Battery charge level is less than 10% when two valid power sources (battery A, battery B or power supply/charger) are detected.	• Connect power supply/charger to MOVES® SLC™ to allow battery or batteries to charge.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
28. Patient monitor failure	Medium			Patient monitors PCB communication failure with MOVES [®] SLC [™] .	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
29. ECG lead off	Low		Inhibitible	ECG lead is disconnected from patient after previously being connected.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect ECG patient connection. 3. Replace ECG patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
30. ECG fault	Medium			ECG module internal fault or communication failure with patient monitors PCB.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
31. Pulse Ox: Cable disconnect	Low		Inhibitible	Pulse oximeter extension cable or patient connection is disconnected from system after previously being connected.	1. Fully disconnect then reconnect pulse oximeter extension cable or patient connection to the system.. 2. Replace pulse oximeter extension cable and/or patient connection. 3. Power cycle MOVES[®] SLC[™]. 4. Return MOVES[®] SLC[™] for servicing.
32. Pulse Ox: Replace cable	Low			Pulse oximeter extension cable has expired and needs to be replaced.	• Replace pulse oximeter extension cable with another that has not expired.
33. Pulse Ox: Invalid cable	Low			Pulse oximeter extension cable is not compatible with system.	• Replace pulse oximeter extension cable with another that is compatible with MOVES[®] SLC[™].

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
34. Pulse Ox: Sensor disconnect	Low		Inhibitible	Pulse oximeter patient connection cable is disconnected from extension cable connected to system after previously being connected.	1. Fully disconnect then reconnect pulse oximeter patient connection cable to extension cable. 2. Fully disconnect then reconnect pulse oximeter extension cable to system. 3. Replace pulse oximeter extension cable and/or patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
35. Pulse Ox: Replace sensor	Low			Pulse oximeter patient connection has expired and needs to be replaced.	• Replace pulse oximeter patient connection with another that has not expired.
36. Pulse Ox: Invalid sensor	Low			Pulse oximeter patient connection is not compatible with system.	• Replace pulse oximeter patient connection with another that is compatible with MOVES® SLC™.
37. PulseOx: Patient disconnect	Low		Inhibitible	Pulse oximeter patient connection is disconnected from patient.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
38. PulseOx: Check sensor	Low			Pulse oximeter connection at system or patient is compromised.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
39. PulseOx: Interference detected	Low			Disturbance of pulse oximeter measurement at patient connection due to ambient or environmental light interference.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to patient. 3. Reduce or eliminate ambient or environmental light interference between pulse oximeter patient connection and possible outside sources. 4. Replace pulse oximeter patient connection. 5. Power cycle MOVES[®] SLC[™]. 6. Return MOVES[®] SLC[™] for servicing.
40. PulseOx: Low perfusion	Low		Inhibitible	Pulse oximeter detects low perfusion related to SpO ₂ measurement.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
41. PulseOx: Low SpCO perfusion	Low		Inhibitible	Pulse oximeter detects low perfusion related to SpCO measurement.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Clinical intervention may be required.
42. PulseOx: Low SpMet perfusion	Low		Inhibitible	Pulse oximeter detects low perfusion related to SpMet measurement.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Clinical intervention may be required.
43. PulseOx: Low SpHb perfusion	Low		Inhibitible	Pulse oximeter detects low perfusion related to SpHb measurement.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Clinical intervention may be required.
44. PulseOx: Low SpOC perfusion	Low		Inhibitible	Pulse oximeter detects low perfusion related to SpOC measurement.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitable	Cause	Action(s)
45. PulseOx: SpO2 reading confidence poor	Low			Pulse oximeter indicating low confidence in SpO ₂ measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
46. PulseOx: HR reading confidence poor	Low			Pulse oximeter indicating low confidence in heart rate measurement value when source of heart rate displayed by the system is the pulse oximeter.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. If available, change heart rate source on MOVES[®] SLC[™]. 5. Power cycle MOVES[®] SLC[™]. 6. Return MOVES[®] SLC[™] for servicing.
47. PulseOx: PI reading confidence poor	Low			Pulse oximeter indicating low confidence in perfusion index measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
48. PulseOx: SpCO reading confidence poor	Low			Pulse oximeter indicating low confidence in SpCO measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
49. PulseOx: SpMet reading confidence poor	Low			Pulse oximeter indicating low confidence in SpMet measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
50. PulseOx: SpHb reading confidence poor	Low			Pulse oximeter indicating low confidence in SpHb measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
51. PulseOx: SpOC reading confidence poor	Low			Pulse oximeter indicating low confidence in SpOC measurement value.	1. Check integrity of interface with patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
52. PulseOx: Only SpO2 available	Low		Inhibitible	SpCO, SpMet, SpHb and SpOC measurements cannot be obtained.	1. Fully disconnect and reconnect pulse oximeter patient connection to patient. 2. Fully disconnect then reconnect pulse oximeter patient connection to system. 3. Replace pulse oximeter patient connection. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
53. PulseOx fault	Low			Pulse oximeter module communication failure with patient monitor PCB.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
54. Gas sensor pressure fault	Medium			Mismatch between O ₂ and CO ₂ cell pressures.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
55. CO2 sensor fault	Medium			CO ₂ sensor internal fault, communication failure with patient monitor PCB or calibration issue.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
56. O2 sensor fault	Medium			O ₂ sensor internal fault, communication failure with patient monitor PCB or calibration issue.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
57. Atmospheric O ₂ out of range	Medium			Ambient O ₂ measurement during patient gas sensor calibration is outside of expected range (20.9 +/- 0.4 %).	1. Ensure there are no obstructions to the gas calibration port on the patient monitor panel. 2. Power cycle MOVES[®] SLC[™] and wait for re-calibration. 3. If unsuccessful, move MOVES[®] SLC[™] into another area, power cycle MOVES[®] SLC[™] and wait for re-calibration. 4. Return MOVES[®] SLC[™] for servicing.
58. Barometer fault	Low			Mismatch between O ₂ and CO ₂ cell pressures and barometric pressure when sample pump is not running.	1. Power cycle MOVES[®] SLC[™]. 2. Return MOVES[®] SLC[™] for servicing.
59. Altitude above system limits	Medium			Barometric pressure is less than 60 kPa.	1. Proceed to a lower altitude. 2. If erroneous, return MOVES[®] SLC[™] for servicing.
60. Altitude below system limits	Medium			Barometric pressure is greater than 110 kPa.	1. Proceed to a higher altitude. 2. If erroneous, return MOVES[®] SLC[™] for servicing.
61. Inspired O ₂ low	High			Inspired O ₂ is less than 19% (while running in O ₂ supplement and ventilation modes).	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required. 3. Return MOVES[®] SLC[™] for servicing.
62. Inspired O ₂ above target	Low			Inspired O ₂ exceeds setting by 15% (ventilation mode only).	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
63. Inspired O ₂ below target	Low			Inspired O ₂ is not reached within 5% of setting while running in ventilation mode.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required. 3. Return MOVES[®] SLC[™] for servicing.
64. Low expired CO ₂	Low			When the apnea alarm is not active, expired CO ₂ partial pressure is (1) less than or equal to the user set low expired CO ₂ alarm limit and greater than or equal to 15 mmHg or (2) less than 15 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
65. High expired CO ₂	Medium			Expired CO ₂ partial pressure is greater than or equal to the high expired PCO ₂ alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
66. High inspired CO ₂ : Change cartridge	Low			While running in ventilate mode, inspired CO ₂ partial pressure is (1) greater than 6 mmHg or (2) greater than or equal to 10 mmHg.	1. Replace ventilator cartridge. 2. Return MOVES[®] SLC[™] for servicing.
67. Low SpO ₂	High			O ₂ saturation level is less than or equal to user-set low SpO ₂ alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
68. Apnea	High			Respiratory rate is less than or equal to 4 breaths per minute and the sample line occlusion alarm is not active.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitable	Cause	Action(s)
69. Low RR	Medium			Respiratory rate is less than or equal to user set low RR alarm limit when apnea alarm is not active and respiratory rate is greater than zero.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
70. High RR	Medium			Respiratory rate is greater than or equal to user set high RR alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
71. High pulse rate	Low			Heart rate is (1) greater than user set high HR alarm limit or (2) greater than or equal to 150 beats per minute and greater than user set high HR alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
72. Low pulse rate	High			Heart rate is less than user set low HR alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
73. Cannot reach set I/E ratio	Low		Inhibitable	Ventilator operational settings and/or patient parameters result in user set I:E ratio not being achieved due to inadequate inspiratory time.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
74. Above target volume	Medium			Exhaled tidal volume is at least 15% greater than the current user set point.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
75. Below target volume	Medium			Exhaled tidal volume is at least 15% less than the current user set point.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
76. High expired tidal volume	Medium			Exhaled tidal volume is 30% greater than the delivered inspired tidal volume.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
77. NIBP: Low systolic pressure	High		Inhibitible	Systolic pressure of NIBP measurement is less than or equal to the user set low systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
78. NIBP: High systolic pressure	Medium		Inhibitible	Systolic pressure of NIBP measurement is greater than or equal to the user set high systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
79. NIBP: Low diastolic pressure	Medium		Inhibitible	Diastolic pressure of NIBP measurement is less than or equal to 40 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
80. NIBP: High diastolic pressure	Medium		Inhibitible	Diastolic pressure of NIBP measurement is greater than or equal to 100 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
81. NIBP inflate failure	Low		Inhibitible	Pressure required for inflation cannot be generated in the NIBP patient arm cuff.	1. Ensure NIBP patient arm cuff is full connected to MOVES[®] SLC[™] and has no leaks. 2. Replace NIBP patient arm cuff. 3. Power cycle MOVES[®] SLC[™]. 4. Return MOVES[®] SLC[™] for servicing.
82. NIBP could not get reading	Low		Inhibitible	NIBP measurement was not obtained unrelated to NIBP patient arm cuff inflation failure.	1. Attempt another NIBP reading. 2. Ensure NIBP patient arm cuff is full connected to MOVES[®] SLC[™] and has no leaks. 3. Replace NIBP patient arm cuff. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
83. NIBP set to: Manual update	Low	Locked & Latching	Inhibitible	NIBP measurement was not obtained three consecutive times when set to take automatic measurements.	1. Attempt another NIBP reading. 2. Ensure NIBP patient arm cuff is full connected to MOVES[®] SLC[™] and has no leaks. 3. Replace NIBP patient arm cuff. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
84. NIBP fault	Low			NIBP module internal fault or communication failure with patient monitors PCB.	1. Power cycle MOVES[®] SLC[™]. 2. Attempt another NIBP reading. 3. Return MOVES[®] SLC[™] for servicing.
85. NIBP Mode changed: STAT -> 5 min	Low	Locked & Latching		NIBP has been set to STAT mode for 15 minutes and has automatically changed to take readings every 5 minutes for long term monitoring.	• No action required. NIBP mode setting can be changed back to STAT at any time if required.
86. IP1: Measurement error	Low			IP measurement contains values that are outside of the measurement range of the sensor.	1. Fully disconnect then reconnect IP patient connection and zero IP transducer as prompted. 2. Replace IP patient connection. 3. Check integrity of interface with patient. 4. Power cycle MOVES[®] SLC[™]. 5. Return MOVES[®] SLC[™] for servicing.
87. IP1 disconnect	Low		Inhibitible	IP patient connection is disconnected from system after previously being connected.	1. Fully disconnect then reconnect IP patient connection. 2. Replace IP patient connection. 3. Power cycle MOVES[®] SLC[™]. 4. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
88. IP1: Zero failed	Low		Inhibitible	IP one-point zero calibration has failed.	1. Fully disconnect then reconnect IP patient connection and re-zero IP transducer as prompted. 2. Replace IP patient connection. 3. If available, use a different IP channel on the MOVES[®] SLC[™]. 5. Power cycle MOVES[®] SLC[™]. 4. Return MOVES[®] SLC[™] for servicing.
89. IP1-ABP: Low systolic pressure	High			Systolic pressure of ABP measurement is less than or equal to the user set low systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
90. IP1-ABP: High systolic pressure	Medium			Systolic pressure of ABP measurement is greater than or equal to the user set high systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
91. IP1-ABP: Low diastolic pressure	Medium			Diastolic pressure of ABP measurement is less than or equal to 40 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
92. IP1-ABP: High diastolic pressure	Medium			Diastolic pressure of ABP measurement is greater than or equal to 100 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
93. IP1: Low CVP	High			CVP measurement is less than or equal to the user set low CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
94. IP1: High CVP	Medium			CVP measurement is greater than or equal to the user set high CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
95. IP1: Low ICP	High			ICP measurement is less than or equal to the user set low ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
96. IP1: High ICP	Medium			ICP measurement is greater than or equal to the user set high ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
97. IP2: Measurement error	Low			IP measurement contains values that are outside of the measurement range of the sensor.	1. Fully disconnect then reconnect IP patient connection and zero IP transducer as prompted. 2. Replace IP patient connection. 3. Check integrity of interface with patient. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
98. IP2 disconnect	Low		Inhibitible	IP patient connection is disconnected from system after previously being connected.	1. Fully disconnect then reconnect IP patient connection. 2. Replace IP patient connection. 3. Power cycle MOVES® SLC™. 4. Return MOVES® SLC™ for servicing.
99. IP2: Zero failed	Low		Inhibitible	IP one-point zero calibration has failed.	1. Fully disconnect then reconnect IP patient connection and re-zero IP transducer as prompted. 2. Replace IP patient connection. 3. If available, use a different IP channel on the MOVES® SLC™. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
100. IP2-ABP: Low systolic pressure	High			Systolic pressure of ABP measurement is less than or equal to the user set low systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitable	Cause	Action(s)
101. IP2-ABP: High systolic pressure	Medium			Systolic pressure of ABP measurement is greater than or equal to the user set high systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
102. IP2-ABP: Low diastolic pressure	Medium			Diastolic pressure of ABP measurement is less than or equal to 40 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
103. IP2-ABP: High diastolic pressure	Medium			Diastolic pressure of ABP measurement is greater than or equal to 100 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
104. IP2: Low CVP	High			CVP measurement is less than or equal to the user set low CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
105. IP2: High CVP	Medium			CVP measurement is greater than or equal to the user set high CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
106. IP2: Low ICP	High			ICP measurement is less than or equal to the user set low ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
107. IP2: High ICP	Medium			ICP measurement is greater than or equal to the user set high ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
108. IP3: Measurement error	Low			IP measurement contains values that are outside of the measurement range of the sensor.	1. Fully disconnect then reconnect IP patient connection and zero IP transducer as prompted. 2. Replace IP patient connection. 3. Check integrity of interface with patient. 4. Power cycle MOVES [®] SLC [™] . 5. Return MOVES [®] SLC [™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
109. IP3 disconnect	Low		Inhibitible	IP patient connection is disconnected from system after previously being connected.	1. Fully disconnect then reconnect IP patient connection. 2. Replace IP patient connection. 3. Power cycle MOVES® SLC™. 4. Return MOVES® SLC™ for servicing.
110. IP3: Zero failed	Low		Inhibitible	IP one-point zero calibration has failed.	1. Fully disconnect then reconnect IP patient connection and re-zero IP transducer as prompted. 2. Replace IP patient connection. 3. If available, use a different IP channel on the MOVES® SLC™. 4. Power cycle MOVES® SLC™. 5. Return MOVES® SLC™ for servicing.
111. IP3-ABP: Low systolic pressure	High			Systolic pressure of ABP measurement is less than or equal to the user set low systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
112. IP3-ABP: High systolic pressure	Medium			Systolic pressure of ABP measurement is greater than or equal to the user set high systolic BP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
113. IP3-ABP: Low diastolic pressure	Medium			Diastolic pressure of ABP measurement is less than or equal to 40 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
114. IP3-ABP: High diastolic pressure	Medium			Diastolic pressure of ABP measurement is greater than or equal to 100 mmHg.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
115. IP3: Low CVP	High			CVP measurement is less than or equal to the user set low CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
116. IP3: High CVP	Medium			CVP measurement is greater than or equal to the user set high CVP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
117. IP3: Low ICP	High			ICP measurement is less than or equal to the user set low ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
118. IP3: High ICP	Medium			ICP measurement is greater than or equal to the user set high ICP alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
119. High airway pressure warning	Medium			Peak inspiratory pressure (PIP) exceeds user set high airway pressure warning alarm limit without reaching user set high airway pressure release alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
120. High airway pressure release	Medium	Latching		Peak inspiratory pressure (PIP) exceeds user set high airway pressure release alarm limit.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
121. Low pressure end of inspire	Medium	Latching		Peak inspiratory pressure (PIP) does not reach 10 cmH ₂ O.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
122. Negative airway pressure	Low	Latching		Airway pressure is less than 0 cmH ₂ O when ventilating in any mode.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
123. Target pressure not reached	Medium			Peak inspiratory pressure (PIP) is $\pm$ 5 cmH ₂ O or greater from target when ventilating in pressure control mode.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
124. Expiratory obstruction	High			Positive end expiratory pressure exceeds user set positive end expiratory pressure (PEEP) by 5 cmH ₂ O or greater.	1. Ensure breathing circuit patient tubing is not kinked or obstructed. 2. Ensure expiratory pathway (ventilator cartridge and valve block) is not obstructed. 3. Replace breathing circuit including ventilator cartridge. 4. Return MOVES[®] SLC[™] for servicing.
125. Insp and exp pressure similar	High			Positive end expiratory pressure is within 10 cmH ₂ O of peak inspiratory pressure (PIP) when ventilating in IMV, SIMV and A/C modes.	1. Check for related alarms and follow actions if required. 2. Clinical intervention may be required.
126. Sample line occlusion	Medium			Pressure at O ₂ /CO ₂ sensor cells is below expected value compared to barometer (below normal resistance when sampling).	1. Ensure sample line has no kinks or obstructions. 2. Replace sample line. 3. Power cycle MOVES[®] SLC[™]. 4. Return MOVES[®] SLC[™] for servicing.
127. Sample line disconnected	Medium			Pressure at O ₂ /CO ₂ sensor cells is above expected value compared to barometer (above normal resistance when sampling).	1. Ensure sample line is fully connected to MOVES[®] SLC[™] and has no leaks. 2. Ensure that Nafion tubing and sample line filters are installed in the sample line 3. Replace sample line and filters. 4. Replace Nafion tubing. 5. Power cycle MOVES[®] SLC[™]. 6. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
128. Patient circuit occlusion	Medium			Inspired tidal volume is at least 15% less than the current user set point and the peak inspiratory pressure (PIP) exceeds user set high airway pressure release alarm limit.	1. Ensure breathing circuit patient tubing is not kinked or obstructed. 2. Ensure inspiratory pathway (ventilator cartridge and valve block) is not obstructed. 3. Replace breathing circuit including ventilator cartridge. 4. Return MOVES[®] SLC[™] for servicing.
129. Patient circuit disconnect	High			While running in volume control ventilator modes, expired CO ₂ partial pressure is less than 5 mmHg or the apnea alarm is active, and the minimum inspired pressure is not reached while the delivered inspired volume is at least 15% below target.	1. Ensure all breathing circuit connections (including patient tubing, Y-piece and sample line) are tight and no leaks are detected. 2. Ensure ventilator cartridge is properly connected to MOVES[®] SLC[™]. 3. Replace breathing circuit including ventilator cartridge. 4. Return MOVES[®] SLC[™] for servicing.
130. Concentrator air intake filter blocked	High			Pressure at intake filter pressure transducer exceeds expected value (above normal resistance when concentrator in use).	1. Check for and remove any obstructions from intake filter. 2. Replace intake filter. 3. Return MOVES[®] SLC[™] for servicing.
131. Concentrator air intake filter missing	Low			Pressure at intake filter pressure transducer is below expected value (below normal resistance when concentrator in use).	1. Ensure intake filter is completely secured to MOVES[®] SLC[™]. 2. Replace intake filter. 3. Return MOVES[®] SLC[™] for servicing.
132. Patient gas sampling system failure	Medium			Patient gas sampling pump flow is insufficient for proper patient gas monitoring.	1. Check for related alarms and follow actions if required. 2. Power cycle MOVES[®] SLC[™]. 3. Return MOVES[®] SLC[™] for servicing.

Alarm / Message	Priority	Locked / Latching	Inhibitible	Cause	Action(s)
133. Unexpected system reset	Medium	Locked	Inhibitible	System suffers unexpected power loss and power is restored within 30 seconds.	Return MOVES® SLC™ for servicing.
134. Safe Gas Mode	Low	Locked		System will revert to "Safe Gas Mode" if any of the following alarms are active while running in ventilation mode: 1. CO₂/O₂ sensors warming up 2. O₂ sensor fault 3. CO₂ sensor fault 4. Gas sensor pressure fault 5. Inspired O₂ low 6. Sample line disconnect 7. Sample line occlusion 8. Patient monitor failure	See actions for relevant alarms.
135. CO ₂ / O ₂ sensors warming up	Message	Locked		O ₂ and CO ₂ sensors are warming up.	Wait up to 5 minutes for system to become ready for use.

15.7 Safe Gas Mode



WARNING! THE MOVES[®] SLC[™] SHOULD NOT BE RUN CONTINUOUSLY IN SAFE GAS MODE. SAFE GAS MODE IS INTENDED FOR SHORT TERM USE ONLY TO COMPLETE TRANSPORTS.



NOTE: The MOVES[®] SLC[™] will enter Safe Gas Mode only when in Ventilate mode.

SGM is used whenever the system cannot rely on the oxygen (O₂) and/or carbon dioxide (CO₂) value that it is sensing. It is also used (when ventilating) if inspired CO₂ becomes too high due to an exhausted CO₂ scrubber cartridge.

SGM can occur for several reasons (each of which generates a unique alarm message):

- Sample line is occluded or disconnected for more than 60 seconds.
- O₂ or CO₂ sensor faults and cannot get a reading for more than 60 seconds.
- Operating the system for more than 60 seconds while the O₂ or CO₂ sensor is still reporting that it is warming up.
- CO₂/O₂ pressure sensor fault for more than 60 seconds.
- Patient monitoring board fault for more than 60 seconds.
- Inspired O₂ is below 19%.
- System is ventilating, and a high PiCO₂ alarm is active at high priority (PiCO₂ ≥ 10 mmHg). System exits SGM when there is no high PiCO₂ warning (PiCO₂ ≤ 6 mmHg), and no other SGM conditions are active.
- O₂ sensor readings are suspected of being biased high..

Standard Safe Gas Mode

If the Ventilator target O₂ setting is ≤ 40%:

- Safe Gas Mode will normally deliver approximately 5.5 to 6 LPM of 40–50% O₂.

High O₂ Safe Gas Mode

If the Ventilator target O₂ setting is > 40%:

- SGM delivers 2.5 LPM of 90% O₂ (nominal).



NOTE: If the concentrator is in a degraded mode, or in a fault mode when operating in Safe Gas Mode, the system will deliver at least 9.5 LPM of air in addition to any oxygen it can produce.

THIS PAGE DELIBERATELY LEFT BLANK.

16.0 Appendix

16.1 System Default Settings

Table 39: System Default Settings

Setting	Default
Alarms	All alarms and audio enabled
System Mode	Ventilator Mode PC-IMV (in Monitor Only)
Vent Mode	IMV
Vent Control	Pressure
Vt (not used by default mode PC-IMV)	500 ml
RR	10 B/M
PEEP	0-3 cmH ₂ O
I/E ratio	1:2
Insp Time (not used by default mode PC-IMV)	1 second
PC	20 cmH ₂ O
Pressure Support (not used by default mode PC-IMV)	disabled
O ₂	40%
NIBP Update	Manual
Suction	325 mmHg
ECG Lead Count	12
ECG Range	3.0 mV
ECG Sweep Speed	Fast
SpO ₂ Avg Time	8s
SpO ₂ Sensitivity Mode	Normal
SpHb Avg Window	Long
High airway pressure warning alarm limit	32 cmH ₂ O
High airway pressure release alarm limit	40 cmH ₂ O
Low expired PCO ₂ alarm limit	25 mmHg
High expired PCO ₂ alarm limit	55 mmHg
Low RR	5 B/M
High RR	30 B/M
Low systolic BP alarm limit	90 mmHg
High systolic BP alarm limit	180 mmHg
Low SpO ₂ alarm limit	85%
Low HR alarm limit	50 BPM
High HR alarm limit	120 BPM
Low CVP alarm limit	0 mmHg

Setting	Default
High CVP alarm limit	20 mmHg
Low ICP alarm limit	0 cmH ₂ O
High ICP alarm limit	27 cmH ₂ O

16.2 System Cleaning

The MOVES[®] SLC[™] is designed for easy maintenance. All exposed parts of MOVES[®] SLC[™] are corrosion resistant. The MOVES[®] SLC[™] device should be serviced after every 6 months of use or after every 12 months of storage.



CAUTION! DO NOT SUBMERGE THE MOVES[®] SLC[™] OR POUR CLEANING LIQUIDS OVER OR INTO THE MOVES[®] SLC[™].

The external body of the system can be cleaned using standard cleaning agents, excluding oxidizing agents. It is recommended that the external metal surfaces be wiped down with isopropyl alcohol during routine maintenance. Cables, NIBP cuffs and tubing can be cleaned with a disinfecting spray.

External surfaces of the MOVES[®] SLC[™] can be wiped clean with one of the following:

- Isopropyl Alcohol
- Chlorine Compounds*
 - Maximum Concentration: 1:10

**These compounds are diluted by volume in water.*

For recommended MOVES[®] SLC[™] accessories not labeled as single use, refer to the cleaning instructions provided by the manufacturer.

16.3 System Maintenance

16.3.1 Replacing Filters

There are two user-serviceable components in the MOVES[®] SLC[™] system:

1. Hydrocarbon Filter
2. Ventilator Filter

The hydrocarbon filter should be replaced when the system alarm indicates it is occluded. The ventilator filter should be replaced every six months. For information on replacing the hydrocarbon filter see *Section 9.6.3 Installing the Hydrocarbon Filter* on page 79.

Replacing the Ventilator Filter

1. Locate the Ventilator Driving Gas Inlet.



2. Locate the cap release.



3. Press down on the cap release while unscrewing the cap counterclockwise.



4. Replace the entire cap and filter and install a replacement cap and filter by pressing down on the cap release again while screwing in the replacement cap and filter clockwise.



16.3.2 General Maintenance



CAUTION! NO LUBRICANTS OTHER THAN THOSE RECOMMENDED BY THE MANUFACTURER SHALL BE USED ON THE MOVES[®] SLC[™].

However, regular maintenance and calibration should be carried out by authorized and qualified service personnel annually or after 100 operational hours of use. Systems should also be checked annually if not in use. In addition, if the MOVES[®] SLC[™] system is subjected to extremely rough handling or environmental stress, or sustains damage, it should be referred to authorized and qualified service personnel for inspection and / or repair.



WARNING! ONLY AUTHORIZED SERVICE AND MAINTENANCE PERSONNEL SHOULD REMOVE ANY COVERS FROM MOVES®. UNAUTHORIZED REMOVAL OF COVERS FROM MOVES® SLC™ MAY RESULT IN ELECTRIC SHOCK AND POSSIBLY DEATH, AND MAY DAMAGE THE SYSTEM COMPONENTS.

16.4 Accessories Maintenance

Reusable accessories should be regularly inspected for wear or damage. Inspect all cables and connections (especially the power cord) for signs of fraying or other damage. Keep accessories clean. Refer to original manufacturer's instructions for cleaning agents and procedures. Ensure that all gels and pastes are removed from electrode cables. Accessories that come in sanitary sealed packages should be inspected for damage to the sealed package before using. If the seal is broken, discard accessory. The label on the *package* of the Ventilator Cartridge contains an expiry date. Always check the expiry date on the Ventilator Cartridge package before using it to make sure that the cartridge has not expired. As well, monitor spare cartridges with regard to their remaining "shelf life".

16.5 Checking the Accuracy of the Temperature Probe

To check the accuracy of the temperature probe:

1. Pour a glass of warm water. Using an external temperature probe, measure the temperature of the water. The water temperature needs to be between 28°C (82.4°F) and 42°C (107.6°F).
2. Place the MOVES® SLC™ temperature probe in the water. Compare the temperature displayed by the MOVES® SLC™ with temperature measured by external temperature probe.
3. The temperature displayed by the MOVES® SLC™ should be +/-0.2F from the temperature measured by the external probe.

16.6 MOVES® SLC™ Specifications

16.6.1 Model Number

The documentation in this manual is for MOVES® SLC™ model number 122752.



CAUTION! OPERATION OF MOVES® SLC™ OUTSIDE OF SPECIFIED LIMITS MAY CAUSE INACCURATE RESULTS.

16.6.2 Physical Properties

Table 40: Physical Properties of MOVES® SLC™

Property	Specification	Notes
Unit Weight (lbs)	37.45	Excludes batteries, accessories, options, cables, etc.
Battery Weight, ea. (lbs)	3.25	
Unit Length (in)	33	
Unit Width (in)	5.5–6.5	
Unit Height (in)	10.25	
Unit Exterior Housing Material	Aluminum	
Operating Sound Level	< 70 dB	At a distance of 1 m level

Property	Specification	Notes
Standards Compliance	IEC 60601-1 IEC 60601-1-1 IEC 60601-1-2 IEC 60601-1-8 ISO 80601-2-13 IEC 60601-2-27 IEC 80601-2-30 IEC 60601-2-34 IEC 60601-2-49 ISO 80601-2-55 ISO 80601-2-61 ANSI/AAMI SP10 ANSI/AAMI EC13 ASTM E1112-00 BS EN 794-3 (2009) MIL-STD-810G JECETS	
Device Classification	Class II, CF Defibrillation Proof	
Screen	115.2 mm (w) x 86.4 mm (h)	Resolution of 640 x 480 pixels

16.6.3 Oxygen Concentrator Specifications

Table 41: Oxygen Concentrator Specifications of MOVES[®] SLC[™]

Property	Specification			Notes
Time to FiO ₂	Temp (°C)	FiO ₂ (%)	Time (min.)	All figures measured with a 1L test lung from cold start. FiO ₂ as measured on MOVES [®] SLC [™] UI with a patient will depend on the patient's FRC (washout time) and their oxygen consumption rate.
	22	50	2	
	22	70	4	
	22	80	6	
	22	85	10	
	0	50	3.5	
	0	70	8	
	0	80	12	
	0	85	18	
	-10	50	5	
	-10	70	11	
	-10	80	15	
	-10	85	22	

16.6.4 Ventilator Specifications

Table 42: Ventilator Specifications of MOVES[®] SLC[™]

Property	Specification
Tidal Volume	100–750 mL
Respiratory Rate	6–40 B/M

Property	Specification
Inspiratory/Expiratory Ratio	1:1 to 1:3
Inspiratory Resistance	6 cmH ₂ O (at 60 LPM)
Expiratory Resistance	6 cmH ₂ O (at 60 LPM)
Inspiratory Time	0.3–3 seconds
Peak Flow	60 LPM
Positive End Expiratory Pressure (PEEP)	0–20 cmH ₂ O
Positive Pressure Relief Valve	70 cmH ₂ O
Pressure Control (PC)	10–55 cmH ₂ O (over PEEP). PC (Pressure Control) = PIP (Peak Inspiratory Pressure) – PEEP (Positive End Expiratory Pressure)
Pressure Support Ventilation	Disabled, 5–40 cmH ₂ O
Trigger Sensitivity (Normal)	10 LPM (flow) or 3 cmH ₂ O (pressure) below PEEP
Trigger Sensitivity (Low)	15 LPM (flow) or 6 cmH ₂ O (pressure) below PEEP
External Oxygen Supply	15 LPM maximum (@ 2 psi minimum)
Modes	PC-IMV (default), VC-IMV, PC-A/C, VC-A/C, PC-SIMV, PC-SIMV+PS, VC-SIMV, VC-SIMV+PS, PSV
Standards Compliance	EN 794–3
Ventilator Circuit Compliance Including Cartridge	Approximately 0.7 mL/cmH ₂ O over the ventilator settable pressure range
Compressible Volume of Ventilator and Cartridge	1350 mL (Note: System compensates for compressible volume)

16.6.5 Suction Specifications

Table 43: Suction Specifications of MOVES[®] SLC[™]

Property	Specification
Suction Vacuum	–100 to –325 mmHg
Free Flow Rate	20 LPM

16.6.6 Electrical Characteristics

Table 44: Electrical Specifications of MOVES[®] SLC[™]

Property	Specification	Notes
External Power	100–240 VAC, 50–60 Hz, 5.5 A max.	
Max Current Output	28 VDC, 14.3 A max.	
Battery Type	25.9 V lithium polymer	
Charge Time (per set 2 batteries)	2.5 hrs when the system is idle	

Property	Specification	Notes
Battery Life (per set 2 batteries)	Minimum 2.5 hrs on two fully charged new batteries at 101 kPa and 21°C. Typical > 4 hours with ventilator and monitors on, running concentrator with ¼ duty cycle.	Battery run time is highly dependent on use of the oxygen concentrator and / or suction.

16.6.7 Environmental Specifications

Table 45: Environmental Specifications of MOVES[®] SLC[™]

Variable	Storage Condition	Operating Condition
Temperature	–14°F to 140°F (–26°C to 60°C) – system & batteries	–14°F to 129°F (–26°C to 54°C). Notes: 1. MOVES [®] SLC [™] can be taken from any operating temperature into extreme cold of –26°C without affecting operation. It can be cold started at –26°C (–14°F) on AC power. However, MOVES [®] SLC [™] can not be cold started below –20°C (–4°F) on batteries. 2. When starting up MOVES [®] SLC [™] in cold temperatures, the system may take a longer time to reach required concentration. See table in Section 13.6.3. 3. NEVER CHARGE BATTERIES IN AMBIENT TEMPERATURES BELOW 32°F (0°C) OR ABOVE 104°F (40°C).
Relative Humidity	15% to 95% non-condensing	Same
Altitude	0–18,000 ft. (5,500m)	Same
Water Resistance	MIL-STD-810G Method 506.4 for blowing rain	Same
Blowing Sand and Dust Resistance	MIL-STD-810G Method 510.4 Procedure II	Same

16.7 Patient Monitoring Specifications

16.7.1 Heart Rate Monitoring Specifications

Table 46: Heart Rate Monitoring Specifications of MOVES[®] SLC[™]

Item	Specification
Source	Auto-detect in the priority sequence of ABP, SpO ₂ , ECG
Range	30–250 bpm for ABP and ECG, 30–239 bpm for SpO ₂
Accuracy	±1% Full Scale (under stationary operation) ±5 BPM (under continuous vibration)

Item	Specification
Pacemaker Pulse Rejection Capability	Pacemaker pulses may be detected by the ECG heart rate monitor and included in its calculation, depending on the type and model of pacemaker detected by the heart rate monitor.
Fixed Delays Due to Signal Processing	ABP: Pulse heart rate is calculated from the previous 6 beats. SpO ₂ : Pulse heart rate is calculated based on the SpO ₂ average time selected on the Advanced screen (2 –16 seconds, default 8 seconds). ECG: Heart rate is calculated from the previous 8 beats.
Alarm Condition Delay (onset of condition to internal realization) NOTE: "Alarm Condition Delay" derives its value from "Fixed Delays Due To Signal Processing" + 100 ms	ABP: 6 beats + 100 ms SpO ₂ : SpO ₂ average time + 100 ms ECG: 8 beats + 100 ms
Alarm Signal Generation Delay (realization to display)	Less than 200 ms
Operating Mode That May Affect Alarm Generation	None

16.7.2 Temperature Monitoring Specifications

Table 47: Temperature Monitoring Specifications of MOVES[®] SLC[™]

Item	Specification
Range	28°C to 42°C (82.4°F to 107.6°F)
Accuracy	±0.1°C

16.7.3 Airflow Monitoring Specifications

Table 48: Airflow Monitoring Specifications of MOVES[®] SLC[™]

Item	Specification
Inspiratory/Expiratory Flow Range	–60 to 60 LPM
Repeatability of Inspiratory/Expiratory Flow Measure	±0.5% (% of reading)
Airway Pressure (Paw) Range	–5 to 100 cmH ₂ O
Accuracy of Airway Pressure Measure	±2 cmH ₂ O + 8% of reading

16.7.4 CO₂ Monitoring Specifications

Table 49: CO₂ Monitoring Specifications of MOVES[®] SLC[™]

Item	Specification
Range	0 to 10% by Volume
Accuracy	± 1.0% Absolute
Rise Time	215 ms (10-90%) at 200 ml/min
Response time of gas sample readings	< 4 seconds
Flow Rate	250 ml/min ± 50 ml/min

Item	Specification
Standards Compliance for CO ₂ Analyzer Used	ISO 80601-2-55: Medical electrical equipment – Particular requirements for the basic safety and essential performance of respiratory gas monitors

16.7.5 O₂ Monitoring Specifications

Table 50: O₂ Monitoring Specifications of MOVES® SLC™

Item	Specification
Range	5 to 100% by Volume
Resolution	0.02%
Accuracy	±4% Absolute
Rise Time	150 ms (10-90%) at 150 ml/min
Response time of gas sample readings	< 4 seconds
Flow Rate	250 ml/min ± 50 ml/min
Operating Mode That May Affect Alarm Generation	None
Standards Compliance for O ₂ Analyzer Used	ISO 80601-2-55: Medical electrical equipment – Particular requirements for the basic safety and essential performance of respiratory gas monitors

16.7.6 ECG Specifications

Table 51: ECG Specifications of MOVES® SLC™

Item	Specification
Number of Leads	12
Lead View	Standard 12 lead with Wilson chest lead placement
Input Impedance	> 10 MOhm
Input Range	> 10 mvpp
Input Range (DC)	> 300 mV
Sensitivity (mV/cm)	1.4, 2.8, 5.5
Filtering	50 Hz , 60 Hz, or no notch filter, optional 15-30 Hz EMG filter
Frequency Response	0.3–70 Hz
Pulse Detection	30 – 250 bpm + - 1%, + - 1 Digit, 8 beat averaging
Defibrillator Protection	Yes
ST Analysis	None
Pacer Detection	None
Standard Complied to	IEC 60601–2–27
Maximum T-Wave Amplitude	0.6 mV
T-Wave Rejection	The tall T-wave rejection is always on. The max T-wave rejection is 0.6 mV when R-wave is 1 mV.

Item	Specification	
Indicated Heart Rate (after a 20-second equipment stabilization period)		Rhythm
	A1 (80 BPM)	78
	A2 (60 BPM)	58
	A3 (120 BPM)	118
	A4 (90 BPM)	88
Response Time to Heart Rate Meter	Longest Time	< 10 seconds
	80 BPM–120 BPM	
	Longest Time	< 15 seconds
	80 BPM–40 BPM	
Time to Alarm for Tachycardia	B1 – Regular Amplitude	6 seconds
	B1 – Half Amplitude	6.2 seconds
	B1 – Double Amplitude	± 6 seconds
	B2 – Regular Amplitude (max)	5 seconds
	B2 – Half Amplitude (max)	4 seconds
	B2 – Double Amplitude (max)	4.2 seconds
Pacemaker Pulse Representation	Pacemaker pulses will be displayed in the MOVES® SLC™ ECG waveform display and can affect heart rate, particularly in Hi-F mode (as an overestimation).	

16.7.7 NIBP Specifications

Table 52: NIBP Specifications of MOVES® SLC™

Item	Specification
Measurement Cycles	Stat, 1, 2, 3, 4, 5, 10, 15 minutes
Max Allowable Cuff Pressure	300 mmHg
Adult Range	Systolic: 40–260 Diastolic: 20–200
Resolution	1 mmHg
Accuracy	± 5 mmHg Average with STD of 8 mmHg
Calibration	The cuff pressure transducer should be verified once every 12 months.

16.7.8 Invasive Pressure Specifications

Table 53: Invasive Pressure Specifications of MOVES® SLC™

Item	Specification
Channels	3
Transducer Sites	ABP, CVP or ICP

Item	Specification
Pressure Range	ABP: -10 to 300 mmHg CVP: -10 to 300 mmHg ICP: -14 to 408 cmH ₂ O
Temperature: Operating Storage	15° to 40°C (57°F to 104°F) -25° to 70°C (13°F to 158°F)
Accuracy	±4 mmHg or 4% of reading whichever is greater

16.7.9 Pulse Oximetry Specifications

Table 54: Pulse Oximetry Specifications of MOVES[®] SLC[™]

Item	Specification
Method	Multiple visible and infrared LEDs (500 to 1400 nm)
Fixed Delays Due to Signal Processing	Current SpO ₂ is calculated from previous 2–16 seconds depending on “SpO ₂ Average Time” setting in Advanced section of Setup Screen. Default: 8 seconds.
Alarm Condition Delay (onset of condition to internal realization)	Fixed delay due to signal processing + 100 ms
Alarm Signal Generation Delay (realization to display)	Less than 200 ms

Equipment Response Time

Table 55: Pulse Oximetry Equipment Response Time

SpO ₂ Values	Average	Latency
Standard/Fast Averaged SpO ₂	Choice of 2–4, 4–6, 8, 10, 12, 14, or 16 seconds. Default: 8 seconds	2 beats
Pulse Rate Values	Average	Latency
Standard/Fast Averaged Pulse Rate	Choice of 2–4, 4–6, 8, 10, 12, 14, or 16 seconds. Default: 8 seconds	2 beats

16.7.10 Drift in Sensing Accuracy

A drift of less than 0.4% in the oxygen reading can be expected over a 6 hour duration when external conditions are held constant.

A drift of less than 0.3% in the carbon dioxide reading can be expected over a 6 hour duration when external conditions are held constant.

16.7.11 Specifications of Masimo Rainbow SET[®] Pulse CO-Oximeter

Measurement Range

Measurement	Display Range
SpO ₂ (Oxygen Saturation)	0% to 100%

Measurement	Display Range
SpMet (Methemoglobin)	0% to 100%
SpCO (Carboxyhemoglobin)	0% to 100%
SpHb (Hemoglobin)	0 g/dL to 25.0 g/dL
SpOC (Oxygen Content)	0 ml of O ₂ /dL to 35 ml of O ₂ /dL of blood
PR (Pulse Rate)	25 bpm to 239 bpm, >239 displayed when PR is 240–260 bpm
PI (Perfusion Index)	0% to 20% (bar graph beside SpO ₂)

Accuracy

Oxygen Saturation Accuracy [1]	
No Motion	60% to 80%
Adults	± 3%
No Motion [2]	70% to 100%
Adults	± 2%
Motion [3]	70% to 100%
Adults	± 3%
Low Perfusion [4]	
Adults	± 2%
Pulse Rate Accuracy	
Pulse rate range	25 bpm to 240 bpm
No Motion	
Adults	± 3 bpm
Motion [4]	
Adults	± 5 bpm
Low Perfusion	
Adults	± 3 bpm
Carboxyhemoglobin Accuracy [1]	
Adults	1% to 40% ± 3%

Methemoglobin Saturation Accuracy [1]	
Adults	1% to 15% $\pm$ 1%
Total Hemoglobin Accuracy [2]	
Adults	8 g/dL to 17 g/dL $\pm$ 1 g/dL

Notes:

1. SpO₂, SpCO, and SpMet accuracy was determined by testing on healthy adult volunteers in the range 60% to 100% SpO₂, 0% to 40% SpCO, and 0% to 15% SpMet against a laboratory CO-Oximeter. SpO₂ and SpMet accuracy was determined on 16 neonatal NICU patients ranging in age from 7 days to 135 days old and weighting between 0.5 kgs and 4.25 kgs. Seventy-nine (79) data samples were collected over a range of 70% to 100% SaO₂ and 0.5% to 2.5% HbMet with a resultant accuracy of 2.9% SpO₂ and 0.9% SpMet. Contact Masimo for testing specifications.

2. SpHb accuracy has been validated on healthy adult male and female volunteers and on surgical patients with light to dark skin pigmentation in the range of 8 g/dL to 17 g/dL SpHb against a laboratory CO-Oximeter. The variation equals plus or minus one standard deviation which encompasses 68% of the population. The SpHb accuracy has not been validated with motion or low perfusion.

Resolution

Parameter	Step Size
%SpO ₂	1%
%SpCO	1%
%SpMet	0.1%
SpHb g/dL	0.1 g/dL
SpOC mL/dL	1 mL O ₂ /dL blood
Pulse Rate	1 beats per minute

16.8 General Accessories Specifications

Table 56: General Accessories Specifications of MOVES[®] SLC[™]

Item	Category	Specification
Batteries	Operating Temperature:	-26°C to 54°C (-14°F to 129°F) (cold start above -20°C)
	Storage Temperature:	-26°C to 60°C (-14°F to 140°F)
	Charging Temperature:	0°C to 40°C (32° to 104°F)

Item	Category	Specification
Hydrocarbon/Particulate Filter (P/N 100915)	Effective Filtration Against:	GME Organic Vapor, Chlorine, Sulfur Dioxide, Chlorine Dioxide, Hydrogen Chloride, Hydrogen Sulfide, Ammonia, Methylamine, Formaldehyde, Hydrogen Fluoride: 99.97% effective against all particulate aerosols
NIBP Cuffs (ALL 4)	Temperature: Operating Storage	0°C to 40°C (32°F to 104°F) -34°C to 70°C (-29.2°F to 158°F)
Temperature Probe (Reusable/Autoclavable)	Temperature:	Interchangeable $\pm 0.1^{\circ}\text{C}$, 25°C to 45°C per EN 12470 Tested to $\pm 0.1^{\circ}\text{C}$, 0°C to 70°C for laboratory use
	Autoclave:	Withstands 100 autoclave cycles to 121°C Withstands 50 autoclave cycles to 134°C
Suction Canister		800 mL capacity

16.8.1 Approved Masimo Pulse Oximeter Accessories

The following Masimo accessories may be used in conjunction with the MOVES[®] SLC[™] pulse oximeter. See the respective sensor instructions for detailed information regarding specified sensor use.

Table 57: Approved Masimo Pulse Oximeter Accessories

Accessory	Description
Rainbow [®] DC-3 SC-360	Adult Reusable Finger Sensor for SpHb, SpO ₂ & SpMet, 3 ft., 60 hours
Rainbow [®] DCI-dc3	Adult Reusable Finger Sensor for SpCO, SpO ₂ & SpMet, 3 ft.
Rainbow [®] RC-4	Rainbow [®] Adhesive Sensor Patient Cable, 4 ft.
Rainbow [®] R1 25	Adult SpO ₂ , SpHb & SpMet Adhesive Sensor NOTE: Requires Rainbow [®] RC-4 Patient Cable
Rainbow [®] R25	Adult SpO ₂ , SpCO & SpMet Adhesive Sensor NOTE: Requires Rainbow [®] RC-4 Patient Cable
Red LNC-10	LNCS Patient Cable for SpO ₂ , 10 ft.
LNCS DCI [®]	Adult Reusable SpO ₂ Finger Sensor, 3 ft. NOTE: Requires Red LNC-10 Patient Cable
LNCS TC-I	Reusable SpO ₂ Ear Sensor, 3 ft NOTE: Requires Red LNC-10 Patient Cable
LNCS Amtx-3	Adult Adhesive Sensor, 3 ft. NOTE: Requires Red LNC-10 Patient Cable

16.8.2 Masimo Pulse Oximeter Accessories Specifications

Table 58: Masimo Pulse Oximeter Accessories Specifications

Sensor	Description	Preferred Application Site	Masimo P/N	Weight Range	SpO2 Accuracy		PR Accuracy		Low Perfusion Accuracy		SpCO/Hb Accuracy	SpMet Accuracy
					No Motion	Motion	No Motion	Motion	SpO2	PR	No Motion	No Motion
Rainbow DC-3 SC-360	Adult Reusable Finger Sensor for SpHb, SpO2, & SpMet, 3 ft., 60 hours	Finger	2241	>30 kg	60-80% ±3% 70-100% ±2%	±3%	±3 bpm	±5 bpm	±2%	±3 bpm	(SpHb) ±1 g/dL	± 1%
Rainbow DCI-dc3	Adult Reusable Finger Sensor for SpO2, SpCO, & SpMet, 3ft.	Finger	2201	>30 kg	60-80% ±3% 70-100% ±2%	±3%	±3 bpm	±5 bpm	±2%	±3 bpm	(SpCO) ±3%	± 1%
LNCS DCI®	Adult Reusable Finger Sensor	Finger	1863	>30 kg	±2%	±3%	±3 bpm	±5 bpm	±2%	±3 bpm	N/A.	N/A.
LNCS TC-I	Reusable Ear Sensor	Ear Lobe	1895	>30 kg	±3.5%	N/A.	±3 bpm	n.a.	±3.5%	±3 bpm	N/A.	N/A.
LNCS Adtx-3	Adult Adhesive Sensor, 3 ft.	Finger	2317	>30 kg	±2%	±3%	±3 bpm	±5 bpm	±2%	±3 bpm	N/A.	N/A.
Rainbow R1 25	Adult Adhesive Sensor for SpHb, SpO2, & SpMet	Finger	2416	>30 kg	60-80% ± 3% 70-100% ± 2%	±3%	±3 bpm	± 5 bpm	± 2%	± 3 bpm	(SpHb) ± 1 g/dL	± 1%

Sensor	Description	Preferred Application Site	Masimo P/N	Weight Range	SpO2 Accuracy		PR Accuracy		Low Perfusion Accuracy		SpCO/Hb Accuracy	SpMet Accuracy
					No Motion	Motion	No Motion	Motion	SpO2	PR	No Motion	No Motion
Rainbow R25	Adult Adhesive Sensor for SpO2, SpCO, & SpMet	Finger	2221	>30 kg	60-80% ± 3% 70-100% ± 2%	±3%	±3 bpm	±5 bpm	± 2%	±3 bpm	(SpCO) ±3%	±1%

16.9 Electromagnetic Conformity Information

16.9.1 IEC 60601-1-2:2007 (Ed 3.0) Table 1 Requirements

MOVES[®] SLC[™] is intended for use in the electromagnetic environment specified below. The customer or the user of MOVES[®] SLC[™] should assure that it is used in such an environment.

Table 59: 15.9.1 IEC 6060-1-2:2001 (Ed 2) Table 201 Requirements

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	MOVES [®] SLC [™] uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	MOVES [®] SLC [™] is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: MOVES [®] SLC [™] is intended for use by healthcare professionals only. MOVES [®] SLC [™] may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating MOVES [®] SLC [™] or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

16.9.2 IEC 60601-1-2:2007 (Ed 3.0) Table 202 Requirements

MOVES[®] SLC[™] is intended for use in the electromagnetic environment specified below. The customer or the user of MOVES[®] SLC[™] should assure that it is used in such an environment.

Table 60: 5.2.2.1f IEC 60601-1-2:2007 (Ed 3.0) Table 2 Requirements

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air with documented necessary	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 610004-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	< 5% UT (> 95% dip in UT) for 0,5 cycle 40 % UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles No anomalies < 5% UT (> 95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. End user shall ensure charged batteries are installed in the Equipment.
Power frequency (50 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

16.9.3 IEC 60601-1-2:200 (Ed 3.0) Table 203 Requirements

MOVES[®] SLC[™] is intended for use in the electromagnetic environment specified below. The customer or the user of MOVES[®] SLC[™] should assure that it is used in such an environment.

Table 61: 5.2.2.2 IEC 60601-1-2:2007 (Ed 3.0) Table 3 Requirements

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	10 Vrms 150 kHz to 80 MHz in ISM Bands	10 V	Portable and mobile RF communications equipment should be used no closer to any part of MOVES[®] SLC[™], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{12}{10} \right] \sqrt{P}$ $d = \left[\frac{12}{10} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{23}{10} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$...where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			

16.9.4 IEC 60601-1-2:2007 (Ed 3.0) Table 5 Requirements

MOVES[®] SLC[™] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of MOVES[®] SLC[™] can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and MOVES[®] SLC[™] as recommended below, according to the maximum output power of the communications equipment.

Table 62: 5.2.2.2 IEC 60601-1-2:2007 (Ed 3.0) Table 5 Requirements

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and MOVES [®] SLC [™]				
Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)			
	150 kHz to 80 MHz outside ISM bands $d = [\frac{3.5}{V_1}] \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d = [\frac{12}{V_2}] \sqrt{P}$	80 MHz to 800 MHz $d = [\frac{12}{E_1}] \sqrt{P}$	800 MHz to 2.5 GHz $d = [\frac{23}{E_1}] \sqrt{P}$
0.01				
0.1				
1				
10				
100				
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 to 40.70 MHz. NOTE 3: An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into the patient area. NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Index

Accessories, General
specifications, 220

Accessories, Patient Monitoring
connecting to MOVES[®] SLC[™], 73

Accessories, Suction
installing, 85

Address
authorized representative, iii
manufacturer, iii
notified body, iii

Advanced Screen
accessing, 122
additional info on, 126
image of, 123
options available, 124
selecting options on, 123

Adverse Effects
associated with respiratory gas monitoring, 22

Airflow Monitoring
specifications, 215

Alarm Audio Pause Button
about, 103

Alarm Limits
list of values & defaults, 139

Alarm Limits List
changing values, 139

Alarm Limits Screen
overview, 138

Alarm On / Off List
toggling value, 141

Alarm On / Off Screen
overview, 140

Alarm Priorities & Characteristics
explained, 175

Alarm Queue
overview, 176

Alarm Status LEDs
explained, 175

Alarm View Button
about, 113

Alarm, High Priority Communication Failure
description, 176

Alarms
about locked, 176
conditions & causes, 179–204
overview, 175
table of conditions & causes, 179
turned off, 178
visual indicators, 104

Alarms, Inhibitible
about, 177

Alarms, Latching
about, 178

Alarms, Standard
types & descriptions, 175

Asterisk Beside Sensor Values
meaning of, 22

Battery Status
table of icons, 114

Batteries
checking charge, 90
disposal of, 51
handling of, 51
inspecting, 91
installing, 92
storage cautions & warnings, 31
storage of, 31

Batteries, Charging
notes, warnings & cautions about, 91

Batteries, Remote Screen
installing, 143

Battery Charger
preparing, 95

Battery Indicator, on Remote Screen
about, 147

Battery Status
icons indicating, 113

Battery Status, Remote Screen
indicators listed & explained, 148

Calibration
of sensors, 21

Carboxyhemoglobin (SpCO)
general description for, 29

Cautions
electrical, 48
general, 45
patient-specific, 48

Cautions, Accessories
MOVES[®] SLC[™], 53

Channel
zeroing, 136

Check Button
about, 104

Circuit, Breathing
overview, 20

Clamps
attaching to MOVES[®] SLC[™], 70

Classification
of equipment, 2

Cleaning
of system, 208

Close Button
about, 104

CO₂ Monitoring
specifications, 215

Compliance, Regulatory
standards complied to, 18
symbols indicating, 17

Concentrator, Oxygen
overview, 20

Connecting Status
shown on remote screen, 155

Consumables
MOVES[®] SLC[™], 65

Contents	Guide, Quick Start Ventilate, 3
of MOVES® SLC™ system, 54	further info on procedures, 10
Copyright	Heart Rate Monitoring
notice of, 1	specifications, 214
Data Views	Hydrocarbon Filter
changing, 105	installing, 79
display of, 105	Icon, Alarm Off
Declaration of Conformity	about, 111
notice about, 2	Icon, Audio Off
Display Validity	about, 112
confirming, 106	Icon, Audio Pause
Disposal	about, 112
of batteries, 51	Icon, Battery
Ear Clip Oximeter	about, 111
attaching, 161	Icon, Battery Status
ECG	about, 113
specifications, 216	Icon, NIBP
ECG Artifacts	about, 112
reducing, 166	Icon, Number of Limits Changed
ECG Cables	about, 111
colour coding & naming conventions, 164	Icon, Plug
ECG Electromyography (EMG) Filter	about, 111
enabling/disabling, 125	Icon, Suction
ECG Screen	about, 112
overview, 136	IEC 60601-1-2
Electrical Specifications	200 (Ed 3.0) Table 203
of MOVES® SLC™, 213	requirements, 225
Electrocardiogram (ECG)	2007 (Ed 3.0) Table 1
overview, 23	requirements, 223
Electrode Package	2007 (Ed 3.0) Table 202
expiry date & expiry time, 166	requirements, 223
Electrodes, ECG	2007 (Ed 3.0) Table 5
correct placement, 164	requirements, 226
expiry date of, 166	Initializing
Electrodes, Heart Rate	message displayed on remote screen, 155
placement of, 163	Interference, Radio
Electromagnetic Conformity	caution about, 51
info about, 223	Introduction
EMC	to operator's manual, 13
notice about, 1	Invasive Pressure
Emergency Operation	overview, 23
preparing for, 50	source buttons explained, 136
Environment, Operating	specifications, 217
of MOVES® SLC™, 15	Main Screen
Environmental Specifications	items displayed & descriptions, 130
of MOVES® SLC™, 214	overview, 126
Equipment	Maintenance
classification of, 2	general, 210
Filter, Ventilator	of accessories, 211
replacing, 209	of system, 208
Filters	Manufacturer
replacing, 208	address, iii
Finger Clip Oximeter	Masimo Pulse Oximeter
attaching, 161	specifications, 218
Glossary	Measurements During Patient Motion
of terms and abbreviations, 13	of SpCO, SpMet & SpHb, 29
Graphs	Messages
about, 136	about locked, 176
Graphs, System	Messages, in Manual
list of available & parameters, 172	about, 33
overview, 171	Methemoglobin (SpMet)
	general description for, 29

- Model Number of MOVES[®] SLC[™], 211
- Monitor Only Mode, on Setup Screen image of, 118
- Monitoring, Patient intended use, 15 types described, 23
- Monitoring, Respiratory Gas overview, 21
- Monitoring, Successful for SpCO, 29 for SpHb, 28 for SpMet, 29 for SpO₂, PR & PI, 26
- MOVES[®] SLC[™]
 - attaching clamps to, 70
 - attaching shoulder strap to, 68
 - connecting patient monitoring accessories, 73
 - connecting remote screen to, 151
 - connecting to power supply / battery charger, 96
 - consumables used with, 65
 - external components, 19
 - intended use, 15
 - model number of, 211
 - physical properties of, 211
 - preparing for activation, 90
 - shutdown of, 168
 - system contents, 54
 - system orientation, 19
 - system overview, 19
- MRI notice about, 1
- Next Button about, 104
- NIBP specifications, 217
- NIBP Control Button about, 103
- Non-Invasive Blood Pressure (NIBP) overview, 23
- Notice
 - copyright, 1
 - declaration of conformity, 2
 - EMC, 1
 - MRI, 1
 - patent, 1
 - regulatory, 1
 - trademark, 1
- Notified Body address, iii
- O₂ Monitoring specifications, 216
- O₂ Supplement, on Setup Screen image of, 118
- Operator, Intended of MOVES[®] SLC[™], 15
- Operator's Manual introduction, 13
- Options, Main Screen selection order for, 128
- Oxygen Concentrator specifications, 212
- Oxygen Saturation functional, 26
- Oxygen Saturation (SpO₂) general description for, 26
- Oxygen, Supplementary delivering, 83 overview, 20
- Panel Buttons on remote screen, 157
- Patent notice of, 1
- Patient connecting, 161
- Patient Monitoring inverted display of values, 135 specifications, 214
- Patient, Intubated connecting, 161
- Patient, Spontaneously Breathing connecting, 161
- Perfusion Index (PI) general description for, 26
- Physical Properties of MOVES[®] SLC[™], 211
- Population, Intended of MOVES[®] SLC[™], 15
- Power Control Button about, 102
- Power Supply preparing, 95
- Power System about MOVES[®] SLC[™], 24
- Pressure Display (PC, PEEP & PIP) explained, 135
- Previous Button about, 104
- Pulse CO-Oximetry vs. drawn whole blood measurements, 28
- Pulse Oximeter
 - about, 25
 - approved accessories, 221
 - indications for use, 25
 - key features, 25
 - overview, 25
 - technology overview, 25
- Pulse Oximeter Accessories specifications, 222
- Pulse Oximeter Sensor avoiding inaccurate readings, 162 warnings, 163
- Pulse Oximetry overview, 23 specifications of, 218
- Pulse Rate (PR) general description for, 26
- Quick Start Ventilate Guide, 3
- Rainbow Pulse CO-Oximetry Technology principles of, 26
- Remote Screen
 - alarm indicators, 158
 - connecting to MOVES[®] SLC[™], 151

- connecting to wall power, 149
- determining software version, 160
- no System Test screen, 155
- panel buttons, 157
- user interface, 154
- using, 143
- when first connecting, 154
- Remote Screen, Graphs
 - independent of main screen, 158
- Remote Screen, Navigating
 - via finger touch, 158
 - via panel buttons, 157
- Representative, Authorized
 - address, iii
- Requirements
 - IEC 60601-1-2
 - 200 (Ed 3.0) Table 203, 225
 - 2007 (Ed 3.0) Table 1, 223
 - 2007 (Ed 3.0) Table 202, 223
- Resources, Loading
 - shown on remote screen, 154
- Respiratory Gas Monitoring
 - adverse effects associated with, 22
- Safe Gas Mode
 - explanation of & reasons for, 205
 - high O₂, 205
 - standard, 205
- Safety
 - electrical, 50
 - general, 49
- Screen Button
 - about, 104
- Screen Dim Button
 - about, 102
- Screen Select Button
 - about, 112
- Screen, Display Orientation
 - adjusting, 99
- Sensing Accuracy
 - drift in, 218
- Sensors
 - calibration of, 21
- Setting
 - modifying, 106
- Settings
 - changing, 105
 - display of, 105
- Settings, Default
 - of system, 207
- Setup Screen
 - about, 116
 - changing settings, 116
 - options listed & described, 119
- Shoulder Strap
 - attaching to MOVES[®] SLC[™], 68
- Shutdown
 - of MOVES[®] SLC[™], 168
- Signal Extraction Technology (SET[®])
 - about, 25
- Specifications
 - airflow monitoring, 215
 - CO₂ monitoring, 215
 - ECG, 216
 - electrical, 213
 - environmental, 214
 - general accessories, 220
 - heart rate monitoring, 214
 - invasive pressure, 217
 - Masimo, 218
 - NIBP, 217
 - O₂ monitoring, 216
 - of oxygen concentrator, 212
 - of suction, 213
 - of ventilator, 212
 - patient monitoring, 214
 - pulse oximetry, 218
 - temperature monitoring, 215
- SpHb Average Window
 - selecting & description, 125
- SpO₂ Avg (Average) Time
 - selecting & description, 124
- SpO₂ Sensitivity Modes
 - selecting and description, 125
- SpOC
 - general description for, 29
- Standards, Regulatory
 - complied to, 18
- Startup Sequence
 - instructions, 106
- Status Bar
 - about items displayed, 111
- Storage
 - of batteries, 31
 - of MOVES[®] SLC[™], 31
- Suction
 - intended use, 15
 - overview, 22
 - specifications, 213
 - using, 168
- Suction Control Button
 - about, 104
- Symbols
 - operational, 34
 - regulatory compliance, 17
- Symbols, in Manual
 - explained, 33
- Symbols, Product
 - on labels, 34
- Symbols, Warning
 - on labels, 34
- System
 - trade name of, 2
- System State Indicator
 - about, 112
- System Status
 - visual indicators, 104
- System Test Screen
 - about, 115
- Temperature Monitoring
 - overview, 23
 - specifications, 215
- Temperature Probe
 - checking accuracy of, 211
- Total Hemoglobin (SpHb)

- general description for, 28
- Trade Name
 - of system, 2
- Trademark
 - notice about, 1
- Transducer
 - extreme warnings, 167
- Transducer Channel
 - zeroing, 167
- Transducer, IP
 - zeroing pressure in, 167
- Trends, System
 - list of & parameters, 173
 - overview, 171
- User Interface (UI)
 - controls & functions, 102
- Ventilate Mode, on Setup Screen
 - image of, 117
- Ventilator
 - overview, 20
 - specifications, 212
- Ventilator Breathing Circuit
 - installing, 80
- Ventilator Cartridge
 - installing, 77
 - using, 75
- Warnings
 - electrical, 39
 - general, 37
 - patient-specific, 39
 - pulse oximeter, 41
- Warnings, Accessories
 - MOVES[®] SLC[™], 53
- Warnings, Extreme
 - transducer, 167